



Wholesale Supplier Cardiovascular_instrument A to Z Hospital Equipment and Surgical Products **B to C**



Cardiovascular

Scissors	328 - 356
Needle Holders	357 - 370
Forceps	371 - 377
Titanium Forceps and Clamps	378 - 389
Forceps and Clamps	390 - 399
Vascular Clamps	400 - 481
Bulldog and Micro Vessel Clamps	482 - 490
Vessel Instruments	491 - 493
Specialty CV/Thoracic Instruments and Retractors ..	494 -
507 Retractors	508 -
517 Rib Spreaders	518 -
	533

8275140877

9822688071

info@amramedtech.com



Sr No.29/1/2, Adarsh Nagar,
Uruli Devachi, Hadapsar, Pune
Branches :- All Maharashtra




amramedtech.com

CARDIOVASCULAR SCISSORS



Vessel Scissors

Catalog Number	Tip	Blade	Length (in.)
BC034R	B/B	1-PROBE-PT	5 3/4





MAYO Operating Scissors


Catalog Number	Tip	Blade	Length (in.)
BC466R	S/S	STR	5 1/2
BC467R	S/S	CVD	5 1/2



MAYO SuperCut Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC924R	B/B	STR	6 1/2
BC925R	B/B	CVD	6 1/2








CARDIOVASCULAR SCISSORS



CHURCH Scissors, Beveled

Catalog Number	Tip	Blade	Length (in.)
MB618R	B/B	STR	7
MB619R	B/B	CVD	7





DEBAKEY Arteriotomy Scissors

Catalog Number	Tip	Blade	Length (in.)
BC629R	B/B	CVD	6 7/8





THOREK Scissors

Catalog Number	Tip	Blade	Length (in.)
BC636R	B/B	HEAVY-CVD	7



SCHMIEDEN-TAYLOR Dura Scissors

Catalog Number	Tip	Blade	Length (in.)
BC660R	1 BALL-PT	ANG	6 1/8



CARDIOVASCULAR SCISSORS

SATINSKY Vena Cava Scissors


Catalog Number	Blade	Length (in.)
BC641R	CVD	9 5/8



BC641R

1/2

1/2



BC641R

1/1


STRULLY Scissors

Catalog Number	Tip	Blade	Length (in.)
BC644R	S/S	CVD	8 3/4



BC644R

1/2




BC644R

1/1


STRULLY Scissors

Catalog Number	Tip	Blade	Length (in.)
BC645R	PROBE	CVD	8 3/4




BC645R

1/2



BC645R

1/1




AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C




CARDIOVASCULAR SCISSORS



BC191R, BC643R, MG381R, US145

POTTS-SMITH Dissecting Scissors


Catalog Number	Tip	Blade	Length (in.)
MG381R	B/B	CVD	5 1/2
BC643R	B/B	CVD	7 1/2
BC191R	B/B	CVD	9
US145	B/B	CVD	11 1/4



MG381R 1/1



BC191R 1/1



BC643R 1/1



US145 1/1



BC196R

MORSE Scissors

Catalog Number	Tip	Blade	Length (in.)
BC196R	S/S	CVD	6 7/8




BC196R 1/1




BC638R-BC639R

HOHENFELLNER Scissors

Catalog Number	Tip	Blade	Length (in.)
BC638R	B/B	CVD	8 7/8
BC639R	B/B	CVD	10



BC638R 1/1





BC639R 1/1



CARDIOVASCULAR SCISSORS

JAMESON-METZENBAUM Scissors

Catalog Number	Tip	Blade	Length (in.)
MG383R	S/S	CVD	9





LILLEHEI-POTTS Tenotomy Scissors

Catalog Number	Tip	Blade	Length (in.)
MG382R	S/S	CVD	7



JAMESON Tenotomy Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
BC177R	B/B	CVD	6



CARDIOVASCULAR SCISSORS

POTTS Tenotomy Scissors

Catalog Number	Tip	Blade	Length (in.)
BC178R	B/B	CVD	6 7/8
BC179R	B/B	CVD	8



BC178R-BC179R



BC178R



BC179R

JAMESON SuperCut Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC911R	B/B	CVD	5 1/8
BC915R	B/B	CVD	6
BC919R	B/B	CVD	7



BC911R, BC915R, BC919R



BC911R



BC915R




BC919R



CARDIOVASCULAR DELICATE SCISSORS


METZENBAUM SuperCut Scissors

Catalog Number	Tip	Blade	Length (in.)
BC936R	B/B	STR	7
BC937R	B/B	CVD	7



METZENBAUM SuperCut Scissors

Catalog Number	Tip	Blade	Length (in.)
BC945R	B/B	CVD	9
BC947R	B/B	CVD	10 1/4
BC949R	B/B	CVD	11 1/4
BC951R	B/B	CVD	11 3/4
BC953R	B/B	CVD	14



CARDIOVASCULAR DELICATE SCISSORS

GORNEY Scissors


Catalog Number	Tip	Blade	Length (in.)
BC188R	B/B	STR	7 3/4
BC189R	B/B	SLIGHT-CVD	7 3/4



BC188R-BC189R



BC188R




BC189R

COOLEY Cardiovascular Scissors

Catalog Number	Tip	Blade	Length (in.)
BC597R	B/B	CVD	7 1/2
BC599R	B/B	CVD	9



BC597R-BC599R




BC597R-BC599R

POTTS-SMITH Lateral Angle Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC646R	S/S	55°	19	7



BC646R



BC646R



CARDIOVASCULAR DELICATE SCISSORS

POTTS-SMITH Reverse Scissors

Catalog Number	Tip	Blade Angle	Handle	Length (in.)
BC655R	S/S	130°	W/ STABILIZING BAR	7



POTTS-SMITH Vessel Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC652R	S/S	60°	10	6 7/8



POTTS-SMITH SuperCut Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC957R	B/B	25°	22	7 1/2
BC958R	B/B	45°	20	7 1/4
BC959R	B/B	60°	19	7





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier


A to Z Hospital Equipment and
Surgical Products B to C



CARDIOVASCULAR DELICATE SCISSORS

POTTS-DE MARTEL Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC647R	S/S	30°	21	7 1/4
BC648R	S/S	45°	21	7 1/4
BC649R	S/S	60°	21	7 1/4
BC650R	S/S	60°	21	8 3/4



BC647R-BC649R



BC647R

1/1



BC648R

1/1




BC649R

1



BC650R


1/2



BC650R

1/1

1/1



BC651R

1/2

MILLS Scissors

Catalog Number	Tip	Blade	Blade Length (mm)	Length (in.)
BC651R	BALL-POINT	60°	23	8 3/4



1/1



BC651R


1/1



CARDIOVASCULAR DELICATE SCISSORS



DIETHRICH Vessel Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC653R	S/S	90°	12	6 3/4
BC654R	S/S	120°	10	6 3/4



DIETHRICH-HEGEMANN Vessel Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC662R	S/S	25°	15	7
BC663R	S/S	45°	15	7
BC664R	S/S	60°	15	7



DIETHRICH Coronary Artery Scissors


Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC676R	S/S	25°	13	7 1/4
BC677R	S/S	45°	14	7 1/4



CARDIOVASCULAR DELICATE SCISSORS

COOLEY Arteriotomy Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC678R	B/B	30°	21.5	8 7/8
BC679R	B/B	30°	21.5	10 3/4



BC678R-BC679R



BC678R



BC679R

DEBAKEY Vessel Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC656R	S/S	45°	18	6 1/8
BC657R	S/S	60°	18	6 1/8



BC656R-BC657R



1/1



1/1

BC656R



1/1


BC657R

DEBAKEY O set Scissors, S-Shaped

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC669R	S/S	45°	24	6 1/8



BC669R



BC669R



CARDIOVASCULAR DELICATE SCISSORS


FAVOLORO Artery Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC670R	S/S	45°	20	5 1/2




POTTS Vascular Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
US200	S/S	45°	13	5 1/4




Valve Lea et Excision Scissors

Catalog Number	Tip	Blade Angle	Length (in.)
MB064R	B/B	45°	5 3/4




POTTS Arteriotomy Scissors

Catalog Number	Tip	Blade Angle	Length (in.)
MB065R	B/B	60°	5 3/8



COOLEY Neonatal Scissors


Catalog Number	Tip	Blade	Length (in.)
MB204R	S/S	CVD	5 3/4



CARDIOVASCULAR SCISSORS

DIETHRICH-POTTS Micro Scissors















Catalog Number	Tip	Blade	Blade Angle	Blade Length (mm)	Length (in.)
BC515R	S/S	X-DEL	25°	8	7
BC516R	S/S	X-DEL	45°	8	7
BC517R	S/S	X-DEL	60°	8	7
BC518R	S/S	X-DEL	90°	8	7
BC519R	S/S	X-DEL	125°	8	7
BC531R	S/S	DEL	25°	10	7
BC532R	S/S	DEL	45°	10	7
BC533R	S/S	DEL	60°	10	7
BC534R	S/S	DEL	90°	10	7
BC535R	S/S	DEL	125°	10	7
BC522R	S/S	FINE	25°	12	7
BC523R	S/S	FINE	45°	12	7
BC524R	S/S	FINE	60°	12	7
BC525R	S/S	FINE	90°	12	7
BC526R	S/S	FINE	125°	12	7



Blade Length: 8 MM

Blade Length: 10 MM

Blade Length: 12 MM



CARDIOVASCULAR SCISSORS

DIETHRICH SuperCut Scissors


Catalog Number	Tip	Blade	Length (in.)
MG346R	S/S	25°	7
MG347R	S/S	45°	7
MG348R	S/S	60°	7
MG349R	S/S	90°	7
MG350R	S/S	125°	7



CARDIOVASCULAR DUROTIP™ TC SCISSORS


DUROTIP TC Tungsten Carbide Scissors

- Excellent cutting accuracy
- Long working life
- Gold-plated rings




DUROTIP TC Wire Cutting Scissors

Catalog Number	Blade	For Wire Size	Length (in.)
DP552R	SERRATED	SOFT, TO 1.2 MM AND HARD, UP TO 0.8 MM	4 3/4



DUROTIP TC DEBAKEY Vascular Scissors

Catalog Number	Tip	Blade Angle	Blade Length (mm)	Length (in.)
BC245R	S/S	60°	21	8 3/4



DUROTIP TC DEMARTEL Vascular Scissors



Catalog Number	Tip	Blade	Length (in.)
MB212R	S/S	CVD	6 3/4



CARDIOVASCULAR DUROTIP™ TC SCISSORS







DUROTIP TC Baby METZENBAUM Dissecting Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC258R	B/B	STR	5 3/4
BC259R	B/B	CVD	5 3/4



DUROTIP TC METZENBAUM Dissecting Scissors



Catalog Number	Tip	Blade	Length (in.)
BC260R	B/B	STR	5 3/4
BC262R	B/B	STR	7
BC261R	B/B	CVD	5 3/4
BC263R	B/B	CVD	7
BC265R	B/B	CVD	7 7/8
BC287R	B/B	CVD	9



CARDIOVASCULAR DUROTIP™ TC SCISSORS




DUROTIP TC METZENBAUM Dissecting Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
MB205R	S/S	STR	9 1/4
MB207R	S/S	CVD	5 3/4
MD496	S/S	CVD	7
MB206R	S/S	CVD	9 1/4



DUROTIP TC MAYO Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC558R	B/B	STR	9
BC559R	B/B	CVD	9




CARDIOVASCULAR DUROTIP™ TC SCISSORS




Wavecut

■ Wavecut scissors have one serrated cutting edge to prevent tissue from slipping.

The W after the item number indicates the scissors are the wavecut style, e.g. BC211W.



BC273W



BC273W

1/1


1/2

DUROTIP TC TOENNIS-ADSON Wavecut Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC273W	B/B	CVD	6 7/8



BC259W



BC259W

1/1


1/2

DUROTIP TC Baby METZENBAUM Wavecut Dissecting Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC259W	B/B	CVD	5 3/4




BC261W-BC265W




BC261W

1/1



BC263W

1/1



BC265W


1/1

DUROTIP TC METZENBAUM Wavecut Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC261W	B/B	CVD	5 3/4
BC263W	B/B	CVD	7
BC265W	B/B	CVD	7 7/8



CARDIOVASCULAR DUROTIP™ TC SCISSORS




BC267W-BC281W

DUROTIP TC NELSON-METZENBAUM Wavecut Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC267W	B/B	CVD	9
BC279W	B/B	CVD	10 1/4
BC281W	B/B	CVD	11



BC267W



BC279W




BC281W




BC271W-BC277W

DUROTIP TC METZENBAUM Wavecut Dissecting Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
BC271W	B/B	CVD	7
BC275W	B/B	CVD	7 7/8
BC277W	B/B	CVD	9



BC271W



BC275W




BC277W



BC287W

DUROTIP TC METZENBAUM Wavecut Dissecting Scissors, Heavy

Catalog Number	Tip	Blade	Length (in.)
BC287W	B/B	CVD	9



BC287W



CARDIOVASCULAR DUROTIP™ TC SCISSORS


DUROTIP TC Ligature Scissors

■ Carbide inserts help ensure sharp cutting edge

(also in the distal area of the blades).






■ Wavecut helps prevent suture materials from slipping.

■ Identified by 1 gold-plated ring.




DUROTIP TC Ligature Wavecut Scissors


Catalog Number	Tip	Blade	Length (in.)
BC294W	B/B	CVD	5 3/4
BC295W	B/B	CVD	7
BC296W	B/B	CVD	9
BC297W	B/B	CVD	10 1/4
BC298W	B/B	CVD	13 3/4



CARDIOVASCULAR NOIR™ SUPREME CUT SCISSORS



Microtome Edge



Wavecut

Noir Supreme Cut Scissors

■ Noir

The combination of the renowned Aesculap Tungsten Carbide technology with the Aesculap Noir coating results in an extremely high surface hardness and therefore high wear resistance in daily hospital use.

■ Microtome edge

The extremely sharp cutting edge of the scissors blade facilitates gentle cutting of tissue with minimal use of force.

■ Wavecut helps prevent tissue from slipping out during the cutting process.

■ The WB after the item number indicates the Noir scissors are the wavecut style, e.g. BC211WB.

For more information about Noir Supreme Cut Scissors see our brochure number DOC637.



BC252B-BC253WB

Noir MAYO Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC252B	B/B	STR	6 3/4
BC253B	B/B	CVD	6 3/4
BC253WB	B/B	CVD	6 3/4



BC252B



BC253B, BC253WB



BC284B-BC284WB

Noir MAYO-LEXER Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC284B	B/B	CVD	6 1/2
BC284WB	B/B	CVD	6 1/2



BC284B-BC284WB



CARDIOVASCULAR NOIR™ SUPREME CUT SCISSORS

Noir LEXER Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC283B	B/B	DEL, CVD	6 1/2
BC283WB	B/B	DEL, CVD	6 1/2

Noir Iris Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC210B	S/S	STR	4 1/4
BC211B	S/S	CVD	4 1/4
BC211WB	S/S	CVD	4 1/4

Noir TOENNIS-ADSON Dissecting Scissors


Catalog Number	Tip	Blade	Length (in.)
BC273B	B/B	DEL, CVD	6 7/8
BC273WB	B/B	DEL, CVD	6 7/8



CARDIOVASCULAR NOIR™ SUPREME CUT SCISSORS

Noir METZENBAUM Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC261B	B/B	CVD	5 3/4
BC261WB	B/B	CVD	5 3/4
BC263B	B/B	CVD	7
BC263WB	B/B	CVD	7
BC271B	B/B	DEL, CVD	7
BC271WB	B/B	DEL, CVD	7
BC265B	B/B	CVD	7 7/8
BC265WB	B/B	CVD	7 7/8
BC275B	B/B	DEL, CVD	7 7/8
BC275WB	B/B	DEL, CVD	7 7/8
BC267B	B/B	CVD	9
BC267WB	B/B	CVD	9
BC277B	B/B	DEL, CVD	9
BC277WB	B/B	DEL, CVD	9



BC261B-BC277WB



BC261B, BC261WB



BC263B, BC263WB



BC271B, BC271WB



BC265B, BC265WB



BC275B, BC275WB



BC267B, BC267WB



BC277B, BC277WB



CARDIOVASCULAR NOIR™ SUPREME CUT SCISSORS

Noir NELSON-METZENBAUM Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC279B	B/B	DEL, CVD	10 1/4
BC279WB	B/B	DEL, CVD	10 1/4
BC281B	B/B	DEL, CVD	11
BC281WB	B/B	DEL, CVD	11



BC279B-BC281WB




BC279B, BC279WB



BC281B, BC281WB




CARDIOVASCULAR MICRO SCISSORS



MILLES Micro Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
FD020R	S/S	STR, SERR	FLAT, SERR	6 1/4
FD021R	S/S	CVD, SERR	FLAT, SERR	6 1/4




BIEMER Vein Scissors

Catalog Number	Blade	Handle	Length (in.)
FD028R	STR	FLAT, SERR	5 1/8

MILLES Micro Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
FD020R	S/S	STR, SERR	FLAT, SERR	6 1/4
FD021R	S/S	CVD, SERR	FLAT, SERR	6 1/4



YASARGIL Micro Vessel Scissors

Catalog Number	Blade	Handle	Length (in.)
FD022R	STR	FLAT, SERR	7 1/4
FD023R	CVD	FLAT, SERR	7 1/4





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



CARDIOVASCULAR MICRO SCISSORS

Micro Scissors, Flat Handle


Catalog Number	Tip	Blade	Blade Angle	Blade Length (mm)	Length (in.)
FM444R	S/S	X-DEL	STR	7	6 1/2
FM430R	S/S	X-DEL	CVD	7	6 1/2
FM431R	S/S	X-DEL	25°	7	6 1/2
FM432R	S/S	X-DEL	45°	7	6 1/2
FM433R	S/S	X-DEL	60°	7	6 1/2
FM434R	S/S	X-DEL	90°	7	6 1/2
FM435R	S/S	X-DEL	125°	7	6 1/2
FM442R	S/S	DEL	STR	10	6 1/2
FM436R	S/S	DEL	CVD	10	6 1/2
FM437R	S/S	DEL	25°	10	6 1/2
FM438R	S/S	DEL	45°	10	6 1/2
FM439R	S/S	DEL	60°	10	6 1/2
FM440R	S/S	DEL	90°	10	6 1/2
FM441R	S/S	DEL	125°	10	6 1/2
FM451R	S/1 BALL-PT	FINE	25°	15	6 1/2
FM452R	S/1 BALL-PT	FINE	45°	15	6 1/2
FM453R	S/1 BALL-PT	FINE	60°	15	6 1/2
FM454R	S/1 BALL-PT	FINE	90°	15	6 1/2
FM455R	S/1 BALL-PT	FINE	125°	15	6 1/2
FM470R	S/S	FINE	25°	15	6 1/2
FM471R	S/S	FINE	45°	15	6 1/2
FM472R	S/S	FINE	60°	15	6 1/2
FM473R	S/S	FINE	90°	15	6 1/2
FM474R	S/S	FINE	125°	15	6 1/2



Blade: STR

Blade: 25°

Blade: 45°



Blade: 60°

Blade: 90°

Blade: 125°






CARDIOVASCULAR MICRO SCISSORS

Micro Scissors, Round Handle



Catalog Number	Tip	Blade	Blade Angle	Blade Length (mm)	Length (in.)
FM503R	S/S	X-DEL	CVD	7	6 1/2
FM504R	S/S	X-DEL	25°	7	6 1/2
FM505R	S/S	X-DEL	45°	7	6 1/2
FM506R	S/S	X-DEL	60°	7	6 1/2
FM507R	S/S	X-DEL	90°	7	6 1/2
FM508R	S/S	X-DEL	125°	7	6 1/2
FM488R	S/S	DEL	STR	10	6 1/2
FM481R	S/S	DEL	CVD	10	6 1/2
FM482R	S/S	DEL	25°	10	6 1/2
FM483R	S/S	DEL	45°	10	6 1/2
FM484R	S/S	DEL	60°	10	6 1/2
FM485R	S/S	DEL	90°	10	6 1/2
FM486R	S/S	DEL	125°	10	6 1/2
FM489R	S/1 BALL-PT	DEL	25°	10	6 1/2
FM490R	S/1 BALL-PT	DEL	45°	10	6 1/2
FM491R	S/1 BALL-PT	DEL	60°	10	6 1/2
FM492R	S/1 BALL-PT	DEL	90°	10	6 1/2
FM493R	S/1 BALL-PT	DEL	125°	10	6 1/2
FM476R	S/S	FINE	25°	15	6 1/2
FM477R	S/S	FINE	45°	15	6 1/2
FM478R	S/S	FINE	60°	15	6 1/2
FM479R	S/S	FINE	90°	15	6 1/2
FM480R	S/S	FINE	125°	15	6 1/2



CARDIOVASCULAR SPECIAL CV SCISSORS


NAGEL Micro Scissors

Catalog Number	Tip	Blade	Length (in.)
FD056R	S/S	45°	7 7/8
FD057R	PROBE	45°	7 7/8




DEBAKEY Circum ex Micro Scissors

Catalog Number	Tip	Blade	Length (in.)
FD059R	S/S	125°	7 1/2





Micro Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
FD100R	S/S	STR	ROUND	5 3/4
FD101R	S/S	CVD	ROUND	5 3/4
FD103R	B/B	CVD	ROUND	5 3/4



Micro Vessel Scissors

Catalog Number	Tip	Blade	Length (in.)
FD105R	S/S	STR	7 1/2
FD106R	S/S	CVD	7 1/2
FD107R	S/S	DEL, STR	7 1/2
FD108R	S/S	DEL, CVD	7 1/2




CARDIOVASCULAR MICRO NEEDLE HOLDERS



BM563R-BM564R

1/2




BM563R



BM564R

CASTROVIEJO Micro Needle Holder, W/ Catch

Catalog Number	Tip	Jaw	Length (in.)
BM563R	STR	SMOOTH	5 3/4
BM564R	CVD	SMOOTH	5 3/4




FD045R-FD046R

1/2



FD045R



FD046R

JACOBSON Micro Needle Holder

Catalog Number	Handle	Jaw	Length (in.)
FD045R	W/O CATCH	SMOOTH	7 1/4
FD046R	W/ CATCH	SMOOTH	7 1/4



FD048R

1/2



FD048R

Micro Needle Holder, W/ Catch

Catalog Number	Jaw	Length (in.)
FD048R	SMOOTH	7 7/8



FD076R-FD079R

1/2



FD076R



FD079R

Micro Needle Holder, W/O Catch

Catalog Number	Tip	Jaw	Length (in.)
FD076R	STR	SMOOTH	6 1/4
FD079R	CVD	SMOOTH	6 1/4



FD080R-FD081R

1/2



FD080R



FD081R

Micro Needle Holder, W/O Catch

Catalog Number	Tip	Jaw	Length (in.)
FD080R	STR	SMOOTH	7 1/4
FD081R	CVD	SMOOTH	7 1/4



CARDIOVASCULAR MICRO NEEDLE HOLDERS

Micro Needle Holder, Bayonet-Shaped






Catalog Number	Tip	Length (in.)
FD120R	CVD-UP	7 3/4

BARRAQUER Micro Needle Holder, Round Knurled Handle

Catalog Number	Tip	Handle	Length (in.)
FD230R	STR	W/OCATCH	4 3/4
MA352R	STR	W/OCATCH	5 1/4
BM562R	CVD,DEL	W/O CATCH	4 3/4
MA360R	CVD	W/O CATCH	4 1/4
BM560R	CVD	W/O CATCH	4 1/2
FD231R	CVD	W/O CATCH	4 3/4
MA357R	CVD	W/O CATCH	5 1/4
FD232R	STR	W/ CATCH	4 3/4
MA351R	STR	W/ CATCH	5 1/4
FD233R	CVD	W/ CATCH	4 3/4
MA358R	CVD	W/ CATCH	5 1/4

REILL Micro Needle Holder, W/O Catch

Catalog Number	Tip	Jaw	Length (in.)
FD240R	STR	SMOOTH	6
FD241R	CVD	SMOOTH	6



CARDIOVASCULAR MICRO NEEDLE HOLDERS

Micro Needle Holder, Round Knurled Handle W/ Catch

Catalog Number	Tip	Jaw	Length (in.)
FD242R	STR	SMOOTH	6
FD243R	CVD	SMOOTH	6

Micro Needle Holder, Round Knurled Handle W/O Catch

Catalog Number	Tip	Jaw	Length (in.)
FD244R	STR	SMOOTH	7 1/4
FD245R	CVD	SMOOTH	7 1/4

Micro Needle Holder, Round Knurled Handle W/ Catch




Catalog Number	Tip	Jaw	Length (in.)
FD246R	STR	SMOOTH	7 1/4
FD247R	CVD	SMOOTH	7 1/4

Round Knurled Handle Micro Needle Holder, W/ Catch

Catalog Number	Tip	Jaw	Length (in.)
FM392R	STR	SMOOTH	9

MUELLER Micro Needle Holder, Plastic Handle W/ Catch

Catalog Number	Tip	Length (in.)
FM021R	CVD	5 3/4



CARDIOVASCULAR DUROGRIP™ TC MICRO NEEDLE HOLDERS

Aesculap features Needle Holders with Tungsten Carbide (TC) inserts.

A smooth jaw and three different serration patterns are available, depending on the size of the suture needle.



smooth



0.2

(16,000 teeth/in.2)



0.4

(4,000 teeth/in.2)



0.5

(2,500 teeth/in.2)

We recommend carbide inserts as follows:

Jaw Surface Smooth 0.2 mm
pitch of serration 0.4 mm pitch
of serration 0.5 mm pitch of
serration


Suture Material

9/0 - 11/0

6/0 - 10/0

4/0 - 6/0

up to 3/0



DUROGRIP TC CASTROVIEJO Micro Needle Holder, Flat Handle

Catalog Number	Jaw	Length (in.)
BM001R	SMOOTH	5 3/4
BM002R	SERR	5 3/4
BM004R	SERR	7
BM005R	SERR	8 1/2

BM001R

1/1

BM001R

BM002R

1/1

0.2

BM002R

BM004R

1/1

0.2

BM004R

BM005R

1/1

0.2

BM005R

US126

1/1

US126

DUROGRIP TC SPENCER Coronary Needle Holder

Catalog Number	Jaw	Length (in.)
US126	SMOOTH	7



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



CARDIOVASCULAR DUROGRIP™ TC MICRO NEEDLE HOLDERS



MB362R


1/2



MB362R


DUROGRIP TC JACOBSON Needle Holder, Flat Handle W/O Catch

Catalog Number	Jaw	Length (in.)
MB362R	SERR	8 3/4



MB214R

1/2



MB214R


DUROGRIP TC Micro Needle Holder, Flat Handle W/ Catch

Catalog Number	Jaw	Length (in.)
MB214R	SMOOTH	7 1/4



MB215R

1/2




MB215R

DUROGRIP TC MILLS Vascular Needle Holder, Round Handle W/ Catch


Catalog Number	Jaw	Length (in.)
MB215R	SERR	8



FM388R-FM390R



FM388R



FM390R


DUROGRIP TC Micro Needle Holder, Round Handle W/ Catch

Catalog Number	Jaw	Length (in.)
FM388R	SERR	7
FM390R	SERR	7 7/8



MG070R

1/2



MG070R

DUROGRIP TC GIANNINI Needle Holder

Catalog Number	Jaw	Length (in.)
MG070R	SERR	10 1/2



CARDIOVASCULAR DUROGRIP™ TC MICRO NEEDLE HOLDERS

■ Dia Dust diamond dust coating ensures an even harder and more wear resistant surface on the inside of the jaws. The structure of the Dia Dust coating enables a reliable grip for the needles, sutures, and tissue.

Dia Dust Micro Needle Holder, Flat Handle W/ Catch

Catalog Number	Tip	Suture Size	Length (in.)
FM540R	STR	7/0 AND SMALLER	6 3/4
FM541R	STR	5/0 AND SMALLER	5 3/4
FM542R	STR	5/0 AND SMALLER	7 1/4
FM543R	STR	5/0 AND SMALLER	8 1/4
FM544R	STR	5/0 AND SMALLER	9

Dia Dust Micro Needle Holder, Flat Handle W/O Catch


Catalog Number	Tip	Suture Size	Length (in.)
FM545R	STR	7/0 AND SMALLER	6 3/4
FM546R	STR	5/0 AND SMALLER	5 3/4
FM547R	STR	5/0 AND SMALLER	7 1/4
FM548R	STR	5/0 AND SMALLER	8 1/4
FM549R	STR	5/0 AND SMALLER	9

Dia Dust Micro Needle Holder, Flat Knurled Handle W/ Catch


Catalog Number	Tip	Suture Size	Length (in.)
FM550R	STR	5/0 AND SMALLER	7 1/4
FM551R	STR	5/0 AND SMALLER	8 1/4
FM552R	STR	5/0 AND SMALLER	9



FM540R-FM544R




FM540R



FM541R-FM544R



FM545R-FM549R



FM545R-FM549R



FM550R-FM552R



FM550R-FM552R

CARDIOVASCULAR DIA DUST™ MICRO NEEDLE HOLDERS



FM553R-FM555R



FM553R-FM555R

Dia Dust Micro Needle Holder, Flat Knurled Handle W/O Catch

Catalog Number	Tip	Suture Size	Length (in.)
FM553R	STR	5/0 AND SMALLER	7 1/4
FM554R	STR	5/0 AND SMALLER	8 1/4
FM555R	STR	5/0 AND SMALLER	9



FD258R-FD259R



FD258R



FD259R

Dia Dust Micro Needle Holder, Golf Ball Handle Design

Catalog Number	Tip	Suture Size	Length (in.)
FD258R	STR	7/0 AND SMALLER	7
FD259R	CVD	7/0 AND SMALLER	7



FM218R-FM230R



FM218R-FM230R



FM221R-FM223R



FM221R-FM223R

Dia Dust Micro Needle Holder, Round Handle W/ Catch

Catalog Number	Tip	Suture Size	Length (in.)
FM218R	STR	4/0 AND SMALLER	7 1/4
FM228R	STR	4/0 AND SMALLER	8 1/4
FM230R	STR	4/0 AND SMALLER	9
FM221R	CVD	4/0 AND SMALLER	8 1/4
FM223R	CVD	4/0 AND SMALLER	9



CARDIOVASCULAR DIA DUST™ MICRO NEEDLE HOLDERS

Dia Dust Micro Needle Holder, Round Handle W/ Catch

Catalog Number	Tip	Suture Size	Length (in.)
FM520R	STR	7/0 AND SMALLER	4 3/4
FM521R	STR	7/0 AND SMALLER	6 1/4
FM563R	STR	7/0 AND SMALLER	7 1/4
FM565R	STR	7/0 AND SMALLER	8 1/4
FM567R	STR	7/0 AND SMALLER	9
FM522R	STR	5/0 AND SMALLER	5 3/4
FM523R	STR	5/0 AND SMALLER	6 1/4
FM564R	STR	5/0 AND SMALLER	7 1/4
FM566R	STR	5/0 AND SMALLER	8 1/4
FM568R	STR	5/0 AND SMALLER	9

Dia Dust Micro Needle Holder, Round Handle W/O Catch

Catalog Number	Tip	Suture Size	Length (in.)
FM524R	STR	7/0 AND SMALLER	4 3/4
FM525R	STR	7/0 AND SMALLER	6 1/4
FM526R	STR	7/0 AND SMALLER	7 1/4
FM527R	STR	7/0 AND SMALLER	8 1/4
FM532R	STR	7/0 AND SMALLER	9
FM528R	STR	5/0 AND SMALLER	5 3/4
FM529R	STR	5/0 AND SMALLER	6 1/4
FM530R	STR	5/0 AND SMALLER	7 1/4
FM531R	STR	5/0 AND SMALLER	8 1/4
FM533R	STR	5/0 AND SMALLER	9

Dia Dust Micro Needle Holder, Round Handle W/ Catch

Catalog Number	Tip	Suture Size	Length (in.)
FM236R	STR	5/0 AND SMALLER	8 1/4



FM520R-FM568R



FM520R, FM521R, FM563R, FM565R, FM567R



FM522R, FM523R, FM564R, FM566R, FM568R



FM524R-FM533R



FM524R-FM527R, FM532R



FM528R-FM531R, FM533R



FM236R





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

CARDIOVASCULAR DIA DUST™ MICRO NEEDLE HOLDERS

Dia Dust Micro Needle Holder, Round Handle W/O Catch

Catalog Number	Tip	Suture Size	Length (in.)
FM276R	STR	5/0 AND SMALLER	8 1/4

Dia Dust Micro Needle Holder, Round Handle, Heavy, W/ Catch

Catalog Number	Tip	Description	Length (in.)
FM538R	STR	2/0 AND SMALLER	7 1/2
FM539R	STR	2/0 AND SMALLER	9
FM534R	STR	2/0 AND SMALLER	9 3/4

Dia Dust Micro Needle Holder, Round Handle, Heavy, W/O Catch

Catalog Number	Tip	Suture Size	Length (in.)
FM535R	STR	2/0 AND SMALLER	7 1/2
FM536R	STR	2/0 AND SMALLER	9
FM537R	STR	2/0 AND SMALLER	9 3/4

Dia Dust GLASSER Micro Needle Holder, Round Handle W/ Catch


Catalog Number	Tip	Suture Size	Length (in.)
FM560R	STR	5/0 AND SMALLER	5
FM561R	STR	5/0 AND SMALLER	7



FM276R



FM276R



FM534R-FM539R



FM534R-FM539R



FM535R-FM537R



FM535R-FM537R



FM560R



FM560R



FM561R



FM561R



CARDIOVASCULAR TITANIUM MICRO NEEDLE HOLDERS

Aesculap Titanium Micro Instruments

■ The extremely light Titanium material is an advantage for

lengthy cardiac and vascular operations. This helps to carry on operating over longer time periods with minimum fatigue.

■ The non-magnetic material makes working with particularly fine needles easier for you.

■ Non-reflecting surface finish.

■ Incomparable operating sensitivity: the Aesculap golf ball handle design.

■ Unmistakable identification: the Aesculap color code.

For more information about Titanium Micro Instruments see our brochure number DOC890.



All needle holder jaws are coated with Aesculap Dia Dust™ diamond dust coating. This coating ensures that the inner surfaces provide the best possible performance in hardness and wear resistance. The texture of the jaws greatly increases its surface area, allowing you to grasp fine needles, sutures and soft tissue safely and securely. The needle quickly aligns itself in an ideal position in the rough surface. This results in a perfect grip that quite simply prevents any needle slippage.

Titanium Micro Needle Holder, Golf Ball Handle Design W/ Catch

Catalog Number	Tip	Suture Size	Handle Color	Length (in.)
FM258T	STR	7/0 AND SMALLER	GREEN	7 1/4
FM259T	STR	7/0 AND SMALLER	GREEN	8 1/4
FM260T	STR	7/0 AND SMALLER	GREEN	9
FM238T	STR	5/0 AND SMALLER	BLUE	7 1/4
FM239T	STR	5/0 AND SMALLER	BLUE	8 1/4
FM240T	STR	5/0 AND SMALLER	BLUE	9

FM238T-FM260T

1/1

FM258T-FM260T

1/1

FM238T-FM240T

FM232T-FM254T

1/1

FM252T-FM254T


1/1

FM232T-FM234T


Titanium Micro Needle Holder, Golf Ball Handle Design W/O Catch

Catalog Number	Tip	Suture Size	Handle Color	Length (in.)
FM252T	STR	7/0 AND SMALLER	GREEN	7 1/4
FM253T	STR	7/0 AND SMALLER	GREEN	8 1/4
FM254T	STR	7/0 AND SMALLER	GREEN	9
FM232T	STR	5/0 AND SMALLER	BLUE	7 1/4
FM233T	STR	5/0 AND SMALLER	BLUE	8 1/4
FM234T	STR	5/0 AND SMALLER	BLUE	9

CARDIOVASCULAR TITANIUM, DIA DUST™, DUROGRIP™ TC NEEDLE HOLDERS




BM301R-BM328R



Dia Dust RYDER Needle Holder

Catalog Number	Tip	Suture Size	Length (in.)
BM301R	X-DEL	5/0 AND SMALLER	5 1/2
BM302R	X-DEL	5/0 AND SMALLER	7
BM325R	X-DEL	5/0 AND SMALLER	6 1/4
BM327R	X-DEL	5/0 AND SMALLER	8 3/4
BM328R	X-DEL	5/0 AND SMALLER	9 3/4




BM385T-BM387T



Titanium Dia Dust RYDER Needle Holder


Catalog Number	Tip	Suture Size	Length (in.)
BM385T	X-DEL	5/0 AND SMALLER	6 1/4
BM386T	X-DEL	5/0 AND SMALLER	7
BM387T	X-DEL	5/0 AND SMALLER	8 1/4



BM051R-BM052R

DUROGRIP TC DEBAKEY Needle Holder, O set

Catalog Number	Jaw	Length (in.)
BM051R	SERR	7 3/4
BM052R	SERR	9 5/8



MG071R



DUROGRIP TC Microvascular Needle Holder











Catalog Number	Jaw	Length (in.)
MG071R	SERR	6



CARDIOVASCULAR DUROGRIP™ TC NEEDLE HOLDERS

DUROGRIP TC DEBAKEY Needle Holder













Catalog Number	Tip	Jaw	Length (in.)
BM032R	DEL	SERR	6
BM033R	DEL	SERR	6 1/2
BM035R	DEL	SERR	7
BM034R	DEL	SERR	8 1/4
BM036R	DEL	SERR	9
BM037R	DEL	SERR	10 1/4
BM038R	DEL	SERR	12
BM063R	X-DEL	SERR	6 1/2
BM064R	X-DEL	SERR	8 1/4
BM040R	X-DEL	SERR	12



CARDIOVASCULAR DUROGRIP™ TC NEEDLE HOLDERS

DUROGRIP TC RYDER Needle Holder


Catalog Number	Tip	Jaw	Length (in.)
BM054R	X-DEL	SERR	5 1/4
BM055R	X-DEL	SERR	6 1/8
BM056R	X-DEL	SERR	6 7/8
BM057R	X-DEL	SERR	8 1/4
BM058R	X-DEL	SERR	9
BM044R	DEL	SERR	5 1/4
BM045R	DEL	SERR	6 7/8
BM046R	DEL	SERR	7 3/4
BM047R	DEL	SERR	10 1/2
BM074R	DEL	SMOOTH	6
BM075R	DEL	SMOOTH	7
BM076R	DEL	SMOOTH	8 7/8





CARDIOVASCULAR DUROGRIPTM TC NEEDLE HOLDERS

DUROGRIP TC STRATTE Needle Holder

Catalog Number	Tip	Jaw	Length (in.)
BM095R	CVD	SERR	8 1/4
BM096R	CVD	SERR	9




BM095R-BM096R




DUROGRIP TC COOLEY Intracardiac Needle Holder

Catalog Number	Jaw	Length (in.)
MB210R	SERR	9



MB210R




1/2

CARDIOVASCULAR FORCEPS


DEBAKEY Tissue Forceps

Catalog Number	Tip Size (mm)	Length (in.)
MB062R	2	14
MB915R	2	16
MB916R	3	16




WHEELER Vessel Forceps

Catalog Number	Tip	Length (in.)
MG440R	SERR	6




GERALD Dressing Forceps, Delicate

Catalog Number	Tip	Jaw	Length (in.)
BD228R	STR	SERR	6 7/8
BD230R	STR	SERR	9
BD229R	CVD	SERR	6 7/8




GERALD Tissue Forceps, Delicate

Catalog Number	Tip	Teeth	Length (in.)
BD662R	STR	1 x 2	6 7/8
BD664R	STR	1 x 2	9
BD663R	CVD	1 x 2	6 7/8



SELMAN Tissue Forceps

Catalog Number	Length (in.)
BD750R	6 1/8
MB795R	7 7/8
MG078R	9 3/4




CARDIOVASCULAR FORCEPS



BD801R

1/2




BD801R



MG428R

1/2



MG428R




BD802R

1/2




BD802R



BD802R



BD710R-BD712R



BD710R



BD711R




FB397R-FB399R



FB397R



FB398R



FB399R

TUTTLE-SINGLEY Lung Forceps

Catalog Number	Length (in.)
BD801R	9
MG428R	14

NELSON Lung Forceps

Catalog Number	Length (in.)
BD802R	9

RESANO Cardiac Valve Forceps

Catalog Number	Tip Size (mm)	Length (in.)
BD710	2.5	7
R	2.5	8
BD711	2.5	9

R
BD712
R

RESANO Cardiac Valve Forceps

Catalog Number	Tip Size (mm)	Length (in.)
FB397R	2	7
FB398R	2	9
FB399R	2	11 3/4


CARDIOVASCULAR FORCEPS

GERALD Micro Forceps, Flat Handle

Catalog Number	Tip	Length (in.)
BD190R	SERR	7
BD197R	SERR	9 3/4



BD190R-BD197R




GERALD Micro Forceps, Flat Handle

Catalog Number	Tip	Length (in.)
BD191R	DEBAKEY	7
BD195R	DEBAKEY	7 7/8




BD191R-BD195R



Micro Suture Forceps, W/ Platform

Catalog Number	Tip Size (mm)	Length (in.)
BD338R	0.3	4 3/4
BD264R	0.5	7 1/4



MUELLER Micro Forceps, Plastic Handle


Catalog Number	Teeth	Length (in.)
FM033R	-	6 1/4
FM034R	1 x 2	6 1/4



FM033R



FM034R



CARDIOVASCULAR FORCEPS



FM035R



FM035R

MUELLER Micro Ring Forceps, Plastic Handle

Catalog Number	Diam. of Ring (mm)	Length (in.)
FM035R	1.2	6 1/4




FM036R



FM036R

MUELLER Micro Forceps, Plastic Handle, W/ Platform

Catalog Number	Length (in.)
FM036R	6 1/4




FM063R




FM063R

Micro Forceps, Round Handle

Catalog Number	Tip	Counter-balance	Length (in.)
FM063R	CROSS-SERR	NO	7 7/8



FM590R-FM591R



FM590R



FM591R

Micro Forceps, Round Handle

Catalog Number	Tip	Counter-balance	Length (in.)
FM590R	DEBAKEY	NO	7 1/4
FM591R	DEBAKEY	NO	8 1/4
FM592R	DEBAKEY	YES	7 1/4
FM600R	DEBAKEY	YES	8 1/4



FM592R-FM600R



FM592R



FM600R



BD154R



BD154R

DUROGRIP TC (POTTS-SMITH) Thumb Forceps

Catalog Number	Tip	Length (in.)
BD154R	X-SERR	5 3/4


CARDIOVASCULAR DIA DUST™ FORCEPS



■ The Dia Dust diamond dust coating ensures an even harder and more wear resistant surface on the inside of the jaws. The structure of the Dia Dust coating enables a reliable grip for the needles, sutures, and tissue.



BD192R-BD193R



BD192R




BD193R

Dia Dust GERALD Micro Ring Forceps, Flat Handle

Catalog Number	Diam. of Ring (mm)	Tip	Length (in.)
BD192R	1.0	STR	7
BD193R	1.0	CVD	7



BD194R-BD196R



BD194R-BD196R

Dia Dust GERALD Micro Forceps, Flat Handle

Catalog Number	Tip	Tip Size (mm)	Length (in.)
BD194R	STR	1.0	7
BD196R	STR	1.0	7 7/8



BD200R-BD201R



BD200R-BD201R

Dia Dust THUMB Forceps, Flat Handle

Catalog Number	Tip	Tip Size (mm)	Length (in.)
BD200R	STR	1.0	7
BD201R	STR	1.0	8 1/4



FD253R



FD253R

Dia Dust Micro Forceps W/ Platform

Catalog Number	Tip	Counter-balance	Length (in.)
FD253R	STR	NO	7



CARDIOVASCULAR DIA DUST™ FORCEPS



FM053R-FM054R




FM053R-FM054R

Dia Dust Micro Ring Forceps, Round Handle

Catalog Number	Diam. of Ring (mm)	Tip	Counter-balance	Length (in.)
FM053R	2.5	STR	NO	7
FM054R	2.5	STR	NO	7 7/8



FM597R-FM598R



FM597R-FM598R



FM587R



FM587R



FM582R-FM584R



FM582R-FM584R

Dia Dust Micro Ring Forceps, Round Handle

Catalog Number	Diam. of Ring (mm)	Tip	Counter-balance	Length (in.)
FM597R	1.0	STR	NO	7 1/4
FM598R	1.0	STR	NO	8 1/4
FM587R	1.0	CVD	NO	7 1/4

Dia Dust Micro Ring Forceps, W/ Platform, Round Handle

Catalog Number	Diam. of Ring (mm)	Tip	Counter-balance	Length (in.)
FM582R	1.0	STR	NO	7 1/4
FM584R	1.0	STR	NO	8 1/4



FM586R-FM588R




FM586R-FM588R


Dia Dust Micro Forceps, W/ Platform, Round Handle

Catalog Number	Tip	Counter-balance	Length (in.)
FM586R	STR	NO	7 1/4
FM588R	STR	NO	8 1/4

CARDIOVASCULAR DIA DUST™ FORCEPS



FM571R-FM595R



Dia Dust Micro Ring Forceps, Round Handle

Catalog Number	Diam. of Ring (mm)	Tip	Counter-balance	Length (in.)
FM570R	1.0	STR	YES	6
FM572R	1.0	STR	YES	7 1/4
FM574R	1.0	STR	YES	8 1/4
FM576R	1.0	STR	YES	9
FM596R	1.0	STR	YES	9 3/4
FM579R	2.0	STR	YES	7 1/4
FM571R	1.0	CVD	YES	6
FM573R	1.0	CVD	YES	7 1/4
FM575R	1.0	CVD	YES	8 1/4
FM577R	1.0	CVD	YES	9
FM595R	1.0	CVD	YES	9 3/4



FM578R-FM599R



FM581R-FM593R



Dia Dust Micro Forceps, W/ Platform, Round Handle

Catalog Number	Tip	Counter-balance	Length (in.)
FM578R	STR	YES	6
FM580R	STR	YES	7 1/4
FM599R	STR	YES	8 1/4
FM594R	STR	YES	9 3/4
FM581R	CVD	YES	6
FM583R	CVD	YES	7 1/4
FM585R	CVD	YES	8 1/4
FM593R	CVD	YES	9 3/4



CARDIOVASCULAR TITANIUM FORCEPS

Aesculap Titanium Micro Instruments

■ The extremely light Titanium material is an advantage for

lengthy cardiac and vascular operations. This helps to carry on operating over longer time periods with minimum fatigue.

■ Non-reflecting surface finish – for optimal glare reduction. ■ The Aesculap golf ball handle design – gives incomparable operating sensitivity.

■ The Aesculap color code – provides unmistakable instrument identification.


For more information about Titanium Micro Instruments see our brochure number DOC890.



All titanium forceps jaws are coated with Aesculap Dia Dust™ diamond dust coating. This coating ensures that the inner surfaces provide the best possible performance in hardness and wear resistance. The texture of the jaws greatly increases its surface area, allowing you to grasp sutures and soft tissue safely and securely.



FM317T-FM319T



FM317T-FM319T




FM327T-FM329T



FM327T-FM329T



FM337T-FM339T



FM337T-FM339T

Titanium Micro Ring Forceps, Golf Ball Handle

Catalog Number	Diam. of Ring (mm)	Tip	Handle Color	Counter-balance	Length (in.)
FM317T	1.0	DIA DUST	BLUE	NO	7 1/4
FM318T	1.0	DIA DUST	BLUE	NO	8 1/4
FM319T	1.0	DIA DUST	BLUE	NO	9

Titanium Micro Forceps, W/ Platform, Golf Ball Handle

Catalog Number	Tip	Handle Color	Counter-balance	Length (in.)
FM327T	DIA DUST	GREEN	NO	7 1/4
FM328T	DIA DUST	GREEN	NO	8 1/4
FM329T	DIA DUST	GREEN	NO	9

Titanium Micro Forceps, Golf Ball Handle

Catalog Number	Tip	Handle Color	Counter-balance	Length (in.)
FM337T	DEBAKEY	PURPLE	NO	7 1/4
FM338T	DEBAKEY	PURPLE	NO	8 1/4
FM339T	DEBAKEY	PURPLE	NO	9

CARDIOVASCULAR TITANIUM FORCEPS

Titanium Standard Forceps

- Non-reflecting surface finish.
- Fatigue-free lightness.
- Blue handle for definite identification of an Aesculap Titanium Instrument.

For more information about Titanium forceps see our brochure number DOC890.



FB387T-FB389T



FB387T-FB389T



FB390T-FB392T



FB390T-FB392T



BD221T



BD221T



BD514T



BD514T



BD516T



BD516T

Titanium DEBAKEY Atraumatic Forceps

Catalog Number	Tip Size (mm)	Length (in.)
FB387T	1.5	5 3/4
FB388T	1.5	7 3/4
FB389T	1.5	9 1/2
FB390T	2.5	5 3/4
FB391T	2.5	7 3/4
FB392T	2.5	9 1/2

Titanium Micro ADSON Dissecting Forceps

Catalog Number	Tip	Length (in.)
BD221T	SERR	4 3/4

Titanium Micro ADSON Tissue Forceps, Delicate

Catalog Number	Tip	Teeth	Length (in.)
BD514T	SERR	1 x 2	4 3/4

Titanium Micro ADSON Tissue Forceps

Catalog Number	Tip	Teeth	Length (in.)
BD516T	SERR	1 x 2	4 3/4

CARDIOVASCULAR TITANIUM FORCEPS

Titanium GERALD Dissecting Forceps, Delicate



Catalog Number	Tip	Length (in.)
BD227T	SERR	6 3/4

Titanium GERALD Tissue Forceps

Catalog Number	Tip	Teeth	Length (in.)
BD657T	SERR	1 x 2	6 3/4

Titanium RESANO Grasping Forceps

Catalog Number	Tip	Length (in.)
BD715T	SERR	9




CARDIOVASCULAR TITANIUM CLAMPS

Titanium Ringhandled Instruments


- Non-reflecting surface finish.
- Fatigue-free lightness.
- Blue rings for definite identification of an Aesculap Titanium Instrument.

For more information about Titanium Ringhandled Instruments see our brochure number DOC890.



Titanium DEBAKEY Atraumatic Vascular Clamp

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
FC320T	STR	70	9



Titanium DEBAKEY Atraumatic Vascular Clamp

Catalog Number	Jaw	Jaw Length (mm)	Shaft	Length (in.)
FC321T	STR	60	ANG	8 1/2




CARDIOVASCULAR TITANIUM CLAMPS

Titanium DEBAKEY Atraumatic Vascular Clamp


Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
FC322T	CVD	70	8 7/8
FC323T	CVD	75	9 5/8



FC322T-FC323T




FC322T




FC323T

Titanium GLOVER Atraumatic Vascular Clamp

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
FC324T	CVD	60	8 3/4
FC325T	CVD	80	10 5/8



FC324T-FC325T



FC324T




FC325T

CARDIOVASCULAR TITANIUM CLAMPS

Titanium LAMBERT-KAY Atraumatic Vascular Clamp




Catalog Number	Jaw Angle	Jaw Length (mm)	Shaft	Length (in.)
FC326T	MULTI	45	ANG	7 3/4



CARDIOVASCULAR TITANIUM CLAMPS

Titanium DEBAKEY Atraumatic Vascular Clamp

Catalog Number	Jaw	Length (in.)
FC327T	SLIGHT-CVD	12 1/4
FC328T	SLIGHT-CVD	12 1/2





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register


CARDIOVASCULAR TITANIUM CLAMPS

Titanium DEBAKEY Atraumatic
Vascular Clamp

Catalog Number	Jaw	Shaft	Length (in.)
FC331T	CVD	ANG	7 3/4
FC330T	ANG	ANG	7




FC330T-FC331T



CARDIOVASCULAR TITANIUM CLAMPS

Titanium DEBAKEY Atraumatic Vascular Clamp


Catalog Number	Jaw	Length (in.)
FC329T	HEAVY-CVD	10 1/4



CARDIOVASCULAR TITANIUM CLAMPS

**Titanium DEBAKEY Atraumatic
Vascular Clamp**

Catalog Number	Jaw	Shaft	Length (in.)
FC332T	ANG	ANG	9 1/2
FC333T	ANG	ANG	8 1/2



CARDIOVASCULAR TITANIUM CLAMPS

Titanium DEBAKEY-SATINSKY Atraumatic Vascular Clamp


Catalog Number	Jaw Angle	Shaft	Length (in.)
FC336T	MULTI	ANG	7 7/8
FC335T	MULTI	ANG	9 5/8
FC334T	MULTI	ANG	9 3/4




FC334T-FC336T



FC336T



FC335T










FC334T

CARDIOVASCULAR TITANIUM CLAMPS

Titanium DEBAKEY Atraumatic Pediatric Vascular Clamp



Catalog Number	Jaw Angle	Length (in.)
FC337T	STR	4 1/4
FC338T	30°	4 1/4
FC339T	45°	4 1/4
FC340T	60°	4 1/4
FC341T	MULTI	3 3/4
FC342T	MULTI	4
FC343T	MULTI	4



CARDIOVASCULAR FORCEPS AND CLAMPS

Chest Tube Passer

Catalog Number	Description	Length (in.)
MG578R	W/O RATCHET	9 3/4
MG589R	W/O RATCHET	14
MG588R	W/ RATCHET	9 3/4






MG588R



CARDIOVASCULAR FORCEPS AND CLAMPS

Tubing Clamp

Catalog Number	Jaw	Length (in.)
AN909R	SMOOTH	6 1/2
MD452R	SERR	6
MD453R	SERR	7
MD454R	SERR	7 7/8



Tubing Clamp, W/ Guard


Catalog Number	Jaw	Length (in.)
MD652	SERR	6
MD653	SERR	7
MD654	SERR	7 7/8



CARDIOVASCULAR FORCEPS AND CLAMPS


SAROT Forceps, Delicate

Catalog Number	Teeth	Length (in.)
BH224R	-	9 1/2
BH234R	1 x 2	9 1/2




CRAFOORD Forceps, Delicate

Catalog Number	Jaw	Length (in.)
BH225R	CVD	9 5/8




HOLLE Forceps, Delicate

Catalog Number	Jaw	Length (in.)
BJ037R	CVD	5 3/4
BJ038R	CVD	7 1/4
BJ039R	CVD	9 1/4




CARDIOVASCULAR FORCEPS AND CLAMPS




COLLER Artery Forceps

Catalog Number	Jaw Angle	Length (in.)
MG421R	90°	9 1/2



SAROT Intrathoracic Artery Forceps

Catalog Number	Jaw	Length (in.)
BH226R	SLIGHT-CVD	9 3/4



JULIAN Thoracic Artery Forceps


Catalog Number	Jaw	Length (in.)
BH233R	CVD	9 1/4



CARDIOVASCULAR FORCEPS AND CLAMPS


LOVELACE Lung Forceps

Catalog Number	Length (in.)
MB056R	8




DUVAL Lung Forceps

Catalog Number	Length (in.)
FB930R	9



BABCOCK Thoracic Tissue Forceps


Catalog Number	Length (in.)
MB057R	9 1/2



CARDIOVASCULAR FORCEPS AND CLAMPS

Thoracoscopy Forceps


Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
MB458R	STR, LONG	158	12 1/2
MB296R	CVD, LONG	158	12 1/4



CARDIOVASCULAR FORCEPS AND CLAMPS


ROCHESTER-OCHSNER Artery Forceps

Catalog Number	Jaw	Teeth	Length (in.)
MB172R	STR	1 x 2	14



CARMALT Artery Forceps, Longitudinal Serrations


Catalog Number	Jaw	Length (in.)
MB182R	SLIGHT-CVD	12



CARDIOVASCULAR FORCEPS AND CLAMPS

COOLEY Bronchus Clamp

Catalog Number	Jaw	Length (in.)
MG449R	ANG	8 1/4




MG449R



MG449R

Bronchus Clamp, Longitudinal Serrations

Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB954R	90°	60	9
FB955R	90°	80	9



FB954R-FB955R



FB954R



FB955R



CARDIOVASCULAR FORCEPS AND CLAMPS

SAROT Bronchus Clamp

Catalog Number	Jaw	Length (in.)
MB610R	CVD-RIGHT	9
MB608R	CVD-LEFT	9



MB610R



MB608R





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

CARDIOVASCULAR FORCEPS AND CLAMPS

FOGARTY-Type Clamp

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
MB091R	STR	33	6 1/2
MB070R	STR	61	11 5/8
MB096R	STR	86	12 1/4
MB092R	ANG	33	4 5/8
MB094R	ANG	61	8 1/2
MB093R	ANG	61	8 3/4
MB095R	ANG	86	9 1/4



MB070R, MB091R-MB096R



MB092R



MB093R



MB094R



MB095R



MB091R



MB070R



MB096R

FOGARTY-Type Spring Clip Applicator

Catalog Number	Teeth	Length (in.)
MB097R	1 x 2	7 1/4




MB097R



MB097R



CARDIOVASCULAR VASCULAR CLAMPS





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier


A to Z Hospital Equipment and
Surgical Products B to C




CARDIOVASCULAR VASCULAR CLAMPS

SATINSKY Vena Cava Clamp

Catalog Number	Jaw	Jaw Size	Shaft	Length (in.)
FB602R	ANG, X-SERR	LARGE	STR	10 1/2
FB604R	ANG, X-SERR	SMALL	LIGHT-CVD	10 1/2
FB605R	ANG, X-SERR	LARGE	LIGHT-CVD	10 1/2



FB602R



FB602R



FB604R-FB605R



FB604R



FB605R



CARDIOVASCULAR VASCULAR CLAMPS

CRAFOORD-SELLORS Coarctation Clamp, Longitudinal Serrations

Catalog Number	Jaw	Length (in.)
FB633R	SLIGHT-CVD	9 1/2



FB633R

1/2




FB633R

1/1

CRAFOORD Auricle Clamp, Longitudinal Serrations

Catalog Number	Jaw	Length (in.)
FB634R	MEDIUM-CVD	9 1/2
FB635R	HEAVY-CVD	9 1/2



FB634R-FB635R

1/2



FB634R

1/1





FB635R

1/1

CARDIOVASCULAR VASCULAR CLAMPS

SEMB Suture Forceps

Catalog Number	Jaw	Length (in.)
FB961R	HEAVY-CVD	9 1/2



HALLMAN Tunnelling Forceps

Catalog Number	Jaw	Length (in.)
MG075R	SERR RING	13 1/2



CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY Mosquito Atraumatic Forceps

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
FB440R	STR	20	5
FB442R	STR	27	6 1/2
FB441R	CVD	20	5
FB443R	CVD	27	6 1/2



CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY-BAINBRIDGE Atraumatic Forceps


Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
FB446R	STR	40	6 1/8
FB448R	STR	65	7
FB447R	CVD	40	6 1/8
FB449R	CVD	65	7



CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY-HARKEN Atraumatic Vascular Clamps


Catalog Number	Tip Size (#)	Jaw	Jaw Length (mm)	Length (in.)
FB521R	1	CVD	86	10
FB522R	2	CVD	86	9 5/8
FB523R	3	CVD	86	9 5/8
FB524R	4	CVD	86	9 3/4



CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY-PEAN Atraumatic Forceps

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
FB526R	STR	60	7 7/8
FB532R	STR	89	10 1/2
FB534R	STR	99	11
FB527R	CVD	60	7 7/8
FB533R	CVD	89	10 1/4
FB535R	CVD	99	10 3/4



FB526R-FB535R



FB526R



FB532R



FB534R



FB527R



FB533R



FB535R




CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY Atraumatic Ring Handle Bulldog Clamp

Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB540R	STR	53	5 1/8
FB541R	CVD	53	5
FB542R	CVD*	53	5
FB543R	45°	53	5
FB544R	90°	53	4 1/8


*Instrument is S-Shaped




FB540R-FB541R, FB543R-FB544R




FB540R



FB543R



FB541R



FB544R



FB542R



1/1


53 mm

FB542R



CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY Atraumatic Multi-Purpose Clamp, Curved

Catalog Number	Jaw Length (mm)	Length (in.)
FB557R	31	6 1/2
FB558R	40	6 3/4




FB557R-FB558R




COOLEY Atraumatic Pediatric Vascular Clamp

Catalog Number	Jaw Angle	Jaw Length	Length (in.)
FB563R	MULTI	24(mm)	6 1/4



FB563R





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier


A to Z Hospital Equipment and
Surgical Products B to C



CARDIOVASCULAR VASCULAR CLAMPS

COOLEY Atraumatic Pediatric Vascular Clamp (continued)



Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB741R	25°	38	5 3/4
FB567R	30°	30	6 1/2
FB568R	60°	30	6 1/4
FB569R	90°	30	6
FB570R	MULTI	90	7 7/8



CARDIOVASCULAR VASCULAR CLAMPS






CASTANEDA Atraumatic Neonatal Clamp, Angled

Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB581R	30°	33	4 3/4
FB582R	45°	31	4 3/4



CASTANEDA Atraumatic Neonatal Clamp


Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB593R	CVD	15	6 1/8
FB590R	MULTI	18	6 1/8
FB591R	MULTI	23	6 1/8
FB592R	MULTI	25	6 1/8



CARDIOVASCULAR VASCULAR CLAMPS

DERRA Atraumatic Pediatric Forceps

Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB559R	MULTI	16	6 3/4
FB560R	MULTI	20	6 3/4
FB561R	MULTI	26	6 7/8




FB559R-FB561R



FB559R



FB560R





FB561R



CARDIOVASCULAR VASCULAR CLAMPS

LELAND-JONES Atraumatic Peripheral Clamp


Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB469R	STR	77	7 3/4
FB468R	15°	65	7 1/2
FB467R	30°	65	7 1/2



CARDIOVASCULAR VASCULAR CLAMPS


DEBAKEY Atraumatic Suture Forceps

Catalog Number	Jaw	Jaw Length (mm)	Shaft	Length (in.)
FB483R	CVD	94	CVD	9 1/2



DIETHRICH Atraumatic Aortic Clamp




Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB525R	MULTI	45	8 1/4



CARDIOVASCULAR VASCULAR CLAMPS



DEBAKEY Atraumatic Forceps

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
FB484R	C-CURVE	77	7 1/2
FB498R	C-CURVE	51	8 1/4
FB485R	D-CURVE	80	9
FB499R	D-CURVE	75	9



HENLY Atraumatic Clamp


Catalog Number	Jaw Length (mm)	Length (in.)
FB617R	75	7 7/8



CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY-RUMEL Atraumatic Forceps

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
FB486R	SIGHT-CVD	81	9 3/4
FB487R	CVD	81	9 5/8




CARDIOVASCULAR VASCULAR CLAMPS

LEES Atraumatic Bronchus Clamp

Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB573R	90°	42	8 3/4



FB573R




FB573R

LEES Atraumatic Bronchus Clamp

Catalog Number	Jaw Angle	Jaw Length (mm)	Shaft	Length (in.)
FB574R	ANG	40	CVD	9 5/8



FB574R




FB574R



CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY Atraumatic Multi-Purpose Clamp

Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB510R	60°	40	9 1/2
FB511R	60°	55	8 3/4
FB512R	60°	65	11
FB513R	60°	100	12




CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY Atraumatic Multi-Purpose Clamp

Catalog Number	Jaw Angle	Jaw Length (mm)	Length (in.)
FB515R	90°	45	8 7/8
FB516R	90°	55	8 1/4
FB517R	90°	72	10
FB518R	90°	108	10 3/4





FB515R-FB518R



CARDIOVASCULAR VASCULAR CLAMPS

**LEMOLE-STRONG Atraumatic Clamp,
Angled Shaft**

Catalog Number	Jaw Size (mm)	Length (in.)
FB618R	40 x 15	7 7/8
FB619R	45 x 20	7 7/8





AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C




Certified Register

CARDIOVASCULAR VASCULAR CLAMPS

STRONG Aorta Clamp, Straight Shaft


Catalog Number	Jaw	Length (in.)
MG074R	CVD	7 7/8



CARDIOVASCULAR VASCULAR CLAMPS

WYLIE Hypogastric Clamp

Catalog Number	Jaw	Jaw Length (mm)	Shaft	Length (in.)
MB241	CVD	70	ANG	8 1/2
MB240	CVD	75	ANG	9 3/4




CARDIOVASCULAR VASCULAR CLAMPS

WYLIE Atraumatic Clamp


Catalog Number	Jaw	Length (in.)
FB669R	HEAVY-CVD	8 3/4
FB670R	CVD-BACK	9



FB669R-FB670R



1/1
FB669R




1/1
FB670R





CARDIOVASCULAR VASCULAR CLAMPS

WYLIE Atraumatic Clamp

Catalog Number	Jaw	Length (in.)
FB672R	CVD-UP	6 1/4
FB673R	CVD-UP	7




FB672R-FB673R





CARDIOVASCULAR VASCULAR CLAMPS

WYLIE Atraumatic Clamp

Catalog Number	Jaw	Length (in.)
FB674R	CVD-DOWN	6
FB675R	CVD-DOWN	6 1/2



FB674R-FB675R



CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY Atraumatic Patent Ductus Clamp

Catalog Number	Jaw Length (mm)	Shaft	Length (in.)
FB450R	51	STR	8 1/4
FB451R	51	ANG	8 1/4




FB450R



FB450R



FB451R



FB451R


CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY Atraumatic Coarctation Clamp

Catalog Number	Jaw Length (mm)	Shaft	Length (in.)
FB453R	41	ANG	7 7/8
FB553R	50	ANG	6 1/2
FB455R	65	ANG	8 7/8
FB457R	80	ANG	9 1/2



FB453R-FB457R



15°




FB453R



FB553R



FB455R



FB457R



CARDIOVASCULAR VASCULAR CLAMPS

DEBAKEY Atraumatic Coarctation Clamp

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
FB452R	STR	41	8
FB552R	STR	50	6 3/4
FB454R	STR	65	8 7/8
FB456R	STR	80	9 1/2



FB452R-FB552R



FB452R



FB552R



FB454R



FB456R



MG423R



MG423R

NOON AV Fistula Clamp

Catalog Number	Jaw Angle	Length (in.)
MG423R	50°	4 3/4



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



CARDIOVASCULAR VASCULAR CLAMPS

GLOVER Atraumatic Clamp


Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
FB458R	CVD	65	8 1/2
FB460R	SPOONSHAPED	60	8 3/4
FB550R	SPOON SHAPED	70	9 5/8
FB459R	CVD	86	9 5/8
FB461R	FULL-CVD	47	7 3/4
FB462R	FULL-CVD	63	8 1/4



FB458R-FB550R




FB458R



FB459R




FB460R



FB461R



FB550R




FB462R



CARDIOVASCULAR VASCULAR CLAMPS

GLOVER Coarctation Clamp

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
MB893R	STR	60	9
MB895R	STR	70	10 1/4
MB894R	ANG	60	8 3/4
MB896R	ANG	70	10




MB893R, MB895R



MB893R



MB895R



MB894R, MB896R



MB894R









MB896R

CARDIOVASCULAR VASCULAR CLAMPS

GREGORY Atraumata Profunda Clamp






Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
MG014R	X-SMALL	29	5
MB750R	SMALL	57	5 3/4
MB751R	MEDIUM	62	7
MB752R	LARGE	75	7 3/4



CARDIOVASCULAR VASCULAR CLAMPS

DALE Femoral Popliteal Clamp

Catalog Number	Jaw	Length (in.)
FB546R	SMALL	6 7/8
FB547R	MEDIUM	7
FB548R	LARGE	7
FB549R	-	7





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



CARDIOVASCULAR VASCULAR CLAMPS

SELMAN Vessel Forceps

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
MG361R	STR, COOLEY SERR	25	5 1/4
MG362R	ANG, COOLEY SERR	25	5 1/4

