




Wholesale Supplier Electrosurgery_Instrument A to Z Hospital Equipment and Surgical Products **B to C**



Electrosurgery

BipoJET® Bipolar Scissors	1168
Non-Irrigating Bipolar Forceps	1169 - 1190
Irrigating Bipolar Forceps	1191 - 1194
Disposable Bipolar Forceps	1195 - 1197
Bipolar Connecting Cords	1198 - 1200

8275140877

9822688071

info@amramedtech.com



Sr No.29/1/2, Adarsh Nagar,
Uruli Devachi, Hadapsar, Pune
Branches :- All Maharashtra



amramedtech.com

ELECTROSURGERY BIPOJET[®] BIPOLAR SCISSORS

The versatility of BipoJET Dissecting Bipolar Scissors enables:

■ Parallel Cutting and
Coagulating

■ Precoagulation

■ Pinpoint Coagulation

■ Surface Coagulation



For more information about Bipolar Scissors see our brochure number DOC456.
See page 1199-1200 for appropriate connecting cords.



BC691R-BC693R



BC691R-BC693R



BC694R-BC696R



BC694R-BC696R



BC685R-BC690R



BC685R



BC686R

METZENBAUM Bipolar Dissecting Scissors

	Catalog Number	Length (in.)
BC691R	X-FINE	4 3/4
BC692R	X-FINE	5 1/2
BC693R	X-FINE	6 1/4

STEVENS Bipolar Dissecting Scissors

Catalog Number	Tip	Length (in.)
BC694R	TENOTOMY	4 3/4
BC695R	TENOTOMY	5 1/2
BC696R	TENOTOMY	6 1/4

METZENBAUM Bipolar Dissecting Scissors

	Catalog Number	Length (in.)
BC685R	STANDARD	7
BC687R	STANDARD	9
BC689R	STANDARD	11
BC686R	FINE	7
BC688R	FINE	9
BC690R	FINE	11

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

YASARGIL Bipolar Forceps Features:

- High precision, tissue preserving preparation instrument
- Stronger tension of the forceps allows for blunt spreading

of the tissue.

Tip Size Reference Chart

0.2 mm		0.8 mm	
0.25 mm		0.9 mm	
0.4 mm		1.0 mm	
0.5 mm		1.3 mm	
0.6 mm		1.5 mm	
0.7 mm		2.0 mm	

YASARGIL Style Bipolar Forceps, Stainless Steel


Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK746R	0.4	BAYONET	INS*	NON-IRRG*	3/4	5 1/2
GK747R	0.7	BAYONET	INS	NON-IRRG	3/4	5 1/2
GK748R	1.0	BAYONET	INS	NON-IRRG	3/4	5 1/2
GK749R	1.3	BAYONET	INS	NON-IRRG	3/4	5 1/2
GK760R	0.4	BAYONET	INS	NON-IRRG	1 3/8	6 1/8
GK761R	0.7	BAYONET	INS	NON-IRRG	1 3/8	6 1/8
GK762R	1.0	BAYONET	INS	NON-IRRG	1 3/8	6 1/8
GK763R	1.3	BAYONET	INS	NON-IRRG	1 3/8	6 1/8
GK764R	0.4	BAYONET	INS	NON-IRRG	2 1/4	6 7/8
GK765R	0.7	BAYONET	INS	NON-IRRG	2 1/4	6 7/8
GK766R	1.0	BAYONET	INS	NON-IRRG	2 1/4	6 7/8
GK767R	1.3	BAYONET	INS	NON-IRRG	2 1/4	6 7/8

* Key:

INS = Insulating

NON-IRRG = Non-Irrigating

IRRG = Irrigating



 1/1
GK746R, GK760R, GK764R

 1/1
GK747R, GK761R, GK765R


 1/1
GK748R, GK762R, GK766R

 1/1
GK749R, GK763R, GK767R

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

YASARGIL Style Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK768R	0.4	BAYONET	INS	NON-IRRIG	3	7 3/4
GK769R	0.7	BAYONET	INS	NON-IRRIG	3	7 3/4
GK770R	1.0	BAYONET	INS	NON-IRRIG	3	7 3/4
GK771R	1.3	BAYONET	INS	NON-IRRIG	3	7 3/4
GK772R	0.4	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK773R	0.7	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK774R	1.0	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK775R	1.3	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2



1/1
GK768R, GK772R

1/1
GK769R, GK773R


1/1
GK770R, GK774R

1/1
GK771R, GK775R

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

YASARGIL Style Bipolar Forceps, Stainless Steel








Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK788R	0.4	BAYONET	INS	NON-IRRIG	4 1/2	9 1/4
GK789R	0.7	BAYONET	INS	NON-IRRIG	4 1/2	9 1/4
GK790R	1.0	BAYONET	INS	NON-IRRIG	4 1/2	9 1/4
GK791R	1.3	BAYONET	INS	NON-IRRIG	4 1/2	9 1/4
GK793R	0.7	BAYONET	INS	NON-IRRIG	5 1/4	10



ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

YASARGIL Style Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK777R	45° UP	0.7	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK780R	90° UP	0.4	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK781R	90° UP	0.7	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK785R	90° DOWN	0.7	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2



YASARGIL Style Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Working Length (in.)	Length (in.)
GK831R	RING	2/3	4 1/2	9 1/4
GK832R	RING	3/4	4 1/2	9 1/4
GK833R	RING	4/5	4 1/2	9 1/4

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK660R	0.6	BAYONET	INS	NON-IRRIG	1 1/2	6 1/4
GK670R	0.6	BAYONET	INS	NON-IRRIG	3 1/8	7 7/8
GK671R	0.9	BAYONET	INS	NON-IRRIG	3 1/8	7 7/8
GK675R	0.6	BAYONET	INS	NON-IRRIG	4	8 3/4
GK677R	0.6	BAYONET	INS	NON-IRRIG	4 3/4	9 1/2
GK679R	0.9	BAYONET	INS	NON-IRRIG	4 3/4	9 1/2

40 mm, 1 1/2"

GK660R

1/2

80 mm, 3 1/8"

GK670R-GK671R

1/2

100 mm, 4"

GK675R

1/2

120 mm, 4 3/4"

GK677R-GK679R

1/2

1/1

GK660R, GK670R, GK675R, GK677R

1/1

GK671R, GK679R

ELECTRO
SURGERY

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

Bayonet Forceps, Rose Gold, Insulated

■ Rose Gold plated tips reduce tissue adhesion



Non-Stick Rose Gold forceps

The Rose Gold forceps enable heat to be quickly discharged away from the tip.

Standard forceps

The heat cannot be discharged fast enough. This leads to overheating and sticking of tissue.

Rose Gold Bipolar Forceps, Bayonet-Shaped

Catalog Number	Tip Size (mm)	Style	Insulation	Irrigation	Working Length (in.)	Length (in.)
US403	0.5	HARDY	INS	NON-IRRG	2	6
US404	0.7	HARDY	INS	NON-IRRG	2	6
US405	1.0	HARDY	INS	NON-IRRG	2	6
US406	1.5	HARDY	INS	NON-IRRG	2	6
US407	2.0	HARDY	INS	NON-IRRG	2	6
US605	0.5	SCOVILLE	INS	NON-IRRG	2 3/4	7 3/4
US606	0.7	SCOVILLE	INS	NON-IRRG	2 3/4	7 3/4
US607	1.0	SCOVILLE	INS	NON-IRRG	2 3/4	7 3/4
US608	1.5	SCOVILLE	INS	NON-IRRG	2 3/4	7 3/4
US609	2.0	SCOVILLE	INS	NON-IRRG	2 3/4	7 3/4
US600	0.5	CUSHING	INS	NON-IRRG	3 1/2	7 1/2
US601	0.7	CUSHING	INS	NON-IRRG	3 1/2	7 1/2
US602	1.0	CUSHING	INS	NON-IRRG	3 1/2	7 1/2
US603	1.5	CUSHING	INS	NON-IRRG	3 1/2	7 1/2
US604	2.0	CUSHING	INS	NON-IRRG	3 1/2	7 1/2
US610	0.5	HARDY	INS	NON-IRRG	4 1/2	8 1/2
US611	0.7	HARDY	INS	NON-IRRG	4 1/2	8 1/2
US612	1.0	HARDY	INS	NON-IRRG	4 1/2	8 1/2
US613	1.5	HARDY	INS	NON-IRRG	4 1/2	8 1/2
US614	2.0	HARDY	INS	NON-IRRG	4 1/2	8 1/2
US615	0.5	HARDY	INS	NON-IRRG	5	9
US616	0.7	HARDY	INS	NON-IRRG	5	9
US617	1.0	HARDY	INS	NON-IRRG	5	9
US618	1.5	HARDY	INS	NON-IRRG	5	9
US619	2.0	HARDY	INS	NON-IRRG	5	9
US530	0.5	HARDY	INS	NON-IRRG	6	10
US620	0.7	HARDY	INS	NON-IRRG	6	10
US621	1.0	HARDY	INS	NON-IRRG	6	10
US173	1.5	HARDY	INS	NON-IRRG	6	10
US531	2.0	HARDY	INS	NON-IRRG	6	10

US173-US621

US403, US530, US600, US610, US615

US404, US601, US611, US616, US621

US405, US602, US612, US617

US173, US406, US603, US613, US618

US407, US531, US604, US614, US619

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

Rose Gold Microform Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US526	0.5	BAYONET	INS	NON-IRRG	1 3/8	6
US527	0.7	BAYONET	INS	NON-IRRG	1 3/8	6
US528	1.0	BAYONET	INS	NON-IRRG	1 3/8	6
US529	1.3	BAYONET	INS	NON-IRRG	1 3/8	6
US570	0.5	BAYONET	INS	NON-IRRG	2 3/8	7
US571	0.7	BAYONET	INS	NON-IRRG	2 3/8	7
US572	1.0	BAYONET	INS	NON-IRRG	2 3/8	7
US573	1.3	BAYONET	INS	NON-IRRG	2 3/8	7
US574	0.5	BAYONET	INS	NON-IRRG	3	7 3/4
US575	0.7	BAYONET	INS	NON-IRRG	3	7 3/4
US576	1.0	BAYONET	INS	NON-IRRG	3	7 3/4
US577	1.3	BAYONET	INS	NON-IRRG	3	7 3/4
US578	0.5	BAYONET	INS	NON-IRRG	3 3/4	8 1/2
US579	0.7	BAYONET	INS	NON-IRRG	3 3/4	8 1/2
US580	1.0	BAYONET	INS	NON-IRRG	3 3/4	8 1/2
US581	1.3	BAYONET	INS	NON-IRRG	3 3/4	8 1/2
US582	0.5	BAYONET	INS	NON-IRRG	4 1/2	9 1/4
US583	0.7	BAYONET	INS	NON-IRRG	4 1/2	9 1/4
US584	1.0	BAYONET	INS	NON-IRRG	4 1/2	9 1/4
US585	1.3	BAYONET	INS	NON-IRRG	4 1/2	9 1/4
US586	0.5	BAYONET	INS	NON-IRRG	5 1/2	10 1/4
US587	0.7	BAYONET	INS	NON-IRRG	5 1/2	10 1/4
US588	1.0	BAYONET	INS	NON-IRRG	5 1/2	10 1/4
US589	1.3	BAYONET	INS	NON-IRRG	5 1/2	10 1/4



ELECTRO
SURGERY

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

Rose Gold Reverse Action Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US547	0.7	BAYONET	INS	NON-IRRIG	6	10
US548	1.0	BAYONET	INS	NON-IRRIG	6	10



Rose Gold Minimally Invasive (MIS) Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US541	STR	0.7	BAYONET	INS	NON-IRRIG	5 1/4	9 1/2
US542	STR	1.2	BAYONET	INS	NON-IRRIG	5 1/4	9 1/2
US543	STR	0.7	BAYONET	INS	NON-IRRIG	6 1/4	10 1/2
US544	STR	1.2	BAYONET	INS	NON-IRRIG	6 1/4	10 1/2
US545	45° UP	0.7	BAYONET	INS	NON-IRRIG	5 1/4	9 1/2
US546	45° UP	0.7	BAYONET	INS	NON-IRRIG	6 1/4	10 1/2



Rose Gold SEMKIN Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US413	1.5	STR	INS	NON-IRRIG	6



ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

SCOVILLE Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US534	1.0	BAYONET	INS	NON-IRRIG	2 3/4	7 3/4
US301	1.5	BAYONET	INS	NON-IRRIG	2 3/4	7 3/4
US300	1.5	BAYONET	NON-INS	NON-IRRIG	2 3/4	7 3/4




GERALD Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US342	1.0	BAYONET	INS	NON-IRRIG	3 1/2	7 1/2
US307	1.0	BAYONET	NON-INS	NON-IRRIG	3 1/2	7 1/2



CUSHING Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US308	0.7	BAYONET	INS	NON-IRRIG	3 1/2	7 1/2
US182	1.0	BAYONET	INS	NON-IRRIG	3 1/2	7 1/2
US183	1.5	BAYONET	INS	NON-IRRIG	3 1/2	7 1/2
US305	2.0	BAYONET	INS	NON-IRRIG	3 1/2	7 1/2



ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

HARDY Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US310	0.5	BAYONET	INS	NON-IRRIG	4 1/2	8 1/2
US536	0.7	BAYONET	INS	NON-IRRIG	4 1/2	8 1/2
US184	1.0	BAYONET	INS	NON-IRRIG	4 1/2	8 1/2
US185	1.5	BAYONET	INS	NON-IRRIG	4 1/2	8 1/2
US186	2.0	BAYONET	INS	NON-IRRIG	4 1/2	8 1/2
US187	0.5	BAYONET	INS	NON-IRRIG	6	10
US188	1.0	BAYONET	INS	NON-IRRIG	6	10
US189	1.5	BAYONET	INS	NON-IRRIG	6	10

US184-US189, US310, US536

US187, US310

US536

US184, US188

US185, US189

US186

HARDY Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US309	0.5	BAYONET	NON-INS	NON-IRRIG	4 1/2	8 1/2

US309

US309

ELECTRO
SURGERY

MILLIMETERS



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

BARRÉ BipoJET® Nerve/Tissue Sparing Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
GK685R	3.0	BAYONET	INS	NON-IRRIG	10 1/4



Strong Tension Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US314	0.5	BAYONET	NON-INS	NON-IRRIG	2 3/8	6 1/4
US315	0.7	BAYONET	NON-INS	NON-IRRIG	2 3/8	6 1/4
US316	1.0	BAYONET	NON-INS	NON-IRRIG	2 3/8	6 1/4
US317	1.3	BAYONET	NON-INS	NON-IRRIG	2 3/8	6 1/4
US318	0.5	BAYONET	NON-INS	NON-IRRIG	3 1/4	7 1/8
US319	0.7	BAYONET	NON-INS	NON-IRRIG	3 1/4	7 1/8
US320	1.0	BAYONET	NON-INS	NON-IRRIG	3 1/4	7 1/8
US321	1.3	BAYONET	NON-INS	NON-IRRIG	3 1/4	7 1/8
US322	0.5	BAYONET	NON-INS	NON-IRRIG	4	7 7/8
US323	0.7	BAYONET	NON-INS	NON-IRRIG	4	7 7/8
US324	1.0	BAYONET	NON-INS	NON-IRRIG	4	7 7/8
US325	1.3	BAYONET	NON-INS	NON-IRRIG	4	7 7/8



ELECTRO
SURGERY





ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

CASPAR Style Bipolar Forceps:

- Developed for precision neurosurgery with a special focus on spine surgery.
- Parallel guiding handle for exact tip alignment.

CASPAR Style Bipolar Forceps, Stainless Steel





Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK899R	0.2	BAYONET	INS	NON-IRRG	2 3/4	6 1/2
GK900R	0.5	BAYONET	INS	NON-IRRG	2 3/4	6 1/2
GK910R	1.0	BAYONET	INS	NON-IRRG	2 3/4	6 1/2
GK920R	2.0	BAYONET	INS	NON-IRRG	2 3/4	6 1/2
GK929R	0.2	BAYONET	INS	NON-IRRG	4	7 3/4
GK930R	0.5	BAYONET	INS	NON-IRRG	4	7 3/4
GK940R	1.0	BAYONET	INS	NON-IRRG	4	7 3/4
GK950R	2.0	BAYONET	INS	NON-IRRG	4	7 3/4
GK959R	0.2	BAYONET	INS	NON-IRRG	5	8 3/4
GK960R	0.5	BAYONET	INS	NON-IRRG	5	8 3/4
GK970R	1.0	BAYONET	INS	NON-IRRG	5	8 3/4
GK980R	2.0	BAYONET	INS	NON-IRRG	5	8 3/4



ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

CASPAR Style Bipolar Forceps, Stainless Steel


Catalog Number	Tip	Shape	Tip Size (mm)	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK976R	90° UP	BAYONET	1.0	INS	NON-IRRIG	4 3/4	8 1/2
GK978R	90° DOWN	BAYONET	1.0	INS	NON-IRRIG	4 3/4	8 1/2
GK972R	45° UP	BAYONET	1.0	INS	NON-IRRIG	5	8 3/4
GK974R	45° DOWN	BAYONET	1.0	INS	NON-IRRIG	5	8 3/4



ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

SINTRAM Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK644R	STR	0.7 x 6	BAYONET	INS	NON-IRRIG	2 3/8	6 1/4
GK645R	STR	1.0 x 6	BAYONET	INS	NON-IRRIG	3 1/4	7 1/4
GK646R	STR	1.0 x 6	BAYONET	INS	NON-IRRIG	4	7 7/8
GK614R	STR	0.7 x 23	BAYONET	INS	NON-IRRIG	2 3/8	6 1/4
GK615R	STR	1.0 x 23	BAYONET	INS	NON-IRRIG	3 1/4	7 1/4
GK616R	STR	1.0 x 23	BAYONET	INS	NON-IRRIG	4	7 7/8
GK604R	STR	0.9 x 23	STR	INS	NON-IRRIG	-	6 1/4
GK605R	STR	0.9 x 23	STR	INS	NON-IRRIG	-	7 1/4
GK606R	STR	0.9 x 23	STR	INS	NON-IRRIG	-	7 7/8
GK624R	ANG 35°	0.9 x 23	STR	INS	NON-IRRIG	-	6 1/4
GK625R	ANG 35°	0.9 x 23	STR	INS	NON-IRRIG	-	7 7/8



ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

MALIS Style Titanium Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US331	STR	0.5	BAYONET	NON-INS	NON-IRRIG	3 3/4	8 3/4
US332	STR	1.0	BAYONET	NON-INS	NON-IRRIG	3 3/4	8 3/4
US330	STR	1.5	BAYONET	NON-INS	NON-IRRIG	3 3/4	8 3/4
US333	ANG-UP	0.5	BAYONET	NON-INS	NON-IRRIG	3 3/4	8 3/4
US175	ANG-DOWN	0.5	BAYONET	NON-INS	NON-IRRIG	3 3/4	8 3/4



US330-US332



US330-US332



US175, US333



US175, US333

MALIS Style Rose Gold Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US161	STR	0.5	BAYONET	INS	NON-IRRIG	3 3/4	8 3/4
US162	STR	1.0	BAYONET	INS	NON-IRRIG	3 3/4	8 3/4
US163	STR	1.5	BAYONET	INS	NON-IRRIG	3 3/4	8 3/4
US164	ANG-UP	0.5	BAYONET	INS	NON-IRRIG	3 3/4	8 3/4
US165	ANG-DOWN	0.5	BAYONET	INS	NON-IRRIG	3 3/4	8 3/4



US161-US165



US161



US162



US163



US164



US165

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

Titanium Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US326	STR	0.5	STR	NON-INS	NON-IRRIG	6
US328	STR	0.5	STR	NON-INS	NON-IRRIG	7 1/8
US327	CVD	0.5	STR	NON-INS	NON-IRRIG	6
US329	CVD	0.5	STR	NON-INS	NON-IRRIG	7 1/8



US326, US328




US327, US329

Micro Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Irrigation	Insulation	Length (in.)
US345	STR	0.7	STR	NON-INS	NON-IRRIG	6
US347	STR	0.7	STR	NON-INS	NON-IRRIG	7 1/8
US346	CVD	0.7	STR	NON-INS	NON-IRRIG	6
US348	CVD	0.7	STR	NON-INS	NON-IRRIG	7 1/8



US345, US347



US346, US348



GK700R-GK703R



GK700R, GK701R, GK703R



GK711R-GK713R





GK711R, GK713R

Classic Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
GK700R	STR	0.6	STR	INS	NON-IRRIG	6 1/4
GK701R	STR	0.9	STR	INS	NON-IRRIG	6 1/4
GK703R	STR	0.9	STR	INS	NON-IRRIG	7 7/8
GK711R	ANG-DOWN	0.9	STR	INS	NON-IRRIG	6 1/4
GK713R	ANG-DOWN	0.9	STR	INS	NON-IRRIG	7 7/8

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS



Classic Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
GK690R	STR	0.4	STR	INS	NON-IRRIG	4 3/4
GK720R	STR	0.6	STR	INS	NON-IRRIG	6 1/4
GK722R	STR	0.9	STR	INS	NON-IRRIG	6 1/4
GK730R	STR	0.9	STR	INS	NON-IRRIG	7 7/8
GK723R	ANG-DOWN	0.9	STR	INS	NON-IRRIG	6 1/4

Jewelers Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US337	STR	0.4	STR	INS	NON-IRRIG	4
US335	STR	0.4	STR	NON-INS	NON-IRRIG	4
US159	CVD	0.4	STR	INS	NON-IRRIG	4
US336	CVD	0.4	STR	NON-INS	NON-IRRIG	4

TENZEL Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US178	CVD	0.4	STR	INS	NON-IRRIG	4
US338	CVD	0.4	STR	NON-INS	NON-IRRIG	4

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS



US339-US340

ADSON Bipolar Forceps, Stainless Steel


Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US340	STR	1.0	STR	INS	NON-IRRIG	5
US339	STR	1.0	STR	NON-INS	NON-IRRIG	5



US179-US344

SEMKIN Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US179	0.7	STR	INS	NON-IRRIG	5 1/2
US311	0.7	STR	NON-INS	NON-IRRIG	5 1/2
US344	2.0	STR	INS	NON-IRRIG	6
US196	2.0	STR	NON-INS	NON-IRRIG	6



US197-US304

CUSHING Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US197	0.7	STR	INS	NON-IRRIG	6
US302	0.7	STR	NON-INS	NON-IRRIG	6
US304	1.5	STR	INS	NON-IRRIG	7
US303	1.5	STR	NON-INS	NON-IRRIG	7

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

GERALD Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Shape	Insulation	Irrigation	Length (in.)
US341	STR	STR	INS	NON-IRRIG	7
US306	STR	STR	NON-INS	NON-IRRIG	7



US306-US341



US171



US171



US172



US172



US176-US177



US176-US177

Rose Gold Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Style	Insulation	Irrigation	Length (in.)
US171	1.0	STR	ADSON	INS	NON-IRRIG	5
US172	0.4	STR	JEWELERS	INS	NON-IRRIG	4

MALIS Style Titanium Bipolar Forceps


Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US176	0.25	STR	NON-INS	NON-IRRIG	6
US177	0.5	STR	NON-INS	NON-IRRIG	6



GK595R-GK596R



GK595R



GK596R



GK597R



GK597R

Micro and Laboratory Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
GK596R	STR	0.3	STR	INS	NON-IRRIG	5 3/4
GK595R	STR	0.6	STR	INS	NON-IRRIG	4 3/4
GK597R	CVD	0.6	STR	INS	NON-IRRIG	4 1/4

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

BiProtect®

Ceramic Insulated Tip, Insulated

■ Ceramic tip insulation helps ensure that heat concentrates on the inner and distal part of the forceps

- Less collateral damage and reduced thermal spread
- Partially insulated tip allows for precise and safe surface coagulation

YASARGIL BiProtect Ceramic Insulated Tip Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GL470R	STR	0.4	BAYONET	INS	NON-IRRIG	3 1/2	7 3/4
GL471R	STR	0.8	BAYONET	INS	NON-IRRIG	3 1/2	7 3/4
GL480R	STR	0.4	BAYONET	INS	NON-IRRIG	4 1/2	8 1/2
GL481R	STR	0.8	BAYONET	INS	NON-IRRIG	4 1/2	8 1/2

GL470R-GL481R

GL470R, GL480R 1/1

GL471R, GL481R 1/1

Transnasal Bipolar Forceps, Stainless Steel


Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK719R	1.0	ANG	INS	NON-IRRIG	5 1/8	8 1/2



ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

MIN Instruments, Pivot-Point Bipolar Forceps, Transnasal

Catalog Number	Tip	Tip Size (mm)	Shape	Working Length (in.)	Length (in.)
GK822R	STR	0.7	BAYONET	3 3/4	8 1/2
GK824R	STR	1.0	BAYONET	3 3/4	8 1/2
GK826R	STR	0.7	BAYONET	5 1/4	10
GK828R	STR	1.0	BAYONET	5 1/4	10
GK823R	CVD	0.7	BAYONET	3 3/4	8 1/2
GK825R	CVD	1.0	BAYONET	3 3/4	8 1/2
GK827R	CVD	0.7	BAYONET	5 1/4	10
GK829R	CVD	1.0	BAYONET	5 1/4	10



Special pin between the branches opens the tip of the forceps by additional compression of the handle, allowing secure coagulation in narrow and deep seated surgical fields.

95 mm, 3 7/8"

GK822R-GK824R

1/2

135 mm, 5 3/8"

GK826R-GK828R

1/2

95 mm, 3 7/8"

GK823R-GK825R

1/2

135 mm, 5 3/8"

GK827R-GK829R

1/2

1/1

GK822R, GK823R, GK826R, GK827R

1/1


GK824R, GK825R, GK828R, GK829R

ELECTRO
SURGERY

ELECTROSURGERY NON-IRRIGATING BIPOLAR FORCEPS

LANDOLT Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Insulation	Irrigation	Description	Working Length (in.)
GK560R	ANG 90°	INS	NON-IRRIG	COAGULATION FORCEPS F/ HYPOPHYSECTOMY	6
GK561R	90°	-	-	SPARE INNER PART F/ GK560R	
GK580R	ANG 120°	INS	NON-IRRIG	COAGULATION FORCEPS F/ HYPOPHYSECTOMY	6
GK581R	120°	-	-	SPARE INNER PART F/ GK580R	




ELECTROSURGERY IRRIGATING BIPOLAR FORCEPS

■ Irrigating coagulation forceps facilitate gentle coagulation without tissue adhesion.

Rose Gold Irrigating Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US622	0.5	BAYONET	INS	IRRG	3 1/2	7 1/2
US623	0.7	BAYONET	INS	IRRG	3 1/2	7 1/2
US624	1.0	BAYONET	INS	IRRG	3 1/2	7 1/2
US625	1.5	BAYONET	INS	IRRG	3 1/2	7 1/2
US626	2.0	BAYONET	INS	IRRG	3 1/2	7 1/2
US627	0.5	BAYONET	INS	IRRG	3	7 3/4
US628	0.7	BAYONET	INS	IRRG	3	7 3/4
US629	1.0	BAYONET	INS	IRRG	3	7 3/4
US630	1.5	BAYONET	INS	IRRG	3	7 3/4
US631	2.0	BAYONET	INS	IRRG	3	7 3/4
US632	0.5	BAYONET	INS	IRRG	4 1/2	8 1/2
US633	0.7	BAYONET	INS	IRRG	4 1/2	8 1/2
US634	1.0	BAYONET	INS	IRRG	4 1/2	8 1/2
US635	1.5	BAYONET	INS	IRRG	4 1/2	8 1/2
US636	2.0	BAYONET	INS	IRRG	4 1/2	8 1/2
US637	0.5	BAYONET	INS	IRRG	5	9
US638	0.7	BAYONET	INS	IRRG	5	9
US639	1.0	BAYONET	INS	IRRG	5	9
US640	1.5	BAYONET	INS	IRRG	5	9
US641	2.0	BAYONET	INS	IRRG	5	9
US532	0.5	BAYONET	INS	IRRG	6	10
US642	0.7	BAYONET	INS	IRRG	6	10
US643	1.0	BAYONET	INS	IRRG	6	10
US174	1.5	BAYONET	INS	IRRG	6	10
US533	2.0	BAYONET	INS	IRRG	6	10



US174-US533

CUSHING Irrigating Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US369	0.7	BAYONET	INS	IRRG	3 1/2	7 1/2
US190	1.0	BAYONET	INS	IRRG	3 1/2	7 1/2
US191	1.5	BAYONET	INS	IRRG	3 1/2	7 1/2
US370	2.0	BAYONET	INS	IRRG	3 1/2	7 1/2



US190-US370



US190-US370

ELECTROSURGERY IRRIGATING BIPOLAR FORCEPS

HARDY Irrigating Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US366	0.5	BAYONET	INS	IRRG	4 1/2	8 1/2
US382	1.0	BAYONET	INS	IRRG	4 1/2	8 1/2
US385	1.5	BAYONET	INS	IRRG	4 1/2	8 1/2
US192	2.0	BAYONET	INS	IRRG	4 1/2	8 1/2
US193	0.5	BAYONET	INS	IRRG	6	10
US194	1.0	BAYONET	INS	IRRG	6	10
US195	1.5	BAYONET	INS	IRRG	6	10
US558	0.5	BAYONET	NON-INS	IRRG	4 1/2	8 1/2

US192-US195, US382, US385, US558

US192-US195, US382, US385, US558

SCOVILLE Irrigating Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US367	1.5	BAYONET	INS	IRRG	2 3/4	7 3/4
US556	1.5	BAYONET	NON-INS	IRRG	2 3/4	7 3/4

US367, US556

GERALD Irrigating Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US368	1.0	BAYONET	INS	IRRG	3 1/2	7 1/2
US557	1.0	BAYONET	NON-INS	IRRG	3 1/2	7 1/2

US368, US557

ELECTROSURGERY IRRIGATING BIPOLAR FORCEPS

MALIS Style Titanium Irrigating Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US360	STR	1.5	BAYONET	INS	IRRG	3 3/4	8 3/4
US361	STR	0.5	BAYONET	INS	IRRG	3 3/4	8 3/4
US362	STR	1.0	BAYONET	INS	IRRG	3 3/4	8 3/4
US363	ANG-UP	0.5	BAYONET	INS	IRRG	3 3/4	8 3/4
US372	ANG-DOWN	0.5	BAYONET	INS	IRRG	3 3/4	8 3/4



US360-US362



US360-US362



US363, US372



US363, US372

MALIS Style Rose Gold Irrigating Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US166	STR	0.5	BAYONET	INS	IRRG	3 3/4	8 3/4
US167	STR	1.0	BAYONET	INS	IRRG	3 3/4	8 3/4
US168	STR	1.5	BAYONET	INS	IRRG	3 3/4	8 3/4
US169	ANG-UP	0.5	BAYONET	INS	IRRG	3 3/4	8 3/4
US170	ANG-DOWN	0.5	BAYONET	INS	IRRG	3 3/4	8 3/4



US166-US168



US169-US170




ELECTRO
SURGERY

ELECTROSURGERY IRRIGATING BIPOLAR FORCEPS

Irrigating Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK838R	0.7	BAYONET	INS	IRRG	2	6 1/2
GK840R	0.4	BAYONET	INS	IRRG	2 3/4	7
GK841R	0.7	BAYONET	INS	IRRG	2 3/4	7
GK842R	1.3	BAYONET	INS	IRRG	2 3/4	7
GK843R	0.4	BAYONET	INS	IRRG	3 1/2	8
GK844R	0.7	BAYONET	INS	IRRG	3 1/2	8
GK845R	1.3	BAYONET	INS	IRRG	3 1/2	8
GK846R	0.4	BAYONET	INS	IRRG	4 1/2	9
GK847R	0.7	BAYONET	INS	IRRG	4 1/2	9
GK848R	1.3	BAYONET	INS	IRRG	4 1/2	9



ELECTROSURGERY DISPOSABLE BIPOLAR FORCEPS

Stainless Steel Disposable Bipolar Forceps, Box of 10

Catalog Number	Tip Size (mm)	Shape	Style	Insulation	Irrigation	Working Length (in.)	Length (in.)
US153SU	1.5	BAYONET	SCOVILLE	INS	NON-IRRIG	2 3/4	7 3/4
US567SU	2.0	BAYONET	CUSHING	INS	NON-IRRIG	2 3/4	7 3/4
US562SU	0.7	BAYONET	CUSHING	INS	NON-IRRIG	3 1/2	7 1/2
US560SU	1.0	BAYONET	GERALD	INS	NON-IRRIG	3 1/2	7 1/2
US154SU	0.5	BAYONET	HARDY	INS	NON-IRRIG	4 1/2	8 1/2
US155SU	1.0	BAYONET	HARDY	INS	NON-IRRIG	4 1/2	8 1/2
US156SU	1.5	BAYONET	HARDY	INS	NON-IRRIG	4 1/2	8 1/2
US157SU	2.0	BAYONET	HARDY	INS	NON-IRRIG	4 1/2	8 1/2



US153SU-US567SU

1/1

US153SU-US567SU



US158SU

1/1

US158SU

ADSON Disposable Bipolar Forceps, Stainless Steel, Box of 10

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US158SU	1.0	STR	INS	NON-IRRIG	5

Disposable Jewelers Bipolar Forceps, Stainless Steel, Box of 10

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US568SU	0.4	STR	NON-INS	NON-IRRIG	4
US569SU	0.4	STR	INS	NON-IRRIG	4



US568SU-US569SU

ELECTROSURGERY DISPOSABLE BIPOLAR FORCEPS

Rose Gold Disposable Jewelers Bipolar Forceps, Box of 5

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US408SU	0.4	STR	INS	NON-IRRG	4
US409SU	1.0	STR	INS	NON-IRRG	5



US408SU-US409SU

Rose Gold Disposable Bipolar Forceps, Box of 5

Catalog Number	Tip Size (mm)	Shape	Style	Insulation	Irrigation	Working Length (in.)	Length (in.)
US565SU	0.7	STR	SEMKIN	INS	NON-IRRG		5 1/2
US559SU	1.5	BAYONET	SCOVILLE	INS	NON-IRRG	2 3/4	7 3/4
US201SU	0.5	BAYONET	CUSHING	INS	NON-IRRG	3 1/2	7 1/2
US561SU	0.7	BAYONET	CUSHING	INS	NON-IRRG	3 1/2	7 1/2
US202SU	1.0	BAYONET	CUSHING	INS	NON-IRRG	3 1/2	7 1/2
US203SU	1.5	BAYONET	CUSHING	INS	NON-IRRG	3 1/2	7 1/2
US566SU	2.0	BAYONET	CUSHING	INS	NON-IRRG	3 1/2	7 1/2
US204SU	0.5	BAYONET	HARDY	INS	NON-IRRG	4 1/2	8 1/2
US563SU	0.7	BAYONET	HARDY	INS	NON-IRRG	4 1/2	8 1/2
US205SU	1.0	BAYONET	HARDY	INS	NON-IRRG	4 1/2	8 1/2
US206SU	1.5	BAYONET	HARDY	INS	NON-IRRG	4 1/2	8 1/2
US564SU	2.0	BAYONET	HARDY		NON-IRRG	4 1/2	8 1/2



US201SU-US564SU



US565SU

ELECTROSURGERY DISPOSABLE BIPOLAR FORCEPS

Rose Gold Disposable Irrigating Bipolar Forceps, Box of 5

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US410SU	0.5	BAYONET	INS	IRRG	4 1/2	8 1/2
US411SU	1.0	BAYONET	INS	IRRG	4 1/2	8 1/2
US412SU	1.5	BAYONET	INS	IRRG	4 1/2	8 1/2



US410SU-US412SU

Stainless Steel Disposable Irrigating Bipolar Forceps, Box of 5

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US523SU	0.5	BAYONET	INS	IRRG	4 1/2	8 1/2
US524SU	1.0	BAYONET	INS	IRRG	4 1/2	8 1/2
US525SU	1.5	BAYONET	INS	IRRG	4 1/2	8 1/2



US523SU-US525SU

Specialty Disposable Non-Stick Bipolar Forceps


Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US501SU	0.5	BAYONET	INS	NON-IRRG	7
US502SU	1.0	BAYONET	INS	NON-IRRG	7
US503SU	1.5	BAYONET	INS	NON-IRRG	7
US504SU	0.5	BAYONET	INS	NON-IRRG	8
US505SU	1.0	BAYONET	INS	NON-IRRG	8
US506SU	1.5	BAYONET	INS	NON-IRRG	8
US507SU	0.5	BAYONET	INS	NON-IRRG	9
US508SU	1.0	BAYONET	INS	NON-IRRG	9
US509SU	1.5	BAYONET	INS	NON-IRRG	9



US501SU-US509SU






US501SU-US509SU



ELECTROSURGERY BIPOLAR CONNECTING CORDS

Bipolar Cords for BipoJET[®]

Catalog Number	New Cord	Pin Distance	Compatible Generators
GN145		22 MM	EMC DOLLEX ERBE Int'l MARTIN
GN085		29 MM	AESCLAP VALLEYLAB ERBE
GN146		38 MM	CODMAN CMC II / III ERBE MARTIN




LANDOLT Bipolar Cord and Adapters

Catalog Number	Pin Distance (mm)	Description
GK195	28.6	REUSABLE CORD F/ INSULATED LANDOLT FORCEPS
GN219	32	ADAPTER F/ CONMED
GN224	38	ADAPTER F/ CODMAN

ELECTROSURGERY BIPOLAR CONNECTING CORDS

Dual Bipolar Cord, Sterile Disposable

Catalog Number	Packaging	Length (ft)
US349SP	10/BOX	12



US349SP

Round Bipolar Cord, Reusable

Catalog Number	Material	Length (ft)
US351	SILICONE	12



US351

Dual Bipolar Cord, Reusable

Catalog Number	Length (ft)
US359	12



US359

ELECTROSURGERY BIPOLAR CONNECTING CORDS

Combined Bipolar and Irrigation Cord

Catalog Number	Pin Type	Packaging	Compatible w/
US352	DUAL	BOX OF 10	-
US353	ROUND	BOX OF 10	-
US355	DUAL	BOX OF 10	MALIS 1000 IRRIGATION SYSTEM
US356	ROUND	BOX OF 10	MALIS 1000 IRRIGATION SYSTEM

