

Linked Page



Certified Register



Wholesale Supplier
Fumigator / O.T.Care / Fogger
A to Z Hospital Equipment and
Surgical Products B to C



 **8275140877**

 **9822688071**

 **info@amramedtech.com**



**Sr No.29/1/2, Adarsh Nagar,
Uruli Devachi, Hadapsar, Pune**
Branches :- All Maharashtra



amramedtech.com



O.T.FUMIGATOR FUEMES

Specification

AMRAMED Make O.T.Fumigator.. WORKING OF FUMIGATOR MACHINE:- PROTECTION :- INDICATION:- M S Powder Quoted with S.S.Top & 100% Copper Bowl Vaporizer Heating Electrode to make Hot Air Fumes of Disinfected Chemical+Water Inbuilt Auto Cut-Off Timer

We Offer Latest O.T.Fumigator Machine Heating Electrode to make Hot Air Fumes Of Disinfected [Formalin] Liquid.Hot Air Fumes will be Spread in Surrounding Area.Here,No Requirements of to add a potassium Permanganate .If we mix a potassium permanganate with formalin liquid,after use its become a Black Coal & Its Difficult to Clean,Maintain Container.It Covers 1000 Sq.Feet Area.

DO's:-Add Formalin 40 ml or Any Type of Chemical to Suitable for Fumigation Purpose 40 ml + 40 ml to 60 ml Water, Just Switch ON then Push Process Switch..

DON'T:-Don't add a potassium permanganate.. DON'T SPRAY CHEMICALS MORE THAN 90 ML.IN TIME INTERVALS OF 24 HOUR'S

WORKING OF FUMIGATOR MACHINE:-

As soon as Fumigator machine is Switched ON, after 5 to 10 Minutes fumigation Liquid starts boiling with help of carbon electrodes inside. It will make hot air fumes and will be Spread in Surrounding Area.

PROTECTION :-

After over of one hour time Fumigator Machine will be switched OFF itself automatically.

INDICATION:-

- Fumigator ON Indcation [Red Led]
- Timer Count ON Indication [Green Led]
- Process ON Fumigation Indication [Blue Led]

NOTE:- After over of 1 hour time above Indications[Green & Blue LED] are Turned OFF Automatically. DON'T SPRAY ANY CHEMICALS MORE THAN 90 ML IN TIME INTERVALS OF 24 HOURS.

ADVANTAGES

- Total Capacity 75/100 Watt Only.
- Power Required 230 VAC/50Hz,+/- 25%.
- Supply Through 3 Pin Mains Cord.
- Inbuilt 1 Hour's Programmed Timer.
- Attractive M.S. Powder Quoted Body with Mirror Finishing S S Plate & 100% Copper Bowl.
- No Countable Electricity Power Consumption.
- Easy to Clean Just add water [Pitambari] & shake the Bowl & remove water
- Light Weight Easy to Carry.
- AREA COVER 1000 SQ.FEET. SINGLE PROCESS.FOR 48 HOUR'S.

AMRA-FF-010



COMPARISON CHART



Indian Fumigation Unit		Existing Model O.T. Care (Doom Type)	Existing Model Fogger Machine
Type of Fumigation	Indian Fumigation Unit	O.T. Care (Doom Type)	Fogger Machine
Rate/Price	Affordable	12000/- to 18000/-	2000/- to 28000/-
Manufacturing	Pune Maharashtra INDIA	Old Technology	China Technology
Working Technology	P.C.B. Base	Electric Motor	Electric Motor With Nasal
Water Use	50 to 75 ML	3 Lit To 5 Lit	1 Lit To 5 Lit
Chemical Name	Pure Formaline 40 ML + 60 ML Water Bailocide, Ecoshield, Envishield, Hydrogen, & Silvermitrate	Potashiyam Pharma Grate & Formaldehyde ETC	Hydrogen Peroxide ETC
Maintance	10% to 20%	10% to 60	60% To 80%
Timer	Inbuildt P.C.B.	Auto Timer Optional (Extra) 2000/- To 3500/-	Auto Timer Optional (Extra) 2500/- To 4500/-
Property Loss	Never (No)	YES 100%	YES 100%
Electricity Consumption (Light Bill)	Not Countable 40 paise	20/- to 30/-	30/- To 50/-
After Using Fumigation	No Physical Work	Furniture & Carpet Swapping Compulsory (Black Spot)	Furniture & Carpet Swapping Compulsory (Black Spot)
Operating (Fumigate)	Any Person System	Trainee or Knowledgeable	Proper Train Person
Working System (Fumigate)	Vaporizing System	Fume's (Drop Spared)	FOG
Cover Area	100 Sq. Ft To 1000 Sq. Ft.	500 Sq. Ft. To 1000 Sq. Ft.	500 Sq. Ft. To 1000Sq. Ft.
Dilution Area Wise	Yes Changes	No Changes	No Changes
Guaranty / Warranty	12 Month Guaranty 12 Month Warranty / 24 Month	12 Month Guaranty	12 Month Guaranty

Uses Full Space: Operation Theater/ICU/NICU/Patient Ward/OPD/Clinic/Laboratory
/ Diagnostic Center / Animal Farm House / Mandir/Loadge/Hotel/Factory/Workshop
/ Godown/College/School/Classes/Transport Taxi/Tempo/Bus/Shopping Mall
/ Marriage & Cinema Hall / Bank/Office/Cabin

O.T.FUMIGATOR FOGGER STAILESS STEEL BODY 3 LTR.

Micro Fogger provides excellent control size, even for near-invisible 7-10 micron droplets needed by concentrated ULV chemicals. The Micro Fogger can treat 2,000-4,000 cubic feet per minute, reducing application time and labor costs and when you need high output and large droplets, the Micro Fogger is just as effective. It can deliver up to 150 ml/min in a mist of large (30 micron) droplets used in humidification, dust abatement and sanitizing work. The advanced Micro Fogger can apply either oil-or water based chemicals and it's easy to calibrate when changing liquids.

AMRA-FSS-060

- Features: a) Max. Area Coverage As Compared To Any Other Conventional Fogger. 304 Grade Stainless Steel
- Bod; b) Chemical Level Indicator.
- Valve c) Comes With Pre-Fitted Liquid Flow Control
- d) 100% Convenient, Portable & Easy To Use



- Body & Tank: 304 Grade Stainless Steel
- Tank Capacity: 4 Litter of Chemical Storage
- Blower Motor: >26000 RPM Motor
- Nozzle: 1 Counter Rotating Vortex Nozzle
- Chemicals: Can atomize both oil-based and water-based liquids.
- Control Valve: On turn. Brass Valve
- Particle Size: 5-30p, adjustable. Liquid viscosity and density can affect particle size.
- Range: Visible fog, 20-30 ft (7.5m)
- Input: 230V 50Hz Single Phase AC Supply
- Dimensions: (33 x 41 x 23.5cm) L x H x Dia
- Net. Weight: 6.50 Kg (Approx)

AMRA-FSS-070



ULV MICRO FOGGER

Capacity

