

Linked Page



Certified Register



Wholesale Supplier
O.T Pendant
A to Z Hospital Equipment and
Surgical Products **B to C**



 **8275140877**

 **9822688071**

 **info@amramedtech.com**



**Sr No.29/1/2, Adarsh Nagar,
Uruli Devachi, Hadapsar, Pune**
Branches :- All Maharashtra



amramedtech.com

O.T PENDANT

Product Description

- The activity range of Horizontal single arm radius 800mm, horizontal arm and functional box can be simultaneously rotated, each horizontal rotation angle 340°.
- Functional box vertical lifting stroke 300 mm. Motor power 100W.

- The gas terminal interface standard configuration: oxygen 2pcs, Vacuum 2pcs; 2pcs Air4 bar

- i.) Different colors and shapes,
- ii.) more than 20,000 times mating cycles;
- iii.) using the secondary seal, with three states (on, off, pull) can carry air service;

Instrument platform: two layers (height adjustable) (Size 470x370x50mm); rounded bumper design, each platform can load 20kg;


- Power supply: AC220V, 50HZ; AC110V, 60HZ
- Electrical outlet: 10 pcs (220V / 6/16A);
- Potential ground terminal: 2pcs;
- The network interface: RJ45 1pc;
- Telecom Interface Rj11 , 1pc;
- Drawer: Size 310x310x125 mm 1pc;
- Weight loading: 150kg;

Syringe pump rack 1pc, with 4 hook IV Stand

- The main material: high-strength aluminium alloy;

The surface treatment using imported high-voltage electrostatic powder spray

AMRA - 01



AMRA - 02



O.T PENDANT

Product Description

- The activity range of Horizontal single arm radius 800mm, horizontal arm and functional box can be simultaneously rotated, each horizontal rotation angle 340°.
- The gas terminal interface standard configuration: oxygen 2pcs, vacuum 2pcs; 2pcs Air4 bar
- I.) Different colors and shapes,
- ii.) more than 20,000 times mating cycles;
- iii.) using the secondary seal, with three states (on, off, pull) can carry air service.
- Instrument platform: two layers (height adjustable) (Size 470x370x50mm); rounded
- bumper design, each platform can load 20kg;
- Power supply: AC220V, 50HZ; AC110V, 60HZ
- Electrical outlet: 10 pcs (220V / 6/16A);
- Potential ground terminal: 2pcs;
- The network interface: RJ45 1pc;
- Drawer: Size 310x310x125 mm 1pc;
- Weight loading: 150kg;
- Syringe pump rack 1pc, with 4 Hook IV Stand.
- The main material: high-strength aluminium alloy;
- The surface treatment using imported high-voltage hdant electrostatic powder spray;

AMRA - 03



AMRA - 04



AMRA - 05



O.T PENDANT ENDOSCOPY/LAPROSCOPY

Technical Parameters

- Working power: AC220V, 50Hz , AC110V, 60Hz
- Cross arm range of motion (radius): 700-1800mm
- Horizontal rotation angle: 0 ~ 340 °, the cross arm and the terminal box can be rotated horizontally or simultaneously.
- The net load weight is 150kg.
- Instrument platform: 3 layers (height adjustable) optional Drawer size 350x350x125mm, Tray size 505x425x50mm rounded bumper design.

AMRA - O6

Standard Configuration

- Gas interface configuration (2 oxygen, 2 vacuum, 2 compressed air):
 - a. The interface color and shape are different, with anti-missing function;
 - b. The number of insertions and removals is more than 20,000 times;
 - c. Use secondary sealing, with three states (on, off, pull), with gas maintenance;
 - Power socket: 10, 220V, 6/16A;
 - Equipotential grounding terminal: 2;
 - One stainless steel adjustable infusion rod holder with 4 hook IV stand.
 - The main material is made of high-strength aluminum alloy profiles;
 - The surface treatment is electrostatic spraying; ceiling-mounted installation, stable and firm.
- One storage item.



O.T Pendant Floor Model

Standard Configuration

- Functional box vertical lifting stroke 300mm. • Motor power 100W.
 - The gas terminal interface standard configuration (Outlets): vaccum 2pcs, oxygen 2 pcs, compressed air 2 pcs.
 - I.) Different colors and shapes,
 - ii.) more than 20,000 times mating cycles.
 - iii.). using the secondary seal, with three states (on, off, pull), can carry compressed air service.
 - Instrument platform: two layers (height adjustable), rounded bumper design, each platform can load 20kg.
 - Weight loading: 150kg.
 - Power supply: AC220V, 50HZ; AC110V, 60HZ
 - Electrical outlet: 10 pcs, 6 Amp+16 Amp
 - Potential ground terminal: 2pcs.
- % The network interface : RJ45 1pc; telecom connector (RJ11 1 PC;)
- 6 Drawer Size: 470x370x50mm 1 pc.
- Syringe pump rack with 4 Hook IV Stand
 - The main material: high-strength aluminium alloy.
 - The surface treatment using imported high-voltage electrostatic powder spray.
 - 4 Double caster Wheel (2 wheels Lockable)

AMRA - 07

