



Wholesale Supplier

Orthopedics_Spine_instrument

A to Z Hospital Equipment and Surgical Products **B to C**

Orthopedics/Spine



Scissors	880 - 885
Plaster/Nail/Wire Instruments	886 - 898
Meniscus Instruments	899
Forceps	900 - 916
Retractors	917 - 956
Rasps	957 - 959
Elevators and Bone Levers	960 - 973
Chisels/Gouges/Osteotomes	974 - 987
Mallets/Impactors/Bone Tamps	988 - 994
Drills/Burr/Reamers	995 - 996
Reamers and Trephines	997
Curettes	998 - 1013
Microdiscectomy Instruments	1014 - 1018
Dissectors and Probes	1019
Rongeurs	1020 - 1049
Kerrison Bone Punches	1050 - 1057
KAIRison Pneumatic Bone Punch	1058 - 1059
Bone Cutters	1060 - 1063
Rulers, Calipers and Gauges	1064 - 1065

8275140877

9822688071

info@amramedtech.com



**Sr No.29/1/2, Adarsh Nagar,
Uruli Devachi, Hadapsar, Pune**
Branches :- All Maharashtra





amramedtech.com

ORTHOPEDICS/SPINE SCISSORS

Universal Bandage Scissors

Catalog Number	Handle Color	Length (in.)
BC881R	BLACK	7
BC980R	RED	7
BC981R	BLUE	7
BC983R	YELLOW	7
BC985R	ORANGE	7
BC987R	PINK	7
BC988R	LIGHT GREEN	7
BC989R	VIOLET	7



LISTER Universal Bandage Scissors

Catalog Number	Handle Color	Length (in.)
BC885R	BLACK	5 3/4
BC970R	RED	5 3/4
BC971R	BLUE	5 3/4
BC973R	YELLOW	5 3/4
BC975R	ORANGE	5 3/4
BC977R	PINK	5 3/4
BC978R	LIGHT GREEN	5 3/4
BC979R	VIOLET	5 3/4


ORTHOPEDICS/SPINE SCISSORS




BC876R, MB146R

LISTER Bandage Scissors, Serrated


Catalog Number	Length (in.)
MB146R	3 1/2
BC876R	6 1/4



BC876R



MB146R





AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C




ORTHOPEDICS/SPINE SCISSORS




BC848R-BC863R

1/2




BC848R

1/1




BC849R

1/1




BC860R

1/1




BC861R

1/1



BC862R

1/1



BC863R

1/1



MB223

1/2



MB223


1/1

LISTER Bandage Scissors

Catalog Number	Blade	Length (in.)
BC848R	LONG	6 1/8
BC849R	LONG	7
BC860R	STANDARD	4 1/2
BC861R	STANDARD	5 1/2
BC862R	STANDARD	7
BC863R	STANDARD	7 7/8

LISTER Bandage Scissors W/ One Large Ring

Catalog Number	Length (in.)
MB223	8



AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



ORTHOPEDICS/SPINE SCISSORS

Plaster and Bandage Scissors

Catalog Number	Length (in.)
LX558R	9 1/4



LX558R

Plaster Scissors

Catalog Number	Length (in.)
LX548R	6 1/4



LX548R


BRUNS Plaster Scissors

Catalog Number	Length (in.)
LX556R	9 1/4



LX556R

LX556R




ORTHOPEDICS/SPINE SCISSORS



BERGMANN Plaster Scissors

Catalog Number	Length (in.)
LX553R	9



ESMARCH Plaster Scissors

Catalog Number	Length (in.)
LX550R	7 7/8
LX552R	9



HERCULES Plaster Shears


Catalog Number	Length (in.)
MB106R	7 1/2



ORTHOPEDICS/SPINE SCISSORS

Plaster Scissors

Catalog Number	Length (in.)
BC888R	8






SEUTIN Plaster and Bandage Scissors

Catalog Number	Length (in.)
LX565R	9



RESANO Cartilage Scissors

Catalog Number	Tip	Blade	Length (in.)
BC595R	B/B	ANG	9 3/4



ORTHOPEDICS/SPINE PLASTER AND NAIL INSTRUMENTS



LX600R

HENNIG Plaster Spreader


Catalog Number	Length (in.)
LX600R	10 5/8



LX605R

Plaster Spreader

Catalog Number	Length (in.)
LX605R	9



LX618R-LX624R

WOLFF Plaster Breaker


Catalog Number	Length (in.)
LX618R	7
LX624R	9 1/2



HF213R

Nail Cutter

Catalog Number	Length (in.)
HF213R	5 1/4



HF213R


ORTHOPEDICS/SPINE NAIL AND WIRE INSTRUMENTS



BC521R

SYSTRUNK Nail Splitting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC521R	S/S	STR	5 1/8



BC521R



HF021R

Nail Scissors

Catalog Number	Tip	Blade	Length (in.)
HF021R	S/S	CVD	4 1/4



HF021R



MB370R

Needle Nose Pliers/Wire Cutters

Catalog Number	Length (in.)
MB370R	7 1/4




MB370R



ORTHOPEDICS/SPINE WIRE INSTRUMENTS



BM178R



BM178R

DUROGRIP™ TC AESCULAP Needle Holder, Double Action

Catalog Number	Handle	Jaw	Length (in.)
BM178R	RATCHET	SERR	7 7/8




LX150R



LX150R

Wire Cutter, W/ Specially Hardened Cutting Edges


Catalog Number	For Wire Size	Length (in.)
LX150R	SOFT, UP TO 2MM AND HARD, UP TO 1.5MM	5 3/4
LX152R	SOFT, UP TO 3.5MM AND HARD, UP TO 2.0MM	8 3/4
LX156R	SOFT, UP TO 3MM AND HARD, UP TO 2.2MM	9 1/4



LX152R



LX152R



LX156R



LX156R



AMRAMED
TECH Pvt. Ltd.


Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

ORTHOPEDICS/SPINE WIRE INSTRUMENTS



LX157R

HARVEY Wire Cutter, W/ Specially Hardened Cutting Edges

Catalog Number	For Wire Size	Tip	Length (in.)
LX157R	SOFT, UP TO 3MM AND HARD, UP TO 2.2MM	ANG	9 1/4




LX158R

Wire Cutter, Front and Lateral, W/ Specially Hardened Cutting Edges

Catalog Number	Wire Size	Length (in.)
LX158R	SOFT, UP TO 2MM AND HARD, UP TO 1.5MM	7



LX158R



ORTHOPEDICS/SPINE WIRE INSTRUMENTS

**REILL Cutter, F/ Wire Close To Bone, W/
Specially Hardened Cutting Edges**

Catalog Number	For Wire Size	Length (in.)
LX159R	SOFT, UP TO 2MM AND HARD, UP TO 1.5MM	6 7/8

**Front Cutter, W/ Specially Hardened
Cutting Edges**

Catalog Number	For Wire Size	Length (in.)
LX163R	SOFT, UP TO 3.5MM AND HARD, UP TO 2.5MM	6 1/2

**TC Wire Cutter W/ Carbide Cutting
Edges**

Catalog Number	Length (in.)
LX162R	9 1/2



LX159R



LX159R



LX163R



LX163R



LX162R



AMRAMED
TECH Pvt. Ltd.


Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

ORTHOPEDICS/SPINE WIRE INSTRUMENTS





MG025R

TC Wire Cutter Double Action

Catalog Number	Packaging	Tip	Length (in.)
MG025R	2.5	ANG	8 3/4




MG025R




MG029R

TC Flush Wire Cutter

Catalog Number	Packaging	Tip	Length (in.)
MG029R	1.6	FRONT AND LATERAL	6 3/4



ORTHOPEDICS/SPINE WIRE INSTRUMENTS




Wire Holding Forceps W/Hard Metal Inserts

Catalog Number	Length (in.)
LX176R	6 3/4



Wire Grasping Forceps W/Hard Metal Inserts

Catalog Number	Tip	Length (in.)
LX179R	FLAT	7 1/2



ORTHOPEDICS/SPINE WIRE INSTRUMENTS



MG088R

HERCULES Pin Cutter

Catalog Number	Packaging	Description	Length (in.)
MG088R	3	TUNGSTEN CARBIDE	8 3/4




MG089R

Pin Cutter, End-Cutting

Catalog Number	Packaging	Length (in.)
MG089R	4.7	14 3/4
MG090R	6.3	22 1/2




MG090R



ORTHOPEDICS/SPINE WIRE INSTRUMENTS

Bolt Cutter


Catalog Number	Item	Length (in.)
LX164R	F/KIRSCHNER WIRES, TIBIA HEAD SCREWS, BOLTS, PINS, UP TO 6MM	18 1/2



LX164R

Bolt Cutter

Catalog Number	Item	Length (in.)
LX165R	F/KIRSCHNER WIRES, TIBIA HEAD SCREWS, BOLTS, PINS, UP TO 6MM	12 1/2





LX165R



LX165R

ORTHOPEDICS/SPINE WIRE INSTRUMENTS



Wire Holding Pliers


Catalog Number	Tip	Length (in.)
LX170R	SERR	6 3/4
LX171R	CROSS-GROOVES	6 3/4

Wire Holding Forceps


Catalog Number	Tip	Length (in.)
LX172R	SERR	7 1/2

Wire Holding Pliers

Catalog Number	Tip	Tip Size (mm)	Length (in.)
MG026R	LOCKING	7	7




ORTHOPEDICS/SPINE WIRE INSTRUMENTS




MG028R

Wire Holding Needle-Nose Pliers

Catalog Number	Tip	Tip Size (mm)	Length (in.)
MG028R	LOCKING, NEEDLE-NOSE	5	9 1/2



MG028R



LX174R

Vise Grip Pliers

Catalog Number	Tip	Length (in.)
LX174R	FLAT NOSE	7 7/8



LX174R


ORTHOPEDICS/SPINE WIRE INSTRUMENTS



LX175S

Bending Iron


Catalog Number	FOR BORING WIRE SIZE	Length (in.)
LX175S	UP TO 2.7MM	4 3/4



LX177R-LX178R


Wire Extraction Pliers

Catalog Number	Length (in.)
LX177R	5 1/8
LX178R	7



LX177R

LX178R



LX181R-LX185R


Wire Extraction Pliers

Catalog Number	Tip	Length (in.)
LX181R	RND	7
LX182R	RND, SHORT JAW	7
LX183R	FLAT	7
LX185R	PARALLEL	5 3/4




LX181R

LX182R




LX183R

LX185R



LX183R

LX185R



ORTHOPEDICS/SPINE WIRE INSTRUMENTS





Utility Pliers

Catalog Number	Length (in.)
MB576R	7 3/4




Flat Nosed Pliers

Catalog Number	Tip	Length (in.)
LX189R	PARALLEL	7 1/4



Wire Bending Pliers

Catalog Number	For Wire Size
LX190R	UP TO 1.8MM





ORTHOPEDICS/SPINE MENISCUS INSTRUMENTS



SMILLIE Meniscotome


Catalog Number	Jaw	Length (in.)
F0857R	STRAIGHT	6 3/4
F0855R	LEFT	6 3/4
F0856R	RIGHT	6 3/4



ORTHOPEDICS/SPINE FORCEPS

BIRCHER Cartilage Forceps


Catalog Number	Jaw	Length (in.)
F0800R	STR	7 7/8
F0801R	CVD	7 7/8



F0800R



F0800R



F0801R





F0801R

ORTHOPEDICS/SPINE FORCEPS



WELLER Cartilage Forceps

Catalog Number	Jaw	Length (in.)
F0807R	CVD	7 7/8



BIRCHER-GANSKE Cartilage Forceps


Catalog Number	Jaw	Length (in.)
F0808R	SLIGHT CVD	7 7/8



ORTHOPEDICS/SPINE FORCEPS


BIRCHER-GANSKE Cartilage Forceps (continued)

Catalog Number	Jaw	Length (in.)
F0810R	CVD	7 7/8



MARTIN Cartilage Forceps


Catalog Number	Jaw	Length (in.)
F0816R	STR	8




ORTHOPEDICS/SPINE BONE HOLDING FORCEPS

MEYER Reposition Forceps, F/ Wires

Catalog Number	For Wire Size	Length (in.)
F0106R	UP TO 1MM	6 3/4
F0108R	UP TO 2MM	8




F0106R




F0106R



F0108R




F0108R



ORTHOPEDICS/SPINE BONE HOLDING FORCEPS



DINGMANN Bone Holding Forceps

Catalog Number	Teeth	Length (in.)
F0110R	2 x 2	7 1/4




DINGMANN Bone Holding Forceps

Catalog Number	Teeth	Shaft	Length (in.)
F0111R	2 x 2	CVD	7 1/4



Bone Holding Forceps


Catalog Number	Length (in.)
F0112R	8 1/4



ORTHOPEDICS/SPINE BONE HOLDING FORCEPS


Finger Reposition Forceps, Stepped Tips

Catalog Number	Length (in.)
F0113R	3 3/4




Finger Reposition Forceps, F/ Small Fragments

Catalog Number	Length (in.)
F0115R	5 1/4



REILL Reposition Forceps

Catalog Number	Length (in.)
F0116R	6 3/4



ORTHOPEDICS/SPINE BONE HOLDING FORCEPS



Forceps for Maxillary Reposition

Catalog Number	Length (in.)
F0117R	5 1/4




Bone Reposition Forceps, Double Ratchet

Catalog Number	Length (in.)
F0118R	8



Bone Holding Forceps

Catalog Number	Length (in.)
F0120R	7 7/8



ORTHOPEDICS/SPINE BONE HOLDING FORCEPS

Bone Reduction Forceps

Catalog Number	Length (in.)
F0123R	5 1/2




F0123R



F0123R

Bone Reduction Forceps, W/ Screw Lock and Thread Fixation


Catalog Number	Length (in.)
F0124R	6 3/4
F0125R	9 1/4



F0124R



F0125R



ORTHOPEDICS/SPINE BONE HOLDING FORCEPS



F0127R



F0127R

LANGENBECK Bone Holding Forceps

Catalog Number	Tip Size (mm)	Length (in.)
F0127R	11	8




F0130R



F0130R

FERGUSON Bone Holding Forceps

Catalog Number	Tip Size (mm)	Length (in.)
F0130R	12	8



F0135R



F0135R

SEMB Rib/Bone Holding Forceps

Catalog Number	Tip Size (mm)	Item	Length (in.)
F0135R	8	W/O RATCHET	7 1/2
F0136R	8	W/ RATCHET	7 1/2



F0136R



F0136R

ORTHOPEDICS/SPINE BONE HOLDING FORCEPS


FARABEUUF Bone Holding Forceps

Catalog Number	Tip Size (mm)	Length (in.)
F0140R	9	9
F0141R	11	10



F0140R

1/2



1/1

Jaw Adustable



F0140R



F0141R

1/2



F0141R

FARABEUUF-LAMBOTTE Forceps, W/ Ratchet

Catalog Number	Tip Size (mm)	Length (in.)
F0143R	11	10




F0143R

1/2



F0143R



ORTHOPEDICS/SPINE BONE HOLDING FORCEPS



**KERN Mini Bone Holding Forceps,
W/ Ratchet**

Catalog Number	Tip Size (mm)	Length (in.)
FO148R	7.5	6 1/8



KERN-LANE Bone Holding Forceps

Catalog Number	Tip Size (mm)	Length (in.)
FO150R	9	9 1/2




**KERN-LANE Bone Holding
Forceps, W/ Ratchet**

Catalog Number	Tip Size (mm)	Length (in.)
FO158R	9	9 1/2

ORTHOPEDICS/SPINE BONE HOLDING FORCEPS

LANE Bone Holding Forceps

Catalog Number	Tip Size (mm)	Description	Length (in.)
F0151R	9	W/O RATCHET	13
F0159R	13	W/ RATCHET	12 1/2



ORTHOPEDICS/SPINE BONE HOLDING FORCEPS

Patella Forceps W/ Screw Lock

Catalog Number	Length (in.)
F0153R	7 1/4



F0153R 1/1

LAMBOTTE Bone Holding Forceps, W/ Movable Jaw and Ratchet

Catalog Number	Tip Size (mm)	Length (in.)
F0160R	10	8 3/4
F0161R	13	10 1/4
F0162R	16	13



ORTHOPEDICS/SPINE BONE HOLDING FORCEPS



F0169R-F0173R

Self-Centering Forceps, W/ Screw Lock
and Thread Fixation


Catalog Number	Tip Size (mm)	Length (in.)
F0169R	3/3	5 3/4
F0170R	7 / 2.5	7 1/2
F0171R	10.5 / 3.5	9 1/2
F0172R	11 / 3.5	10 1/4
F0173R	11.5 / 4.5	10 3/4




F0169R




F0170R




F0171R



F0172R





F0173R



ORTHOPEDICS/SPINE BONE HOLDING FORCEPS

VERBRUGGE Bone Forceps

Catalog Number	Tip Size (mm)	Length (in.)
F0180R	5	6 7/8
F0181R	10 / 4.0	9 3/4
F0182R	11 / 4.5	10 1/2
F0183R	12 / 4.5	10 5/8



ORTHOPEDICS/SPINE BONE HOLDING FORCEPS



F0190R-F0191R

LOWMAN Bone Clamp

Catalog Number	Length (in.)
F0190R	6 7/8
F0191R	8



F0192R

LOWMAN-GERSTER Bone Clamp


Catalog Number	Length (in.)
F0192R	8 1/4




F0193R

LOWMAN-HOGLUND Bone Clamp

Catalog Number	Length (in.)
F0193R	6 7/8




ORTHOPEDICS/SPINE BONE HOLDING FORCEPS



GERSTER Traction Bar
F/FO190R-F0193R

Catalog Number
F0199R



LEWIN Bone Clamp

Catalog Number	Tip	Length (in.)
MD455	SERR	7

ORTHOPEDICS/SPINE RETRACTORS

25 mm

38 mm

Dimensions for blade size of retractors

1st number = depth

Example: 25 x 38 mm


2nd number = width

Knee Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT419R	39 x 9	8 7/8




BT419R



BT419R




BT420R



BT420R

BLOUNT Knee Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT420R	40 x 4	6 7/8




ORTHOPEDICS/SPINE RETRACTORS



Knee Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT430R	56 x 10	6 3/4



DOANE "Z" Knee Retractor

Catalog Number	Blade Size (mm)	Size
MB976R	55 x 75	MEDIUM
MB977R	66 x 80	LARGE





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier **Certified Register**


A to Z Hospital Equipment and
Surgical Products B to C



ORTHOPEDICS/SPINE RETRACTORS

T-Handle Retractor

Catalog Number	Blade Size (mm)	Length (in.)
MD512R	120 x 32	7
MD513R	120 x 32	11
MD514R	120 x 32	13



MD512R-MD514R

HIBBS Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT503R	25 x 75	9 1/8
BT504R	20 x 50	9 1/8
BT505R	13 x 25	9 1/8



BT503R-BT505R




BT503R



BT504R



BT505R



ORTHOPEDICS/SPINE RETRACTORS


TAYLOR Spinal Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT501R	30 x 70	7 1/4
BT502R	30 x 100	7 7/8

BT501R-BT502R

BT501R






BT502R



ORTHOPEDICS/SPINE RETRACTORS

SMILLIE Retractor


Catalog Number	Blade Size (mm)	Length (in.)
BT421R	32 x 19	5 1/2
BT422R	48 x 19	5 1/2
BT423R	57 x 19	5 1/2
BT424R	52 x 13	5 1/2
BT425R	12 x 13	5 1/2




ORTHOPEDICS/SPINE RETRACTORS

MEYERDING Retractor


Catalog Number	Blade Size (mm)	Length (in.)
BT485R	52 x 16	9
BT486R	76 x 25	9
BT487R	89 x 50	10




BT485R-BT487R



BT485R



BT486R



BT487R



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C




Certified Register

ORTHOPEDICS/SPINE RETRACTORS

AUFRANC Cobra Retractor


Catalog Number	Tip	Length (in.)
MB448R	SMOOTH, POINTED	10
MB514R	SMOOTH, BLUNT	10



ORTHOPEDICS/SPINE RETRACTORS

AUFRANC Cobra Retractor


Catalog Number	Tip	Length (in.)
MB379R	SERR, BLUNT	10



MB379R

1/2

1/1



MB379R

1/1

ORTHOPEDICS/SPINE RETRACTORS



MF058R-MF060R

**CLOWARD Style Cervical Retractor
W/ Lip**


Catalog Number	Width (mm)	Length (in.)
MF058R	18	7
MF059R	20	7
MF060R	23	7



MN073R-MN075R

**CLOWARD Style Cervical Retractor
W/O Lip**


Catalog Number	Width (mm)	Length (in.)
MN073R	18	7
MN074R	20	7
MN075R	23	7



ORTHOPEDICS/SPINE RETRACTORS

HAYES Bone Retractor

Catalog Number	Blade Size (mm)	Length (in.)
FK110R	5.5	6 1/8
FK111R	2.8	5 1/2
FK113R	4.5	6
FK114R	7	6



ORTHO
SPINE

MILLIMETERS

ORTHOPEDICS/SPINE RETRACTORS

BUCK-GRAMCKO Bone Retractor

Catalog Number	Blade Size (mm)	Length (in.)
FK116R	2/7.5	6

HOHMANN Bone Retractor


Catalog Number	Blade Size (mm)	Length (in.)
FK117R	2/6	6 1/4
FK118R	2/8	6 1/4
MB449R	2/15	6

Baby HOHMANN Retractor

Catalog Number	Blade Size (mm)	Length (in.)
FK119R	2/15	5




FK116R 1/1



FK117R 1/1



FK118R 1/1




MB449R




MB449R



FK119R



FK119R





ORTHOPEDICS/SPINE RETRACTORS



HOHMANN Bone Retractor


Catalog Number	Blade Size (mm)	Length (in.)
FH170R	9	8 3/4
FH166R	17	9 1/2



ORTHOPEDICS/SPINE RETRACTORS

HOHMANN Bone Retractor (continued)

Catalog Number	Blade Size (mm)	Length (in.)
FK171R	18	9 1/2
FK176R	25	10 5/8



ORTHOPEDICS/SPINE RETRACTORS

HOHMANN Bone Retractor (continued)


Catalog Number	Blade Size (mm)	Length (in.)
FK172R	24	10 1/4
FK173R	34	11 1/4



FK172R-FK173R



1/1




1/1

FK172R



1/1



1/1

FK173R



ORTHOPEDICS/SPINE RETRACTORS

LANGE-HOHMANN Bone Retractor

Catalog Number	Blade Size (mm)	Length (in.)
FK160R	24	9 3/4



FK160R




ORTHOPEDICS/SPINE RETRACTORS

LANGE-HOHMANN Bone Retractor (continued)

Catalog Number	Blade Size (mm)	Length (in.)
FK161R	30	10 1/2
FK162R	33	11




FK161R-FK162R



FK161R

FK162R



FK161R
30 mm


FK162R
33 mm

ORTHOPEDICS/SPINE RETRACTORS





HOHMANN (WAGNER) Humoral Head Retractor

Catalog Number	Blade Size (mm)	Length (in.)
FK167R	17	8 3/4



HOHMANN (ALDINGER) Bone Retractor

Catalog Number	Blade Size (mm)	Length (in.)
FK169R	24	10 5/8




ORTHOPEDICS/SPINE RETRACTORS



HOHMANN Bone Retractor

Catalog Number	Blade Size (mm)	Length (in.)
FK185R	8	8 3/4
FK186R	15	9 3/4
FK187R	22	11 3/4





ORTHOPEDICS/SPINE RETRACTORS

HOHMANN Bone Retractor (continued)

Catalog Number	Blade Size (mm)	Length (in.)
FK188R	43	9




FK188R




ORTHOPEDICS/SPINE RETRACTORS

HOHMANN Bone Retractor (continued)

Catalog Number	Blade Size (mm)	Length (in.)
FK189R	40	11 1/2



FK189R



FK189R



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

ORTHOPEDICS/SPINE RETRACTORS

SCHUMACHER Modified WAGNER Bone Retractor

Catalog Number	Blade Size (mm)	Length (in.)
FK168R	37	14



FK168R


FK168R

SCHUMACHER Bone Retractor


Catalog Number	Blade Size (mm)	Length (in.)
FK177R	35	11 1/2



FK177R



FK177R




ORTHOPEDICS/SPINE RETRACTORS

VERBRUGGE-MÜLLER Bone Retractor

Catalog Number	Blade Size (mm)	Length (in.)
FK174R	44	9 1/2
FK175R	65	10



FK174R-FK175R




ORTHOPEDICS/SPINE RETRACTORS



FK180R, MD563

BENNETT Bone Retractor

Catalog Number	Blade Size (mm)	Length (in.)
FK180R	65	9 1/2
MD563	44	9 1/2




FK193R

MURPHY-LANE Bone Skid

Catalog Number	Blade Size (mm)	Length (in.)
FK193R	18 / 30	11 3/4



FK193R




ORTHOPEDICS/SPINE RETRACTORS

Bone Hook


Catalog Number	Tip	Blade End	Length (in.)
BT232R	SHARP	1 PRONG	9 5/8



BT232R



BT234R




KOCHER Bone Hook


Catalog Number	Tip	Blade End	Length (in.)
BT234R	SHARP	1 PRONG	8
BT238R	SHARP	1 PRONG	8 3/4
BT239R	BLUNT	1 PRONG	8 3/4



BT238R




BT239R




ORTHOPEDICS/SPINE RETRACTORS



BT240R




BT241R




VOLKMANN Bone Hook


Catalog Number	Tip	Blade End	Length (in.)
BT240R	SHARP	1 PRONG	8 3/4
BT241R	BLUNT	1 PRONG	8 3/4



MB418R, MB974R




MB975R



Bone Hooks

Catalog Number	Tip	Diameter (mm)	Length (in.)
MB974R	SHARP	19	9
MB975R	BLUNT	19	9
MB418R	SHARP	30	9





ORTHOPEDICS/SPINE RETRACTORS


ANDERSON-ADSON Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV250R	SHARP	4 x 4 PRONGS	7 1/2




NORFOLK-NORWICH Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV248R	BLUNT	4 x 5 PRONGS	8 3/4



ADSON Self Retaining Retractor


Catalog Number	Tip	Blade End	Length (in.)
BV252R	SHARP	3 x 4 PRONGS	10 1/2




ORTHOPEDICS/SPINE RETRACTORS

Laminectomy Self Retaining Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV254R	SHARP	6 x 6 PRONGS	11



BV254R




BV254R

ADSON Self Retaining Retractor, Hinged


Catalog Number	Tip	Blade End	Length (in.)
BV251R	SEMI-S	4 x 4 PRONGS	8 1/2



BV251R




BV251R




ORTHOPEDICS/SPINE RETRACTORS

HARVEY-JACKSON Retractor, Hinged

Catalog Number	Tip	Blade End	Length (in.)
BV258R	SHARP	3 x 3 PRONGS	11 3/4



BV258R




BV258R

BECKMANN-ADSON Retractor, Hinged

Catalog Number	Tip	Blade End	Length (in.)
BV271R	SEMI-S	4 x 4 PRONGS	12 1/4
BV272R	SHARP	4 x 4 PRONGS	12 1/4



BV271R-BV272R




BV271R

BV272R

ORTHOPEDICS/SPINE RETRACTORS



ADSON Retractor, Hinged

Catalog Number	Tip	Blade End	Length (in.)
BV273R	SEMI-S	5 x 4 PRONGS	12 1/2



BECKMANN Retractor, Hinged




Catalog Number	Tip	Blade End	Length (in.)
BV276R	SHARP	7 x 7 PRONGS	12 1/2



ORTHOPEDICS/SPINE RETRACTORS

BECKMANN Retractor, Hinged (continued)

Catalog Number	Tip	Blade End	Length (in.)
MG064R	SHARP	3 x 4 PRONGS	6 1/2
MG371R	BLUNT	3 x 4 PRONGS	6 1/4



ORTHOPEDICS/SPINE RETRACTORS

MISKIMON Cerebellum Retractor


Catalog Number	Tip	Blade End	Length (in.)
MF061R	SEMI-S	7 x 7 PRONGS	7 1/2

MEYERDING Retractor


Catalog Number	Blade Size (mm)	Length (in.)
MF065R	25 x 57	6 1/2
MF066R	25 x 51	6 1/2
MF067R	25 x 31	6 1/2

WILLIAMS Retractor

Catalog Number	Tip	Blade Size (mm)	Length (in.)
MF068R	BLADE LEFT	10 x 50	7
MF069R	BLADE RIGHT	10 x 50	7
MF070R	BLADE LEFT	10 x 60	7
MF071R	BLADE RIGHT	10 x 60	7
MF072R	BLADE LEFT	10 x 70	7
MF073R	BLADE RIGHT	10 x 70	7
MF074R	BLADE LEFT	20 x 50	7
MF075R	BLADE RIGHT	20 x 50	7
MF076R	BLADE LEFT	20 x 60	7
MF077R	BLADE RIGHT	20 x 60	7
MF078R	BLADE LEFT	20 x 70	7
MF079R	BLADE RIGHT	20 x 70	7




MF061R




MF061R




MF065R-MF067R




MF065R-MF067R



MF068R-MF079R





MF068R-MF079R



ORTHOPEDICS/SPINE RETRACTORS



DUCKER Retractor

Catalog Number	Tip	Blade Size (mm)	Length (in.)
MF329R	HOOK LEFT	50 x 20	6 3/4
MF330R	HOOK RIGHT	50 x 20	6 3/4
MF331R	HOOK LEFT	70 x 20	6 3/4
MF332R	HOOK RIGHT	70 x 20	6 3/4



MARKHAM-MEYERDING Retractor


Catalog Number	Tip	Blade Size (mm)	Length (in.)
MF442R	BLADE RIGHT	46 x 26	7 1/2
MF443R	BLADE LEFT	46 x 26	7 1/2








ORTHOPEDICS/SPINE RETRACTORS

Zelpi-Type Retractor, Deep Blade

Catalog Number	Tip	Depth (in.)
MD529R	ANG 45°	3
MD527R	ANG 90°	3
MD528R	ANG 90°	4
MD530R	ANG 90°	5



MD529R-MD530R



ORTHOPEDICS/SPINE RETRACTORS AND SPREADERS



BV296R

CLOWARD Lamina Spreader

Catalog Number	Tip	Tip Size (mm)	Length (in.)
BV296R	SERR OUTER	5.0	6 1/8



BV299R

CLOWARD Lamina Spreader

Catalog Number	Tip	Tip Size (mm)	Length (in.)
BV299R	SERR	4.5	5 1/4



MB989R

INGE Lamina Spreader

Catalog Number	Tip	Length (in.)
MB989R	TOOTHED	6 1/2



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



ORTHOPEDICS/SPINE RETRACTORS AND SPREADERS

INGE Lamina Spreader

Catalog Number	Tip	Tip Size (mm)	Length (in.)
BV297R	SERR	6.0	6 3/4




BV297R

INGE Lamina Spreader

Catalog Number	Tip	Tip Size (mm)	Length (in.)
BV298R	TOOTHED	12.0	10 3/4



BV298R




ORTHOPEDICS/SPINE RETRACTORS



Aesculap Ball Snap System on Spreaders and Blades

Advantages:

- Blades are safely locked without additional fixation, preventing them from detaching when the instrument is passed to the surgeon. In addition, the smooth surface minimizes the risk of catching a suture.
- Blades can be quickly exchanged during the operation.
- Different retractor patterns can be combined and exchanged at random, which allows for optimum selection of the blade size according to surgical needs.



BV510R

SCOVILLE-RICHTER Laminectomy Retractor, Reversible

Catalog Number	Description
BV510R	COMPLETE

BV510R Includes:

- 1 Retractor BV509R and
 - 1 Pair of blades BV516R
- The retractor bar and shanks rotate 180°



BV516R



BV517R



BV518R



BV519R

SCOVILLE-RICHTER Retractor Blade

Catalog Number	Blade Size (mm)
BV516R	62 x 52 AND 62 x 37, PAIR
BV517R	47 x 51
BV518R	30 x 40
BV519R	37 x 62

ORTHOPEDICS/SPINE RETRACTORS

- Complete retraction system consisting of a retractor body, various blade configurations, and storage tray. Designed to enable less invasive spine procedures such as discectomies and decompressions.
 - Unique organization system enables easy storage and retrieval of desired blade.
 - Sturdy design provides rigid retraction even at large opening widths.
- All retractor components are titanium. Blade designs improve visualization and intra-operative imaging.

For more information see our brochure number **DOC658**.

MF393 Includes:

- 1 Hinged Titanium Retractor Frame MF365T
- 2 Wide Blades, 40 mm MF368T
- 2 Wide Blades, 50 mm MF369T
- 2 Wide Blades, 60 mm MF370T
- 2 Wide Blades, 70 mm MF371T
- 2 Narrow Blades, 40 mm MF375T
- 2 Narrow Blades, 50 mm MF376T
- 2 Narrow Blades, 60 mm MF377T
- 2 Narrow Blades, 70 mm MF378T
- 1 Storage Tray MF392
- 1 Hook Blade, 30 mm MF381T
- 1 Hook Blade, 40 mm MF382T
- 1 Hook Blade, 50 mm MF383T
- 1 Hook Blade, 60 mm MF384T
- 1 Toothed Blade, 30 mm MF387T
- 1 Toothed Blade, 40 mm MF388T
- 1 Toothed Blade, 50 mm MF389T
- 1 Toothed Blade, 60 mm MF390T

Spine Retractor System


Catalog Number
MF393

Storage Tray for Spine Retractor System

Catalog Number
MF392


The storage tray will fit into Aesculap's full size container for sterilization.

For more information about Aesculap's SterilContainter System see our brochure number **DOC132**.




MF392

1/4




ORTHOPEDICS/SPINE RETRACTORS



Titanium Retractor Frame, Hinged

Catalog Number
MF365T



Titanium Retractor Frame, Rigid

Catalog Number
MF366T



ORTHOPEDICS/SPINE RETRACTORS



MF367T-MF373T

Titanium Retractor Blade, Wide

Catalog Number	Blade Size (mm)
MF367T	27 x 30
MF368T	27 x 40
MF369T	27 x 50
MF370T	27 x 60
MF371T	27 x 70
MF372T	27 x 80
MF373T	27 x 90



MF374T-MF380T

Titanium Retractor Blade, Narrow

Catalog Number	Blade Size (mm)
MF374T	20 x 30
MF375T	20 x 40
MF376T	20 x 50
MF377T	20 x 60
MF378T	20 x 70
MF379T	20 x 80
MF380T	20 x 90



MF381T-MF386T

Titanium Retractor Hook

Catalog Number	Length (mm)
MF381T	30
MF382T	40
MF383T	50
MF384T	60
MF385T	70
MF386T	80



MF387T-MF391T

Titanium Toothed Blades

Catalog Number	Length (mm)
MF387T	30
MF388T	40
MF389T	50
MF390T	60
MF391T	70




ORTHOPEDICS/SPINE RETRACTORS

CHARNLEY Initial Incision Retractor

Catalog Number	Description	Blade	Size	Blade Size (mm)
MB191R	SET COMPLETE	-	-	-
MB156R	INDIVIDUAL BLADE	SMALL	STANDARD DEPTH	25 x 38
MB177R	INDIVIDUAL BLADE	SMALL	MEDIUM DEPTH	25 x 75
MB192R	INDIVIDUAL BLADE	SMALL	LONG DEPTH	25 x 127
MB157R	INDIVIDUAL BLADE	LARGE	STANDARD DEPTH	25 x 44
MB178R	INDIVIDUAL BLADE	LARGE	MEDIUM DEPTH	25 x 75
MB193R	INDIVIDUAL BLADE	LARGE	LONG DEPTH	25 x 127
MG081R	WEIGHT AND CHAIN REPLACEMENT FOR MB191R	-	-	-

MB191R includes:

- 1 frame
- 1 small retractor blade (MB156R)
- 1 large retractor blade (MB157R)
- 1 weight w/ chain (3 lbs.)



MB191R-MG081R

ORTHOPEDICS/SPINE RASPS



FK507R

Bone File/Rasp

Catalog Number	Tip Size (mm)	Length (in.)
FK507R	20	8 3/4

One Side File Cut
One Side Rasp Cut



FK507R

Bone File

Catalog Number	Tip Size (mm)	Length (in.)
FK511R	13	9 1/2
FK510R	20	11 5/8




FK511R



FK511R



FK510R





ORTHOPEDICS/SPINE RASPS



DUROGRIP™ TC FOMON Rasp, Double Ended

Catalog Number	Size	Teeth	Length (in.)
OL423R	3 + 4	16/20	8 1/4
OL425R	5 + 6	24/30	8 1/4
OL427R	7 + 8	36/44	8 1/4
OL429R	9 + 10	52/44	8 1/4




FOMON Bone Rasp, Double Ended

Catalog Number	Tip Size (mm)	Tip	Length (in.)
OL405R	6.5 / 7.5	X-SERR Fine/Coarse	8 1/2
OL406R	6.5 / 7.5	Fine/Coarse	8 1/2
MG584R	8/8	FLAT, CONVEX FINE, COARSE	8 1/2
MA418R	7.5 / 7.5	Coarse/ Xtra-Coarse	8 3/8

ORTHOPEDICS/SPINE RASPS



OL407R



OL407R



MA419R



MA420R



OL408R-OL409R



OL408R




OL409R



FK301R-FK302R



1/1




1/1

FK301R



1/1



1/1

FK302R



FK505R



FK505R

MALTZ Bone Rasp

Catalog Number	Tip Size (mm)	Tip	Length (in.)
OL407R	9	DOWN-STROKE	6 1/2
MA419R	33	DOWN-STROKE	6 1/2
MA420R	33	UP STROKE	6 1/2

LEWIS Bone Rasp

Catalog Number	Tip	Tip Size (mm)	Length (in.)
OL408R	FINE-SERR DOWN-STROKE	8	6 7/8
OL409R	COARSE-SERR DOWN-STROKE	8	6 7/8

Bone Raspatory

Catalog Number	Blade Size (mm)	Length (in.)
FK301R	6	6 1/2
FK302R	8	6 1/2

PUTTI Bone File, Double Ended

Catalog Number	Tip Size (mm)	Length (in.)
FK505R	18/20	10 5/8

ORTHO
SPINE

INCHES

ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS



MB816R-MB817R



1/1



1/1



1/1



1/1

MB816R-MB817R

SAYRE Periosteal Elevator

Catalog Number	Tip	Length (in.)
MB816R	SHARP	6 3/4
MB817R	BLUNT	6 3/4



FK300R

WILLIGER Bone Elevator

Catalog Number	Blade Size (mm)	Length (in.)
FK300R	6	6 1/4



1/1



1/1

FK300R



FK306R-FK309R

ADSON Elevator

Catalog Number	Tip	Blade Size (mm)	Blade	Length (in.)
FK306R	SHARP	8	STR	6 3/4
FK308R	SEMI-S	7	STR	6 3/4
FK307R	BLUNT	7	CVD	6 3/4
FK309R	SHARP	7	CVD	6 3/4



1/1



1/1



1/1


FK306R



1/1



1/1



1/1

FK307R



1/1



1/1



1/1

FK308R



1/1



1/1





1/1

FK309R

ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS

MCGLAMRY Metatarsal Elevators


Catalog Number	Blade Size (mm)	Length (in.)
MB978R	11	6 1/2
MB979R	13	6 1/2
MB980R	15	6 1/2
MB981R	17	6 1/2



ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS



FK323R




Elevator W/ PEEK Handle

Catalog Number	Tip	Blade Size (mm)	Blade	Length (in.)
FK323R	SHARP	13	CVD	6 7/8



FK324R




Elevator W/ PEEK Handle

Catalog Number	Tip	Blade Size (mm)	Blade	Length (in.)
FK324R	SHARP	14	STR	7 1/2



FK325R



Elevator W/ PEEK Handle

Catalog Number	Tip	Blade Size (mm)	Blade	Length (in.)
FK325R	SHARP	14	CVD	7

ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS


LANGENBECK Periosteal Elevator

Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK360R	16	STR	7 1/2
MB279R	17	STR	7 3/4



FK360R

16 mm




FK360R



MB279R

17 mm




MB279R



FK350R

18 mm




FK350R



FK351R

15 mm



FK351R



MB536R




MB536R

SEDILLOT Elevator

Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK350R	18	STR	7 1/4
FK351R	15	CVD	8 1/2

KEY Periosteal Elevator

Catalog Number	Blade Size (mm)	Length (in.)
MB536R	3.1	6 3/4



ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS




KEY Elevator

Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK410R	6.4	STR	7
FK411R	9.5	STR	7 1/2
FK412R	12.7	STR	7 1/2
FK413R	19	STR	7 7/8
FK414R	25.4	STR	8 1/4



COBB Elevator, Round

Catalog Number	Blade Size (mm)	Blade	Length (in.)
MG024R	19	CVD	9 1/2



ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS

COBB Elevator, Sharp

Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK146R	10	CVD	11
FK147R	13	CVD	11
FK148R	19	CVD	11



FK146R-FK148R



10 mm

1/1



1/1

FK146R



13 mm

1/1



1/1

FK147R



19 mm

1/1



1/1

FK148R

COBB Periosteal Elevator, Round

Catalog Number	Blade Size (mm)	Blade	Length (in.)
MB518R	7	CVD	12



MB518R




7 mm

1/1



1/1

MB518R




ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS

COBB Elevator


Catalog Number	Blade Size (mm)	Blade	Length (in.)
MD645	10	ROUND	11
MD646	13	ROUND	11
MD647	19	ROUND	11
MD648	26	ROUND	11
MD649	32	ROUND	11



32 mm




26 mm



Elevator W/ PEEK Handle

Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK348R	20	CVD	11 1/2





FK348R

ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS

WAGNER Elevator W/ PEEK Handle

Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK352R	8	CVD	13
FK353R	17	CVD	13




ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS

Black Handled Elevator

Catalog Number	Tip Size (mm)	Handle	Length (in.)
MN001R	10	8 1/2	13 1/4
MN002R	15	8 1/2	13 1/4
MN003R	20	8 1/2	13 1/4
MN004R	25	8 1/2	13 1/4
MN005R	30	8 1/2	13 1/4
MN006R	35	8 1/2	13 1/4
MN007R	9	9	15
MN008R	13	9	15
MN009R	19	9	15
MN010R	25	9	15
MN011R	32	9	15




MN001R-MN006R




MN001R




MN002R




MN003R



MN005R



MN004R



MN006R




MN007R-MN011R




MN007R




MN008R



MN009R



MN010R



MN011R

ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS




WILLIGER Elevator

Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK120R	5	CVD	6 1/2
FK121R	5.5	CVD	6 3/4

SEBILEAU Elevator

Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK124R	5	CVD	6 7/8





ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS



LANGENBECK Elevator

Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK130R	8	CVD	7 3/4
FK131R	10	CVD	7 3/4
FK139R	11	CVD	9




QUERVAIN Elevator, Joker

Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK133R	6	CVD	7 3/4



LANE TRETHOWAN Bone Lever

Catalog Number	Blade	Length (in.)
FK135R	CVD	8 1/2



LANE Bone Lever

Catalog Number	Blade	Length (in.)
FK136R	CVD	10

ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS




CREGO Elevator

Catalog Number	Tip	Blade Size (mm)	Blade	Length (in.)
MB818R	BLUNT	6.5	CVD	8
MB819R	BLUNT	7.1	CVD	8
MB820R	BLUNT	14.6	CVD	8

BRISTOW Bone Lever

Catalog Number	Length (in.)
FK140R	5 3/4




ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS



TAEGER Bone Lever F/ Angle Joint
Operations


Catalog Number	Blade Size (mm)	Blade	Length (in.)
FK158R	12	CVD	6 7/8



ORTHOPEDICS/SPINE ELEVATORS AND BONE LEVERS

CHANDLER Bone Elevator

Catalog Number	Blade Size (mm)	Blade	Length (in.)
MB275R	14	SMALL	7 1/4
MB276R	19	MEDIUM	8
MB277R	25	LARGE	9 1/4
MB278R	30	EXTRA LARGE	10



ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES



HIBBS Chisel

Catalog Number	Tip Size (mm)	Blade	Length (in.)
MB850R	7	STR	9 1/2
MB991R	10	STR	9 1/2
MB851R	15	STR	9 1/2
MB852R	19	STR	9 1/2
MB853R	22	STR	9 1/2
MB854R	25	STR	9 1/2
MB855R	38	STR	9 1/2
MB856R	7	CVD	9 1/2
MB857R	9	CVD	9 1/2
MB992R	13	CVD	9 1/2
MB858R	15	CVD	9 1/2
MB859R	19	CVD	9 1/2
MB993R	22	CVD	9 1/2
MB994R	25	CVD	9 1/2
MB995R	32	CVD	9 1/2
MB996R	38	CVD	9 1/2




STILLE Chisel

Catalog Number	Tip Size (mm)	Length (in.)
FL503R	10	8
FL504R	12	8
FL505R	15	8
FL506R	20	8
FL507R	25	8



Chisel Handle F/ Interchangeable Blades


Catalog Number	Item
FL665R	HANDLE W/KEY
GB020R	REPLACEMENT KEY



ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES




FL666R-FL672R



FL666R-FL672R




FL682R




FL682R




FL685R-FL687R



FL685R



FL685R



FL687R

Chisel Blade, F/FL665R


Catalog Number	Tip Size (mm)	Tip Size (in.)	Working Length (mm)
FL666R	5	3/16	45
FL668R	10	3/8	45
FL670R	16	5/8	45
FL672R	25	1	45

Osteotome Blade, F/FL665R

Catalog Number	Tip Size (mm)	Tip Size (in.)	Working Length (mm)
FL682R	25	1	78

Concave Chisel Blade, F/FL665R


Catalog Number	Radius (mm)	Tip Size (mm)	Tip Size (in.)
FL685R	50	25	1
FL687R	70	25	1



ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES

Lamina Gouge W/ HARTPRESS Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
FL437R	6	1/4	9 1/4
FL438R	8	5/16	9 1/4
FL439R	11	7/16	9 1/4
FL440R	16	5/8	9 1/4



Chisel W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
FL430R	20	3/4	9 5/8




LEXER-MINI Gouge W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
FL390R	4	5/32	7
FL391R	6	1/4	7
FL392R	8	5/16	7
FL393R	10	3/8	7
FL394R	12	1/2	7




ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES



LEXER Gouge W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
FL411R	5	3/16	8 7/8
FL412R	7	9/32	8 7/8
FL413R	10	3/8	8 7/8
FL415R	15	5/8	8 7/8
FL416R	20	3/4	8 7/8
FL417R	25	1	8 7/8
FL418R	30	1 3/16	8 7/8



Gouge, F/ Removing Broken Screws, W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
FL420R	9	3/8	8 7/8



LEXER Spongiosa Gouge W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
FL423R	10	3/8	11

ORTHO
SPINE

INCHES

ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES



ZIELKE Gouge W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
FL701R	5.5	7/32	9
FL702R	9	3/8	9




FL701R-FL702R

Spongiosa Gouge W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Blade	Length (in.)
FL441R	5	3/16	CVD	9 3/4
FL443R	10	3/8	CVD	9 3/4
FL445R	15	5/8	CVD	9 3/4



FL441R-FL445R




FL441R-FL445R



FL573R-FL577R

STILLE Gouge

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
FL573R	10	3/8	8
FL574R	12	1/2	8
FL575R	15	5/8	8
FL576R	20	3/4	8
FL577R	25	1	8



FL573R-FL577R


ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES

WAGNER Gouge W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Blade Shape	Length (in.)
ND120R	9	3/8	ANG FWD	13 3/8
ND121R	9	3/8	ANG BACK	13 3/8
ND122R	15	5/8	ANG FWD	13 3/8
ND123R	15	5/8	ANG BACK	13 3/8
ND124R	9	3/8	STR	13 3/8
ND125R	15	5/8	STR	13 3/8




ND120R-ND125R




9 mm

ND120R, ND121R, ND124R




ND121R
ND123R



15 mm

ND122R, ND123R, ND125R



ND124R
ND125R



MF465R-MF478R


SMITH-PETERSON Gouge

Catalog Number	Tip	Tip Size (mm)	Length (in.)
MF465R	STR	6	8
MF466R	STR	10	8
MF467R	STR	13	8
MF468R	STR	16	8
MF469R	STR	19	8
MF470R	STR	25	8
MF471R	STR	32	8
MF472R	CVD	6	8
MF473R	CVD	10	8
MF474R	CVD	13	8
MF475R	CVD	16	8
MF476R	CVD	19	8
MF477R	CVD	25	8
MF478R	CVD	32	8

ORTHO
SPINE

INCHES

ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES




LEXER-MINI Osteotome W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
FL380R	4	5/32	7
FL381R	6	1/4	7
FL382R	8	5/16	7
FL383R	10	3/8	7
FL384R	12	1/2	7



LEXER Osteotome W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
FL401R	5	3/16	8 7/8
FL402R	7	9/32	8 7/8
FL403R	10	3/8	8 7/8
FL405R	15	5/8	8 7/8
FL406R	20	3/4	8 7/8
FL407R	25	1	8 7/8
FL408R	30	1 3/16	8 7/8
ND108R	10	3/8	10 3/4
ND109R	15	5/8	10 3/4



Chisel W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
ND105R	16	5/8	10 1/2





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES

COBB Gouge

Catalog Number	Blade	Length (in.)
MB707R	STR	11 1/2



MB707R



MB707R

Gouge, Double Curved, W/ HARTPRESS Handle

Catalog Number	Tip Size (mm)	Tip Size (in.)	Length (in.)
ND118R	25	1	11 1/4



ND118R




MB709R-MB715R



MB716R-MB723R

HIBBS Gouge

Catalog Number	Tip Size (mm)	Tip Size (in.)	Blade	Length (in.)
MB709R	16	5/8	STR	9 1/2
MB710R	19	3/4	STR	9 1/2
MB711R	22	7/8	STR	9 1/2
MB712R	25	1	STR	9 1/2
MB713R	28	1 1/8	STR	9 1/2
MB714R	32	1 1/4	STR	9 1/2
MB715R	38	1 1/2	STR	9 1/2
MB716R	6	1/4	CVD	9 1/2
MB717R	9	3/8	CVD	9 1/2
MB718R	13	1/2	CVD	9 1/2
MB719R	13	5/8	CVD	9 1/2
MB720R	19	3/4	CVD	9 1/2
MB721R	25	1	CVD	9 1/2
MB722R	28	1 1/8	CVD	9 1/2
MB723R	38	1 1/2	CVD	9 1/2



ORTHO
SPINE

ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES



FL156R-FL170R



FI156R-FL160R



FL166R-FL170R

HOKE Osteotome

Catalog Number	Tip Size (mm)	Tip Size (in.)	Blade	Length (in.)
FL156R	3	1/8	STR	6 1/8
FL157R	6	1/4	STR	6 1/8
FL158R	10	3/8	STR	6 1/8
FL159R	13	1/2	STR	6 1/8
FL160R	16	5/8	STR	6 1/8
FL166R	3	1/8	CVD	6 1/8
FL167R	6	1/4	CVD	6 1/8
FL168R	10	3/8	CVD	6 1/8
FL169R	13	1/2	CVD	6 1/8
FL170R	16	5/8	CVD	6 1/8



FL512R-FL517R



FL512R-FL517R



FL523R-FL527R



FL523R-FL527R



FL553R-FL557R




FL553R-FL557R

STILLE Osteotome

Catalog Number	Tip Size (mm)	Tip Size (in.)	Blade	Length (in.)
FL512R	8	5/16	STR, HVY	8
FL513R	10	3/8	STR, HVY	8
FL514R	12	1/2	STR, HVY	8
FL515R	15	5/8	STR, HVY	8
FL516R	20	3/4	STR, HVY	8
FL517R	25	1	STR, HVY	8
FL523R	10	3/8	STR	8
FL524R	12	1/2	STR	8
FL525R	15	5/8	STR	8
FL526R	20	3/4	STR	8
FL527R	25	1	STR	8
FL553R	10	3/8	CVD	8
FL554R	12	1/2	CVD	8
FL555R	15	5/8	CVD	8
FL556R	20	3/4	CVD	8
FL557R	25	1	CVD	8

ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES




Gouge

Catalog Number	Tip Size (mm)	Tip	Length (in.)
MF458R	4	CVD	7
MF459R	4	STR	7
MF460R	6	CVD	7 1/2
MF461R	6	STR	7 1/2



SMITH-PETERSON Osteotome

Catalog Number	Tip Size (mm)	Tip Size (in.)	Blade	Length (in.)
MB657R	7	9/32	STR	8
MB658R	11	7/16	STR	8
MB659R	14	9/16	STR	8
MB660R	17	21/32	STR	8
MB661R	19	3/4	STR	8
MB662R	25	1	STR	8
MB663R	6	1/4	CVD	8
MB664R	9.5	3/8	CVD	8
MB665R	13	1/2	CVD	8
MB666R	16	5/8	CVD	8
MB667R	19	3/4	CVD	8
MB668R	25	1	CVD	8
MB669R	32	1 1/4	CVD	8



HIBBS Osteotome

Catalog Number	Tip Size (mm)	Tip Size (in.)	Blade	Length (in.)
MB675R	6	1/4	STR	9 5/8
MB676R	10	3/8	STR	9 5/8
MB677R	13	1/2	STR	9 5/8
MB678R	16	5/8	STR	9 5/8
MB679R	19	3/4	STR	9 5/8
MB680R	22	7/8	STR	9 5/8
MB681R	25	1	STR	9 5/8
MB682R	28	1 1/8	STR	9 5/8
MB683R	32	1 1/4	STR	9 5/8
MB684R	38	1 1/2	STR	9 5/8
MB685R	6	1/4	CVD	9 5/8
MB686R	10	3/8	CVD	9 5/8
MB687R	13	1/2	CVD	9 5/8
MB688R	16	5/8	CVD	9 5/8
MB689R	19	3/4	CVD	9 5/8
MB690R	22	7/8	CVD	9 5/8
MB691R	25	1	CVD	9 5/8
MB692R	28	1 1/8	CVD	9 5/8
MB693R	32	1 1/4	CVD	9 5/8
MB694R	38	1 1/2	CVD	9 5/8



ORTHO
SPINE

INCHES

ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES

COBB Osteotomes

Catalog Number	Tip Size (mm)	Blade	Length (in.)
MB084R	6	STR	11
MB085R	10	STR	11
MB742R	13	STR	11
MB087R	13	STR	11
MB088R	16	STR	11
MB089R	19	STR	11
MB090R	25	STR	11
MG030R	6	CVD	11
MG031R	10	CVD	11
MG032R	13	CVD	11
MG033R	16	CVD	11
MG036R	19	CVD	11
MG034R	25	CVD	11
MG035R	32	CVD	11



ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES



MB404R-MB373R



MB404R-MB373R



MB398R-MB403R



MB398R-MB403R



MB388R-MB389R

Osteotome W/ Hexagonal Handle

Catalog Number	Tip Size (mm)	Length (in.)
MB404R	2	6 1/4
MB405R	3	6 1/4
MB406R	4	6 1/4
MB407R	5	6 1/4
MB408R	6	6 1/4
MB409R	7	6 1/4
MB410R	8	6 1/4
MB411R	9	6 1/4
MB412R	10	6 1/4
MB413R	11	6 1/4
MB373R	12	6 1/4

CONVERSE Osteotome

Catalog Number	Tip Size (mm)	Length (in.)
MB398R	2	7
MB399R	4	7
MB400R	6	7
MB401R	8	7
MB402R	10	7
MB403R	12	7

RUBIN Osteotome, W/ Finger Grip

Catalog Number	Tip Size (mm)	Length (in.)
MB388R	10	6 1/2
MB389R	14	6 1/2

ORTHO
SPINE

INCHES

ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES

FL088R-MG049R 1/2

FL088R-MG049R

FL088R-MG049R 1/1

FL088R-MG049R

MG037R-MG044R 1/2

MG037R-MG044R

MG037R-MG044R 1/1

MG037R-MG044R

MB621R-MB628R 1/2

MB621R-MB628R

MB620R-MB628R 1/1

MB620R-MB628R

LAMBOTTE Osteotome

Catalog Number	Tip Size (mm)	Blade	Length (in.)
MG050R	SET OF 7	STR	-
FL088R	2	STR	5
FL101R	4	STR	5
FL102R	6	STR	5
FL104R	8	STR	5
FL106R	10	STR	5
FL108R	12	STR	5
MG046R	14	STR	5
FL109R	15	STR	5
MG047R	18	STR	5
MG048R	19	STR	5
MG049R	20	STR	5
MG037R	2	CVD	5
MG038R	4	CVD	5
MG039R	6	CVD	5
MG040R	8	CVD	5
MG041R	10	CVD	5
MG042R	12	CVD	5
MG043R	14	CVD	5
MG044R	15	CVD	5

LAMBOTTE MINI OSTEOTOME SET:

CONTAINS SEVEN MG050R STRAIGHT 5" OSTEOTOMES:

2MM, 4MM, 6MM, 8MM, 10MM, 12MM, 14MM PLUS MALLET (FL036R)

LAMBOTTE Osteotome with Calibrations

Catalog Number	Tip Size (mm)	Blade	Length (in.)
MB620R	COMPLETE SET	-	-
MB621R	4	STR	7 1/4
MB622R	6	STR	7 1/4
MB623R	7	STR	7 1/4
MB624R	8	STR	7 1/4
MB625R	10	STR	7 1/4
MB626R	12	STR	7 1/4
MB627R	16	STR	7 1/4
MB628R	22	STR	7 1/4

MB620R includes:

Osteotome MB621R-MB628R and Mallet MB629R

ORTHOPEDICS/SPINE CHISELS, GOUGES AND OSTEOTOMES



MB640R-MB646R



MB640R-MB646R



MB650R-MB656R



MB650R-MB656R

LAMBOTTE Osteotome

Catalog Number	Tip Size (mm)	Tip Size (in.)	Blade	Length (in.)
MB640R	COMPLETE SET	STR		
MB641R	6 1/4 STR			9
MB642R	13 1/2 STR			9
MB643R	19 3/4 STR			9
MB644R	25 1 STR			9
MB645R	32 1/4 STR			9
MB646R	38 1/2 STR			9
MB650R	COMPLETE SET	CVD		
MB651R	6 1/4 CVD			9
MB652R	13 1/2 CVD			9
MB653R	19 3/4 CVD			9
MB654R	25 1 CVD			9
MB655R	32 1/4 CVD			9
MB656R	38 1/2 CVD			9

LAMBOTTE Osteotome

Catalog Number	Tip Size (mm)	Blade	Length (in.)
FL650R	4	STR	9 5/8
FL651R	8	STR	9 5/8
FL652R	10	STR	9 5/8
FL653R	13	STR	9 5/8
FL654R	15	STR	9 5/8
FL655R	18	STR	9 5/8
FL656R	20	STR	9 5/8
FL658R	25	STR	9 5/8
FL659R	30	STR	9 5/8
FL660R	38	STR	9 5/8
FL661R	44	STR	9 5/8
FL662R	50	STR	9 5/8
FL781R	8	CVD	9 5/8
FL783R	13	CVD	9 5/8
FL785R	18	CVD	9 5/8
FL788R	25	CVD	9 5/8
FL789R	30	CVD	9 5/8
FL790R	38	CVD	9 5/8



FL650R-FL662R



FL650R-FL662R



FL781R-FL790R



FL781R-FL790R



ORTHOPEDICS/SPINE MALLET

Mallet, PEEK

Catalog Number	Handle	Weight (gr)	Diameter (mm)	Length (in.)
FL010	HARTPRESS	230	42	9 3/4
FL011	HARTPRESS	315	50	9 3/4



FL010-FL011

MEAD Mallet, Removable Discs, 180gr


Catalog Number	Description	Diameter (mm)	Length (in.)
FL035R	COMPLETE	20	7 1/2
FL034	SPARE DISCS, F/FL035R, 2/SET	20	



FL035R

Mallet W/ Removable Nylon Discs, 135gr

Catalog Number	Description	Diameter (mm)	Length (in.)
FL036R	COMPLETE	25	7 1/4
FL037	NYLON SPARE DISCS, F/FL036R, 2/SET	25	



FL036R

Mallet, 226gr

Catalog Number	Description	Diameter (mm)	Length (in.)
MB629R	ONE SIDE NYLON, ONE SIDE METAL	26	7 1/4
MB629200	REPLACEMENT NYLON DISC FOR MB629R MALLET HEAD		



MB629R

ORTHOPEDICS/SPINE MALLET



DF114R

Mallet, 155gr

Catalog Number	Diameter (mm)	Length (in.)
DF114R	25	7 1/2



US142

Mallet, Copper


Catalog Number	Diameter (mm)	Length (in.)
US142	20	6 1/2



MG583R

Kirk Mallet, 560gr

Catalog Number	Diameter (mm)	Length (in.)
MG583R	35	7 1/4



ORTHOPEDICS/SPINE MALLET



FL038R

COTTLE Mallet, Flat/Round, 235gr

Catalog Number	Diameter (mm)	Length (in.)
FL038R	30	7 1/4



FL044R

HAJEK Mallet, 140gr

Catalog Number	Diameter (mm)	Length (in.)
FL044R	27	8 3/4



FL055R

BERGMANN Mallet, 260gr

Catalog Number	Diameter (mm)	Length (in.)
FL055R	30	9 1/4

PEEK

1/2

ORTHOPEDICS/SPINE MALLET

Mallet, 285gr, W/ Silicone Handle

Catalog Number	Diameter (mm)	Length (in.)
FL059R	30 / 42	9 1/2



BERGMANN Mallet, 300gr, W/ Silicone Handle

Catalog Number	Diameter (mm)	Length (in.)
FL060R	45	9 5/8




Orthopedic Mallet

Catalog Number	Diameter (mm)	Length (in.)
MB630R	30	7 1/4
MB631R	38	7 1/4

MB630R: 1 LB Weight (485gr)

MB631R: 2 LB Weight (907gr)




ORTHOPEDICS/SPINE MALLET

WAGNER Mallet Set, 130gr, W/ Silicone Handle

Catalog Number	Diameter (mm)	Length (in.)
FL062	35	9 1/4


Set includes:

- 1 PR Metal Discs FL064R
- 1 PR Plastic Discs FL065
- 1 Wrench FL062/206



OMBREDANNE Mallet, 520gr, W/ Silicone Handle

Catalog Number	Diameter (mm)	Length (in.)
FL066R	40	9 1/2



HEATH Copper Head Mallet, 725gr, W/ Silicone Handle

Catalog Number	Diameter (mm)	Length (in.)
FL068R	40	9 1/4





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

ORTHOPEDICS/SPINE MALLETS, IMPACTORS AND BONE TAMPS



MG095R

HEATH Mallet, 823gr

Catalog Number	Diameter (mm)	Size	Length (in.)
MG095R	40	1.8 lbs	7 1/4



FL090R

Bone Graft Tamper

Catalog Number	Diameter (mm)	Length (in.)
FL090R	12	7 7/8



FL090R



MB796R-MB800R

Bone Tamp

Catalog Number	Diameter (mm)	Length (in.)
MB796R	2	6 1/4
MB798R	4	6 1/4
MB800R	6	6 1/4
MB801R	8	6 1/4
MB802R	10	6 1/4



MB796R




MB798R



MB800R




MB801R



MB801R



MB802R



MB802R



FG093R-FG096R

SYPERT Spinal Impactor, W/ HART-PRESS Handle

Catalog Number	Tip Size (mm)	Length (in.)
FG093R	3 x 8	11 1/2
FG094R	4 x 10	11 1/2
FG095R	5 x 15	11 1/2
FG096R	6 x 20	11 1/2



FG093R



FG094R



FG095R



FG096R



ORTHOPEDICS/SPINE BONE IMPACTOR

COBB Bone Impactor


Catalog Number	Tip Size (mm)	Length (in.)
MB545R	3 x 12	11 1/4



MB545R

1/2

1/1
MB545R



ORTHO
SPINE

MILLIMETERS




AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C





ORTHOPEDICS/SPINE HAND DRILL



MG021R

RALKS Hand Drill, Cannulated

Catalog Number	Length (in.)
MG021R	6 1/4



ORTHOPEDICS/SPINE DRILLS, BURRS AND REAMERS



CHAMPIONNIÈRE Bone Drill

Catalog Number	Length (in.)
FR110R	6 3/4



Awl/Reamer

Catalog Number	Length (in.)
FR118R	6



Reamer

Catalog Number	Length (in.)
FR120R	7

ORTHOPEDICS/SPINE REAMERS AND TREPHINES



FR125R

PERTHES Reamer

Catalog Number	Length (in.)
FR125R	8 1/2



FR140R-FR141R

Bone Biopsy Trephine


Catalog Number	Diameter (mm)	Length (in.)
FR140R	3.9 / 5	6 1/4
FR141R	6.4 / 8	6 1/4



FR140R



FR141R



ORTHOPEDICS/SPINE CURETTES

Bone Curette, Double Ended

Catalog Number	Tip Size (mm)	Length (in.)
FK845R	6 / 7	5 1/8
FK846R	7 / 8	6 1/4
FK848R	8 / 11	7 7/8
FK850R	7 / 6	8 1/4

FK845R 1/2



FK845R 1/1

FK846R 1/2



FK846R 1/1

FK848R 1/2




FK848R 1/1

FK850R 1/2




FK850R 1/1

FK658R-FK659R 1/2



FK658R 1/1



FK659R 1/1

WILLIGER Bone Curette

Catalog Number	Tip Size (#)	Tip Size (mm)	Length (in.)
FK658R	00	3.4	6 3/4
FK659R	1	5.2	6 3/4

ORTHO
SPINE

MILLIMETERS

ORTHOPEDICS/SPINE CURETTES

BRUNS Curette, Round

Catalog Number	Tip Size (#)	Tip Size (mm)	Length (in.)
FK611R	000	3	6 1/2
FK612R	00	4	6 1/2
FK613R	0	5	6 1/2
FK614R	1	6	6 1/2
FK615R	2	8	6 1/2
FK616R	3	10	6 1/2
FK617R	4	12	6 1/2

SCHEDE Curette


Catalog Number	Tip Size (#)	Tip Size (mm)	Length (in.)
FK621R	000	2.5	6 3/4
FK622R	00	3.4	6 3/4
FK623R	0	4.2	6 3/4
FK624R	1	5.2	6 3/4
FK625R	2	6.3	6 3/4
FK626R	3	7.4	6 3/4
FK627R	4	8.0	6 3/4

VOLKMANN Curette


Catalog Number	Tip Size (#)	Tip Size (mm)	Length (in.)
FK630R	0000	2.8	6 3/4
FK631R	000	3.6	6 3/4
FK632R	00	4.4	6 3/4
FK633R	0	5.2	6 3/4
FK634R	1	6.8	6 3/4
FK635R	2	8.5	6 3/4
FK636R	3	10.0	6 3/4
FK637R	4	11.5	6 3/4
FK638R	5	13.0	6 3/4
FK639R	6	14.5	6 3/4



3 mm
FK611R




4 mm
FK612R




5 mm
FK613R




6 mm
FK614R




8 mm
FK615R




10 mm
FK616R




12 mm
FK617R




2.5 mm
FK621R




3.4 mm
FK622R




4.2 mm
FK623R




5.2 mm
FK624R




6.3 mm
FK625R




7.4 mm
FK626R




8 mm
FK627R




2.8 mm
FK630R




3.6 mm
FK631R




4.4 mm
FK632R




5.2 mm
FK633R




6.8 mm
FK634R




8.5 mm
FK635R




10 mm
FK636R



11.5 mm
FK637R



13 mm
FK638R



14.5 mm
FK639R







ORTHOPEDICS/SPINE CURETTES








DAUBENSPECK Bone Curette

Catalog Number	Tip Size (#)	Tip Size (mm)	Length (in.)
FK662R	0000	2.8	7 7/8
FK663R	000	3.6	7 7/8
FK664R	00	4.4	7 7/8
FK665R	0	5.2	7 7/8



LEMPERT Bone Curette

Catalog Number	Tip Size (mm)	Length (in.)
FK675R	1.8	8 1/2
FK676R	2.0	8 1/2
FK677R	2.4	8 1/2
FK678R	2.6	8 1/2
FK679R	2.8	8 1/2



EPSTEIN Spinal Curette


Catalog Number	Tip Size (mm)	Size (#)	Length (in.)
MD565	2.0	000	8
MD566	2.5	00	8
MD567	3.0	0	8



ORTHOPEDICS/SPINE CURETTES


SCOVILLE Bone Curette

Catalog Number	Tip Size (#)	Tip Size (mm)	Tip	Length (in.)
FK702R	000	3.6	FWD STR	9 3/4
FK703R	00	4.4	FWD STR	9 3/4
FK704R	0	5.2	FWD STR	9 3/4
FK705R	1	6.8	FWD STR	9 3/4
FK712R	000	3.6	FWD ANG	9 3/4
FK713R	00	4.4	FWD ANG	9 3/4
FK714R	0	5.2	FWD ANG	9 3/4
FK715R	1	6.8	FWD ANG	9 3/4



BUSHE Bone Curette

Catalog Number	Tip Size (mm)	Tip	Length (in.)
FK723R	3.0	REV ANG	9 3/4
FK724R	4.0	REV ANG	9 3/4
FK725R	5.0	REV ANG	9 3/4





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



ORTHOPEDICS/SPINE CURETTES




MF519R-MD490, MF534R-MF544R, MF567R-MF577R




MF521R-MD481, MF545R-MF555R, MF578R-MF588R




MF523R-MF533R, MF556R-MF566R, MF589R-MF599R




MD473




MD474




MD475



MD476




MD477



MD478



MD479



MD480



MD481

BRUNS Spinal Curette

Catalog Number	Tip Size (#)	Tip	Length (in.)
MF519R	00000	STR	8
MF520R	0000	STR	8
MD482	000	STR	8
MD483	00	STR	8
MD484	0	STR	8
MD485	1	STR	8
MD486	2	STR	8
MD487	3	STR	8
MD488	4	STR	8
MD489	5	STR	8
MD490	6	STR	8
MF521R	00000	FWD ANG	8
MF522R	0000	FWD ANG	8
MD473	000	FWD ANG	8
MD474	00	FWD ANG	8
MD475	0	FWD ANG	8
MD476	1	FWD ANG	8
MD477	2	FWD ANG	8
MD478	3	FWD ANG	8
MD479	4	FWD ANG	8
MD480	5	FWD ANG	8
MD481	6	FWD ANG	8
MF523R	00000	REV ANG	8
MF524R	0000	REV ANG	8
MF525R	000	REV ANG	8
MF526R	00	REV ANG	8
MF527R	0	REV ANG	8
MF528R	1	REV ANG	8
MF529R	2	REV ANG	8
MF530R	3	REV ANG	8
MF531R	4	REV ANG	8
MF532R	5	REV ANG	8
MF533R	6	REV ANG	8

Catalog Number	Tip Size (#)	Tip	Length (in.)
MF534R	00000	STR	9
MF535R	0000	STR	9
MF536R	000	STR	9
MF537R	00	STR	9
MF538R	0	STR	9
MF539R	1	STR	9
MF540R	2	STR	9
MF541R	3	STR	9
MF542R	4	STR	9
MF543R	5	STR	9
MF544R	6	STR	9
MF545R	00000	FWD ANG	9
MF546R	0000	FWD ANG	9
MF547R	000	FWD ANG	9
MF548R	00	FWD ANG	9
MF549R	0	FWD ANG	9
MF550R	1	FWD ANG	9
MF551R	2	FWD ANG	9
MF552R	3	FWD ANG	9
MF553R	4	FWD ANG	9
MF554R	5	FWD ANG	9
MF555R	6	FWD ANG	9
MF556R	00000	REV ANG	9
MF557R	0000	REV ANG	9
MF558R	000	REV ANG	9
MF559R	00	REV ANG	9
MF560R	0	REV ANG	9
MF561R	1	REV ANG	9
MF562R	2	REV ANG	9
MF563R	3	REV ANG	9
MF564R	4	REV ANG	9
MF565R	5	REV ANG	9
MF566R	6	REV ANG	9


Catalog Number	Tip Size (#)	Tip	Length (in.)
MF567R	00000	STR	12
MF568R	0000	STR	12
MF569R	000	STR	12
MF570R	00	STR	12
MF571R	0	STR	12
MF572R	1	STR	12
MF573R	2	STR	12
MF574R	3	STR	12
MF575R	4	STR	12
MF576R	5	STR	12
MF577R	6	STR	12
MF578R	00000	FWD ANG	12
MF579R	0000	FWD ANG	12
MF580R	000	FWD ANG	12
MF581R	00	FWD ANG	12
MF582R	0	FWD ANG	12
MF583R	1	FWD ANG	12
MF584R	2	FWD ANG	12
MF585R	3	FWD ANG	12
MF586R	4	FWD ANG	12
MF587R	5	FWD ANG	12
MF588R	6	FWD ANG	12
MF589R	00000	REV ANG	12
MF590R	0000	REV ANG	12
MF591R	000	REV ANG	12
MF592R	00	REV ANG	12
MF593R	0	REV ANG	12
MF594R	1	REV ANG	12
MF595R	2	REV ANG	12
MF596R	3	REV ANG	12
MF597R	4	REV ANG	12
MF598R	5	REV ANG	12
MF599R	6	REV ANG	12




ORTHOPEDICS/SPINE CURETTES




MB872R-MB880R




MB872R




MB873R




MB874R




MB875R




MB876R




MB877R



MB878R



MB879R




MB880R

BRUNS Curette


Catalog Number	Tip Size (#)	Tip Size (mm)	Tip	Length (in.)
MB872R	000	2.0	STR	9
MB873R	00	2.5	STR	9
MB874R	0	3.3	STR	9
MB875R	1	3.8	STR	9
MB876R	2	4.8	STR	9
MB877R	3	5.3	STR	9
MB878R	4	6.0	STR	9
MB879R	5	6.7	STR	9
MB880R	6	7.3	STR	9




FK688R-FK690R



FK688R



FK689R



FK690R

SIMON Bone Curette

Catalog Number	Tip Size (#)	Tip Size (mm)	Length (in.)
FK688R	000	3.6	9 1/4
FK689R	00	4.4	9 1/4
FK690R	0	5.2	9 1/4



MB415R-MB417R



MB415R




MB416R



MB417R

HIBBS Curette

Catalog Number	Tip	Length (in.)
MB415R	SMALL	8
MB416R	MEDIUM	8
MB417R	LARGE	8







ORTHOPEDICS/SPINE CURETTES

Bone Curette, Long, W/ PEEK Handle

Catalog Number	Tip Size (#)	Tip Size (mm)	Length (in.)
FK666R	00	4.4	10 5/8
FK668R	1	6.8	10 5/8
FK670R	2	8.5	10 5/8
FK672R	3	10.0	10 5/8



FK666R-FK672R






CASPAR Curette, W/ PEEK Handle




Catalog Number	Tip Size (#)	Tip Size (mm)	Blade	Length (in.)
FK773R	000	3.6	STR	8 3/4
FK774R	00	4.4	STR	8 3/4
FK775R	0	5.2	STR	8 3/4
FK783R	000	3.6	CVD	8 3/4
FK784R	00	4.4	CVD	8 3/4
FK785R	0	5.2	CVD	8 3/4



FK773R-FK775R



FK783R-FK785R







ORTHOPEDICS/SPINE CURETTES



CASPAR Curette, W/ PEEK Handle

Catalog Number	Tip Size (mm)	Tip	Length (in.)
FK834R	4	SQUARE	8 3/4
FK835R	5	SQUARE	8 3/4
FK836R	6	SQUARE	8 3/4




ORTHOPEDICS/SPINE CURETTES

CASPAR Bone Curette, Toothed

Catalog Number	Tip Size (mm)	Tip	Length (in.)
FK854R	4.0	SQUARE	9 3/4
FK855R	5.0	SQUARE	9 3/4



FK854R-FK855R



CASPAR Bone Curette, Toothed

Catalog Number	Tip Size (mm)	Tip	Length (in.)
FK863R	3.0	SQUARE, ANG	10 5/8
FK864R	4.0	SQUARE, ANG	10 5/8




FK863R-FK864R



ORTHOPEDICS/SPINE CURETTES

COBB Curette

Catalog Number	Tip Size (#)	Tip Size (mm)	Tip	Length (in.)
FK581R	000	2.0	STR	11
FK582R	00	2.5	STR	11
FK583R	0	3.0	STR	11
FK584R	1	3.5	STR	11
FK585R	2	4.5	STR	11
FK586R	3	5.5	STR	11
FK587R	4	7.5	STR	11
FK588R	5	8.5	STR	11
FK589R	6	10.0	STR	11
FK591R	000	2.0	ANG	11
FK592R	00	2.5	ANG	11
FK593R	0	3.0	ANG	11
FK594R	1	3.5	ANG	11
FK595R	2	4.5	ANG	11
FK596R	3	5.5	ANG	11
FK597R	4	7.5	ANG	11
FK598R	5	8.5	ANG	11
FK599R	6	10.0	ANG	11



$\frac{1}{1}$
○ 2.0 mm
FK581R, FK591R

$\frac{1}{1}$
○ 2.5 mm
FK582R, FK592R

$\frac{1}{1}$
○ 3.0 mm
FK583R, FK593R

$\frac{1}{1}$
○ 3.5 mm
FK584R, FK594R

$\frac{1}{1}$
○ 4.5 mm
FK585R, FK595R

$\frac{1}{1}$
○ 5.5 mm
FK586R, FK596R

$\frac{1}{1}$
○ 7.5 mm
FK587R, FK597R

$\frac{1}{1}$
○ 8.5 mm
FK588R, FK598R

$\frac{1}{1}$
○ 10.0 mm
FK589R, FK599R

ORTHO
SPINE

INCHES

ORTHOPEDICS/SPINE CURETTES

Colored Cervical Curette


Catalog Number	Tip Size (#)	Tip	Color	Length (in.)
MD216R	000000	STR	RED	10
MD217R	00000	STR	SILVER	10
MD218R	0000	STR	PURPLE	10
MD219R	000	STR	TURQUOISE	10
MD220R	00	STR	GOLD	10
MD221R	0	STR	GREEN	10
MD222R	1	STR	ORANGE	10
MD223R	2	STR	BLUE	10
MD224R	3	STR	GRAY	10
MD225R	4	STR	BLACK	10
MD226R	5	STR	BLACK	10
MD227R	6	STR	BLACK	10
MD234R	000000	ANG	RED	10
MD235R	00000	ANG	SILVER	10
MD236R	0000	ANG	PURPLE	10
MD237R	000	ANG	TURQUOISE	10
MD238R	00	ANG	GOLD	10
MD239R	0	ANG	GREEN	10
MD240R	1	ANG	ORANGE	10
MD241R	2	ANG	BLUE	10
MD242R	3	ANG	GRAY	10
MD243R	4	ANG	BLACK	10
MD244R	5	ANG	BLACK	10
MD245R	6	ANG	BLACK	10



ORTHOPEDICS/SPINE CURETTES

Colored Cervical Curette

Catalog Number	Tip Size (#)	Tip	Color	Length (in.)
MD309R	000000	REVERSE ANG	RED	10
MD310R	00000	REVERSE ANG	SILVER	10
MD311R	0000	REVERSE ANG	PURPLE	10
MD312R	000	REVERSE ANG	TURQUOISE	10
MD313R	00	REVERSE ANG	GOLD	10
MD314R	0	REVERSE ANG	GREEN	10
MD315R	1	REVERSE ANG	ORANGE	10
MD316R	2	REVERSE ANG	BLUE	10
MD317R	3	REVERSE ANG	GRAY	10
MD318R	4	REVERSE ANG	BLACK	10
MD319R	5	REVERSE ANG	BLACK	10
MD320R	6	REVERSE ANG	BLACK	10
MD383R	000000	ANG 90°	RED	10
MD384R	00000	ANG 90°	SILVER	10
MD385R	0000	ANG 90°	PURPLE	10
MD386R	000	ANG 90°	TURQUOISE	10
MD387R	00	ANG 90°	GOLD	10
MD388R	0	ANG 90°	GREEN	10
MD389R	1	ANG 90°	ORANGE	10
MD390R	2	ANG 90°	BLUE	10
MD391R	3	ANG 90°	GRAY	10
MD392R	4	ANG 90°	BLACK	10
MD393R	5	ANG 90°	BLACK	10
MD394R	6	ANG 90°	BLACK	10



ORTHOPEDICS/SPINE CURETTES

Colored Lumbar Curette




Catalog Number	Tip Size (#)	Tip	Color	Length (in.)
MD500R	000000	STR	RED	13 3/4
MD501R	00000	STR	SILVER	13 3/4
MD502R	0000	STR	PURPLE	13 3/4
MD503R	000	STR	TURQUOISE	13 3/4
MD504R	00	STR	GOLD	13 3/4
MD505R	0	STR	GREEN	13 3/4
MD506R	1	STR	ORANGE	13 3/4
MD507R	2	STR	BLUE	13 3/4
MD508R	3	STR	GRAY	13 3/4
MD509R	4	STR	BLACK	13 3/4
MD510R	5	STR	BLACK	13 3/4
MD511R	6	STR	BLACK	13 3/4
MD620R	000000	ANG	RED	13 3/4
MD621R	00000	ANG	SILVER	13 3/4
MD622R	0000	ANG	PURPLE	13 3/4
MD623R	000	ANG	TURQUOISE	13 3/4
MD624R	00	ANG	GOLD	13 3/4
MD625R	0	ANG	GREEN	13 3/4
MD626R	1	ANG	ORANGE	13 3/4
MD627R	2	ANG	BLUE	13 3/4
MD628R	3	ANG	GRAY	13 3/4
MD629R	4	ANG	BLACK	13 3/4
MD630R	5	ANG	BLACK	13 3/4
MD631R	6	ANG	BLACK	13 3/4
MD632R	000000	REVERSE ANG	RED	13 3/4
MD633R	00000	REVERSE ANG	SILVER	13 3/4
MD634R	0000	REVERSE ANG	PURPLE	13 3/4
MD635R	000	REVERSE ANG	TURQUOISE	13 3/4
MD636R	00	REVERSE ANG	GOLD	13 3/4
MD637R	0	REVERSE ANG	GREEN	13 3/4
MD638R	1	REVERSE ANG	ORANGE	13 3/4
MD639R	2	REVERSE ANG	BLUE	13 3/4
MD640R	3	REVERSE ANG	GRAY	13 3/4
MD641R	4	REVERSE ANG	BLACK	13 3/4
MD642R	5	REVERSE ANG	BLACK	13 3/4
MD643R	6	REVERSE ANG	BLACK	13 3/4



ORTHOPEDICS/SPINE CURETTES

Bayonet Curette

Catalog Number	Tip	Tip Size (#)	Length (in.)
MD660R	STR	0	9
MD661R	STR	00	9
MD662R	STR	000	9
MD663R	STR	0000	9
MD664R	STR	00000	9
MD655R	REVERSE ANG UP	0	9
MD656R	REVERSE ANG UP	00	9
MD657R	REVERSE ANG UP	000	9
MD658R	REVERSE ANG UP	0000	9
MD659R	REVERSE ANG UP	00000	9
MD665R	REVERSE ANG DOWN	0	9
MD666R	REVERSE ANG DOWN	00	9
MD667R	REVERSE ANG DOWN	000	9
MD668R	REVERSE ANG DOWN	0000	9
MD669R	REVERSE ANG DOWN	00000	9



ORTHOPEDICS/SPINE CURETTES

Black Handle Curette

Catalog Number	Tip	Tip Size (#)	Handle	Length (in.)
MN012R	STR	000	8 1/2	13 1/4
MN013R	STR	00	8 1/2	13 1/4
MN014R	STR	0	8 1/2	13 1/4
MN015R	STR	1	8 1/2	13 1/4
MN016R	STR	2	8 1/2	13 1/4
MN017R	STR	3	8 1/2	13 1/4
MN018R	STR	4	8 1/2	13 1/4
MN019R	STR	5	8 1/2	13 1/4
MN020R	STR	6	8 1/2	13 1/4
MN021R	STR	7	8 1/2	13 1/4
MN022R	ANG	000	8 1/2	13 1/4
MN023R	ANG	00	8 1/2	13 1/4
MN024R	ANG	0	8 1/2	13 1/4
MN025R	ANG	1	8 1/2	13 1/4
MN026R	ANG	2	8 1/2	13 1/4
MN027R	ANG	3	8 1/2	13 1/4
MN028R	ANG	4	8 1/2	13 1/4
MN029R	ANG	5	8 1/2	13 1/4
MN030R	ANG	6	8 1/2	13 1/4
MN031R	ANG	7	8 1/2	13 1/4
MN032R	REVERSE ANG	000	8 1/2	13 1/4
MN033R	REVERSE ANG	00	8 1/2	13 1/4
MN034R	REVERSE ANG	0	8 1/2	13 1/4
MN035R	REVERSE ANG	1	8 1/2	13 1/4
MN036R	REVERSE ANG	2	8 1/2	13 1/4
MN037R	REVERSE ANG	3	8 1/2	13 1/4
MN038R	REVERSE ANG	4	8 1/2	13 1/4
MN039R	REVERSE ANG	5	8 1/2	13 1/4
MN040R	REVERSE ANG	6	8 1/2	13 1/4
MN041R	REVERSE ANG	7	8 1/2	13 1/4



MN012R-MN021R



MN022R-MN031R



MN032R-MN041R



ORTHOPEDICS/SPINE CURETTES

XL Black Handled Curette

Catalog Number	Tip	Tip Size (#)	Handle	Length (in.)
MN042R	STR	000	9	15
MN043R	STR	00	9	15
MN044R	STR	0	9	15
MN045R	STR	1	9	15
MN046R	STR	2	9	15
MN047R	STR	3	9	15
MN048R	STR	4	9	15
MN049R	STR	5	9	15
MN050R	STR	6	9	15
MN051R	STR	7	9	15
MN052R	ANG	000	9	15
MN053R	ANG	00	0	15
MN054R	ANG	0	0	15
MN055R	ANG	1	9	15
MN056R	ANG	2	9	15
MN057R	ANG	3	9	15
MN058R	ANG	4	9	15
MN059R	ANG	5	9	15
MN060R	ANG	6	9	15
MN061R	ANG	7	9	15
MN062R	REVERSE ANG	000	9	15
MN063R	REVERSE ANG	00	9	15
MN064R	REVERSE ANG	0	9	15
MN065R	REVERSE ANG	1	9	15
MN066R	REVERSE ANG	2	9	15
MN067R	REVERSE ANG	3	9	15
MN068R	REVERSE ANG	4	9	15
MN069R	REVERSE ANG	5	9	15
MN070R	REVERSE ANG	6	9	15
MN071R	REVERSE ANG	7	9	15



1/3




MN042R-MN051R



MN052R-MN061R




MN062R-MN071R




ORTHOPEDICS/SPINE CURETTES AND MICRODISSECTING INSTRUMENTS

Micro Pro le Instruments

Catalog Number	Description	Tip Size (#)	Tip	Length (in.)
MN087R	DISSECTOR	3 MM	STR	10 1/4
MN088R	DISSECTOR	5 MM	STR	10 1/4
MN077R	CURETTE	1	ANG	10 1/4
MN079R	CURETTE	2	ANG	10 1/4
MN081R	CURETTE	3	ANG	10 1/4
MN083R	CURETTE	4	ANG	10 1/4
MN085R	CURETTE	5	ANG	10 1/4
MN078R	CURETTE	1	REV-ANG	10 1/4
MN080R	CURETTE	2	REV-ANG	10 1/4
MN082R	CURETTE	3	REV-ANG	10 1/4
MN084R	CURETTE	4	REV-ANG	10 1/4
MN086R	CURETTE	5	REV-ANG	10 1/4



ORTHOPEDICS/SPINE MICRODISCECTOMY INSTRUMENTS




ORTHOPEDICS/SPINE MICRODISCECTOMY INSTRUMENTS

Microdiscectomy Cervical Instruments

Catalog Number	Tip Size (#)	Tip	Length (in.)
MF250	COMPLETE SET (Includes One Each of All Items Listed)		
MF251R	0	FWD STR	8 1/2
MF252R	00	FWD STR	8 1/2
MF253R	000	FWD STR	8 1/2
MF254R	0000	FWD STR	8 1/2
MF255R	Bone Cutter	FWD STR	8 1/2
MF256R	0	FWD ANG	8 1/2
MF257R	00	FWD ANG	8 1/2
MF258R	000	FWD ANG	8 1/2
MF259R	0000	FWD ANG	8 1/2
MF260R	Bone Cutter	FWD ANG	8 1/2
MF261R	0	BACK STR	8 1/2
MF262R	00	BACK STR	8 1/2
MF263R	000	BACK STR	8 1/2
MF264R	0000	BACK STR	8 1/2
MF265R	Bone Cutter	BACK STR	8 1/2
MF266R	0	BACK ANG	8 1/2
MF267R	00	BACK ANG	8 1/2
MF268R	000	BACK ANG	8 1/2
MF269R	0000	BACK ANG	8 1/2
MF270R	Bone Cutter	BACK ANG	8 1/2
MF321R		MICRO NERVE HOOK	8 1/2
MF322R		ULTRA MICRO NERVE HOOK	8 1/2
MF323R		BONE WAXER	8 1/2
MF320	Case - 12 Instruments (Set includes 2 cases)		

Microdiscectomy Cervical Long Instruments

Catalog Number	Tip Size (#)	Tip	Length (in.)
MF271	COMPLETE SET (Includes One Each of All Items Listed)		
MF272R	0	FWD STR	9 1/2
MF273R	00	FWD STR	9 1/2
MF274R	000	FWD STR	9 1/2
MF275R	0000	FWD STR	9 1/2
MF276R	Bone Cutter	FWD STR	9 1/2
MF277R	0	FWD ANG	9 1/2
MF278R	00	FWD ANG	9 1/2
MF279R	000	FWD ANG	9 1/2
MF280R	0000	FWD ANG	9 1/2
MF281R	Bone Cutter	FWD ANG	9 1/2
MF282R	0	BACK STR	9 1/2
MF283R	00	BACK STR	9 1/2
MF284R	000	BACK STR	9 1/2
MF285R	0000	BACK STR	9 1/2
MF286R	Bone Cutter	BACK STR	9 1/2
MF287R	0	BACK ANG	9 1/2
MF288R	00	BACK ANG	9 1/2
MF289R	000	BACK ANG	9 1/2
MF290R	0000	BACK ANG	9 1/2
MF291R	Bone Cutter	BACK ANG	9 1/2
MF292R		NERVE HOOK	9 1/2
MF293R		ULTRA MICRO	9 1/2
MF294R		NERVE HOOK	9 1/2
MF320	Case - 12 Instruments (Set includes 2 cases)		



Tip Size 1:1


Tip Size (#)



ORTHOPEDICS/SPINE MICRODISCECTOMY INSTRUMENTS



Microdiscectomy Lumbar Instruments

Catalog Number	Tip Size (#)	Tip	Length (in.)
MF295	COMPLETE SET (Includes One Each of All Items Listed)		
MF296R	0	FWD STR	9 1/2
MF297R	1	FWD STR	9 1/2
MF298R	2	FWD STR	9 1/2
MF299R	0	FWD ANG	9 1/2
MF300R	1	FWD ANG	9 1/2
MF301R	2	FWD ANG	9 1/2
MF302R	0	FWD DOWN	9 1/2
MF303R	1	FWD DOWN	9 1/2
MF304R	2	FWD DOWN	9 1/2
MF305R	0	BACK STR	9 1/2
MF306R	1	BACK STR	9 1/2
MF307R	2	BACK STR	9 1/2
MF308R	0	BACK ANG	9 1/2
MF309R	1	BACK ANG	9 1/2
MF310R	2	BACK ANG	9 1/2
MF311R	0	BACK DOWN	9 1/2
MF312R	1	BACK DOWN	9 1/2
MF313R	2	BACK DOWN	9 1/2
MF314R		#1 PENFIELD	9 1/2
MF315R		#2 PENFIELD	9 1/2
MF316R		NERVE HOOK	9 1/2
MF317R		#4 PENFIELD	9 1/2
MF318R		BONE WAXER	9 1/2
MF319R		SPINE PROBE	9 1/2
MF320	Case - 12 Instruments (Set includes 2 cases)		




Tip Size 1:1

Tip Size (#)



ORTHOPEDICS/SPINE MICRODISCECTOMY INSTRUMENTS



Storage Case for Microdiscectomy
Instruments

Catalog Number
MF320

ORTHOPEDICS/SPINE DISSECTORS AND PROBES



LEMMON Intima Dissector

Catalog Number	Tip	Length (in.)
FB141R	SLIGHT-CVD	6 1/2
FB142R	CVD	6 1/2
FB143R	HEAVY-CVD	6 1/2






BUNNELL Dissecting Probe

Catalog Number	Tip	Length (in.)
MB456R	ANG	5 1/2




WATSON-CHEYNE Dissector, Double Ended

Catalog Number	Blade Size (mm)	Length (in.)
FK292R	2 / 2.5	5 1/8
FK293R	2 / 2.7	7 1/2



ORTHOPEDICS/SPINE RONGEURS

Rongeur jaw measurement guide



F0399R



STELLBRINK Synovectomy Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0399R	2.0 x 39.5	SLIGHT-CVD	6 7/8
F0400R	2 x 33	CVD	6 3/4



F0400R



Synovectomy Rongeur




Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0599R	1.8 x 24	CVD	5 1/8



F0599R



ORTHOPEDICS/SPINE RONGEURS



FRIEDMANN Micro Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
FO406R	1.3 x 19	STR	5 1/2
FO407R	1.3 x 20	CVD	5 1/2

LUER-FRIEDMANN Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
FO409R	2.2 x 22	CVD	5 3/4

Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
FO410R	3.8 x 19	CVD	5 1/2



ORTHOPEDICS/SPINE RONGEURS

HAKANSSON Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0428R	4.5 x 24	CVD	6 3/4



LUER Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0460R	7.5 x 25	STR	6 7/8
F0462R	9.5 x 33	STR	7
F0411R	4.2 x 20	CVD	6
F0412R	4.2 x 26	STRONG CVD	6 1/8
F0461R	7.5 x 24.5	CVD	6 7/8
F0463R	9.0 x 30	CVD	7



ORTHOPEDICS/SPINE RONGEURS

BLUMENTHAL Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0415R	3.5 x 25	30° CVD	6 1/8
F0416R	3.5 x 22	45° CVD	6



JANSEN Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0430R	4.2 x 35	STR	6 7/8
F0432R	4.5 x 30	CVD	6 7/8
F0442R	3.8 x 37	ANG	7 1/2



ORTHOPEDICS/SPINE RONGEURS




ALEXANDER Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
FO440R	3.2 x 27	ANG	5 1/2



LEMPERT Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
MD617	2.0 x 23	DEL, STR	6 1/2
FO450R	2.2 x 27	STR	7 7/8
MG059R	2.5 x 30	CVD	6 1/4
FO451R	2.2 x 27	CVD	7 7/8





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

ORTHOPEDICS/SPINE RONGEURS

BEYER Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0453R	3.2 x 32	CVD	7





MEAD Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0455R	3.5 x 30	CVD	6 1/2



BANE-HARTMANN Bone Rongeur


Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0457R	4.0 x 32	CVD	6 7/8



ORTHOPEDICS/SPINE RONGEURS

TROTTER Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
FO467R	6.5 x 24	CVD	8



FO467R

FO467R

OLIVECRONA Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
FO470R	5.0 x 36.5	CVD	8



FO470R

FO470R

GULEKE Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
FO472R	6.5 x 27	CVD	8




FO472R

FO472R

ORTHOPEDICS/SPINE RONGEURS

PENNYBACKER Bone Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
FO475R	4.0 x 26	STR	8




FO475R



FO475R

ADSON Cranial Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
MD615	21.0 x 31.5	STR	7 7/8



MD615



MD615

KLEINERT-KUTZ Rongeur, Double Action


Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
MG060R	3.0 x 15.0	SLIGHT, CVD	6 1/8



MG060R




MG060R



ORTHOPEDICS/SPINE RONGEURS

BOEHLER Bone Rongeur, Double Action

Catalog Number	Jaw Size (mm)	Jaw	Jaw End	Length (in.)
F0500R	3.0 x 27.0	CVD	NARR	6 1/8
F0502R	3.5 x 30.0	CVD	WIDE	6 1/8



BEYER Bone Rongeur, Double Action

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0514R	4.0 x 30.0	CVD	7



MAYFIELD Bone Rongeur, Double Action

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0511R	3.5 x 23	CVD	6 7/8





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

ORTHOPEDICS/SPINE RONGEURS




F0515R



F0515R

ZAUFAL-JANSEN Bone Rongeur, Double Action


Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0515R	4.5 x 29.0	CVD	7




MG056R-MG058R

RUSKIN Bone Rongeur, Double Action


Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
MG056R	2	STR, MINI	6 1/4
MG057R	2	CVD, MINI	6 1/4
MG058R	3	CVD, MINI	6 1/4
F0518R	5	CVD	7 1/2
MG579R	6	CVD	7 1/4



MG056R




MG057R




MG058R



MG579R, F0518R



MG579R, F0518R




ORTHOPEDICS/SPINE RONGEURS



MARQUARDT Bone Rongeur, Double Action

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0520R	3.5 x 15.5	CVD	7 7/8



ROETTGEN-RUSKIN Bone Rongeur, Double Action

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0521R	6.3 x 30	CVD	9 1/2



SMITH-PETERSON Rongeur

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
MG577R	3	STR	9 1/2

ORTHOPEDICS/SPINE RONGEURS



F0522R-F0523R

ECHLIN Bone Rongeur, Double Action, Duckbill


Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0522R	2 x 10	ANG	9
F0523R	3 x 10	ANG	9



F0522R




F0523R



MG581R

LEKSELL Bone Rongeur, Double Action, Duckbill

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
MG581R	5	FULL CVD	9
F0526R	5 x 16	ANG	9 5/8
F0527R	3 x 16	ANG	9 5/8



MG581R




F0526R-F0527R



F0526R



F0527R




ORTHOPEDICS/SPINE RONGEURS




FRYKHOLM Bone Rongeur, Double Action

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0525R	4.5 x 25	ANG	9 1/2



STILLE Bone Rongeur, Double Action, Duckbill

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0530R	7.0 x 25	ANG	9



LEKSELL-STILLE Bone Rongeur, Double Action, Duckbill

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0532R	7.5 x 29.5	ANG	9 1/2

ORTHOPEDICS/SPINE RONGEURS


SYPERT Bone Rongeur, Double Action

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
FG052R	8.0 x 35.0	ANG	14 1/4






STILLE-RUSKIN Bone Rongeur, Double Action

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0534R	8.0 x 34.0	CVD	9 1/2



SAUERBRUCH-CORYLLOS Rongeur, Double Action

Catalog Number	Length (in.)
F0570R	12 1/4
F0572R	12 1/4



ORTHOPEDICS/SPINE RONGEURS



F0550R-F0551R

LUER-STILLE Bone Rongeur, Double Action


Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0550R	10.5 x 27.5	STR	8 7/8
F0551R	10 x 26.5	CVD	8 7/8



F0550R



F0551R



F0555R

SEMB Bone Rongeur, Double Action

Catalog Number	Jaw Size (mm)	Jaw	Length (in.)
F0555R	16.5 x 38.5	CVD	9 1/2



1/1



F0555R