



Wholesale Supplier Plastic_Surgery_Instrument A to Z Hospital Equipment and Surgical Products **B to C**



Plastic Surgery

Scissors	748 - 768
Needle Holders	768 - 772
Forceps	773 - 782
Hooks and Retractors	783 - 793
Aluminum Hand	794
Knives	795 - 797
Suction Instruments	798
Bone Instruments	799 - 807
Micro Vessel Clips and Clamps	808 - 810

8275140877

9822688071

info@amramedtech.com



**Sr No.29/1/2, Adarsh Nagar,
Uruli Devachi, Hadapsar, Pune
Branches :- All Maharashtra**











amramedtech.com



PLASTIC SURGERY SCISSORS

Iris Scissors, Delicate

Catalog Number	Tip	Blade	Handle	Length (in.)
BC004R	S/S	STR	SQUARE	4 1/2
BC005R	S/S	CVD	SQUARE	4 1/2
BC056R	S/S	STR	ROUND	3 3/4
BC060R	S/S	STR	ROUND	4 1/8
BC064R	S/S	STR	ROUND	4 3/4
BC057R	S/S	CVD	ROUND	3 3/4
BC061R	S/S	CVD	ROUND	4 1/8
BC065R	S/S	CVD	ROUND	4 3/4



PLASTIC SURGERY SCISSORS








Iris Scissors, Ribbon Handle

Catalog Number	Tip	Blade	Length (in.)
MB605R	S/S	STR	4
BC190R	S/S	CVD	4



Suture Scissors

Catalog Number	Tip	Blade	Length (in.)
BC130R	S/S	STR	4 3/4
BC140R	S/S	STR	5 1/8
BC144R	B/B	STR	5 1/8
BC135R	B/B	CVD, DEL	4 3/4
BC131R	S/S	CVD	4 3/4
BC141R	S/S	CVD	5 1/8
BC145R	B/B	CVD	5 1/8






PLASTIC SURGERY SCISSORS

COTTLE-MASING Plastic Scissors

Catalog Number	Tip	Blade	Length (in.)
OK365R	S/S	CVD	4 1/8
OK367R	B/B	CVD	4 1/8



OK365R-OK367R




OK365R




OK367R

STEVENS Tenotomy Scissors


Catalog Number	Tip	Blade	Handle	Weight	length (in.)
MA470R	B/B	STR	RIBBON LIGHT		3 3/4
MA468R	B/B	STR	RIBBON HEAVY		3 7/8
MA469R	B/B	CVD	RIBBON HEAVY		3 3/4



MA470R-MA468R




MA470R



MA468R



MA469R



MA469R



PLASTIC SURGERY











MILLIMETERS



PLASTIC SURGERY SCISSORS

**STEVENS Tenotomy Scissors
(continued)**

Catalog Number	Tip	Blade	Handle	Length (in.)
BC008R	S/S	STR	SQUARE	4 1/2
BC010R	B/B	STR	SQUARE	4 1/2
BC009R	S/S	CVD	SQUARE	4 1/2
BC011R	B/B	CVD	SQUARE	4 1/2
BC170R	S/S	STR	ROUND	4 1/4
BC172R	B/B	STR	ROUND	4 1/4
MA467R	B/B	STR	ROUND	5 3/4
BC171R	S/S	CVD	ROUND	4 1/4
BC173R	B/B	CVD	ROUND	4 1/4
MA466R	B/B	CVD	ROUND	5 3/4







PLASTIC SURGERY SCISSORS



STEVENS Tenotomy Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC194R	B/B	STR	5
BC195R	B/B	CVD	5





REYNOLDS-JAMESON Tenotomy Scissors, Delicate

Catalog Number	Tip	Blade	Handle	Length (in.)
BC015R	B/B	CVD	SQUARE	5 1/2





REYNOLDS Tenotomy Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC180R	B/B	CVD	5 3/4
BC181R	B/B	CVD	6 7/8




JAMESON Tenotomy Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC177R	B/B	CVD	6



PLASTIC SURGERY SCISSORS




JAMESON SuperCut Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC911R	B/B	CVD	5 1/8
BC915R	B/B	CVD	6
BC919R	B/B	CVD	7




POTTS Tenotomy Scissors

Catalog Number	Tip	Blade	Length (in.)
BC178R	B/B	CVD	6 7/8
BC179R	B/B	CVD	8



Nerve Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC628R	S/S	CVD	6












PLASTIC SURGERY SCISSORS

JOSEPH Scissors, Delicate

Catalog Number	Tip	Blade	Handle	Length (in.)
BC006R	S/S	STR	SQUARE	6
BC007R	S/S	CVD	SQUARE	6
BC156R	S/S	STR	ROUND	6
BC157R	S/S	CVD	ROUND	6





KILNER Scissors

Catalog Number	Tip	Blade	Length (in.)
BC162R	S/S	STR	4 1/2
BC166R	FLAT	CVD	5
BC168R	FLAT	CVD	6






PLASTIC SURGERY SCISSORS



Strabismus Scissors

Catalog Number	Tip	Blade	Length (in.)
BC164R	B/B	STR	4 1/2
BC165R	B/B	CVD	4 1/2



LITTLER Scissors

Catalog Number	Tip	Blade	Length (in.)
BC163R	W/EYE	CVD	4 1/2










PLASTIC SURGERY SCISSORS


METZENBAUM SuperCut Scissors

Catalog Number	Tip	Blade	Length (in.)
BC945R	B/B	CVD	9
BC947R	B/B	CVD	10 1/4
BC949R	B/B	CVD	11 1/4
BC951R	B/B	CVD	11 3/4
BC953R	B/B	CVD	14




GRAEFE Scissors

Catalog Number	Tip	Blade	Length (in.)
OC415R	B/B	LAT-ANG	4







PLASTIC SURGERY SCISSORS



OK373R-OK372R



OK373R




OK372R

FOMON Lower Lateral Scissors

Catalog Number	Tip	Blade	Length (in.)
OK373R	B/B	HEAVY-CVD	4 1/2
OK372R	B/B	HEAVY-CVD	5 1/4



BC801R-BC802R



BC801R




BC802R

SPENCER Suture Scissors

Catalog Number	Blade	Length (in.)
BC801R	STR	3 1/2
BC802R	STR	4 1/2



BC803R



BC803R

LITTAUER Ligature Scissors

Catalog Number	Blade	Length (in.)
BC803R	STR	5 1/4





PLASTIC SURGERY SCISSORS

NORTHBENT Ligature Scissors

Catalog Number	Tip	Blade	Length (in.)
BC805R	B/B	LAT-CVD	4 3/4

Cartilage Scissors






Catalog Number	Tip	Blade	Length (in.)
BC672R	S/S	SERR	7 1/4

KAYE Facelift Scissors

Catalog Number	Tip	Blade	Length (in.)
MA458R	SEMI-SHARP	CVD, SERR	6
MB109R	SEMI-SHARP	CVD, SERR	7 1/2

METZENBAUM-LIPSETT Scissors

Catalog Number	Tip	Blade	Length (in.)
MB150R	B/B	CVD	5 1/2






PLASTIC SURGERY SCISSORS



BC184R-BC189R

GORNEY Scissors

Catalog Number	Tip	Blade	Length (in.)
BC184R	B/B	STR	5
BC188R	B/B	STR	7 3/4
BC189R	B/B	SLIGHT-CVD	7 3/4



BC184R



BC188R



BC189R



MB632R

GORNEY-FREEMAN SuperCut Scissors


Catalog Number	Tip	Blade	Length (in.)
MB632R	BEVELED BLUNT	STR	7 1/4
MB633R	BEVELED BLUNT	CVD	7 1/4



MB632R



MB633R




MB633R




PLASTIC SURGERY



PLASTIC SURGERY DUROTIP™ TC SCISSORS



DO208R-DO209R



DO208R




DO209R

DUROTIP TC Delicate Scissors


Catalog Number	Tip	Blade	Length (in.)
DO208R	S/S	STR	4 3/4
DO209R	S/S	CVD	4 3/4



BC210R-BC211R



BC210R




BC211R

DUROTIP TC Iris Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC210	S/S	STR	4 1/4
R	S/S	CVD	4 1/4

BC211
R



DO250R



DO250R

DUROTIP TC GOLDMAN-FOX Tenotomy Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
DO250R	S/S	CVD, 1/SERR	5 1/8



MB198R-MB197R



MB198R



MB197R

DUROTIP TC Par Scissors




Catalog Number	Tip	Length (in.)
MB198R	CVD, BLUNT	4 1/2
MB197R	CVD, BLUNT	5 3/4




PLASTIC SURGERY DUROTIP™ TC SCISSORS

**DUROTIP TC GORNEY-FREEMAN
Facelift Scissors**

Catalog Number	Blade	Length (in.)
MB648R	STR, SERR	7 3/4
MB649R	CVD, SERR	7 3/4



Wavecut
Wavecut scissors have one serrated cutting edge to help prevent tissue from slipping.
The W after the item number indicates the scissors are the wavecut style, e.g. BC211W.



DUROTIP TC Iris Scissors, Wavecut

Catalog Number	Tip	Blade	Length (in.)
BC211W	S/S	CVD	4 1/4






PLASTIC SURGERY DUROTIP™ TC SCISSORS




**DUROTIP TC Baby METZENBAUM
Wavecut Dissecting Scissors, Delicate**

Catalog Number	Tip	Blade	Length (in.)
BC259W	B/B	CVD	5 3/4




**DUROTIP TC METZENBAUM
Wavecut Dissecting Scissors**

Catalog Number	Tip	Blade	Length (in.)
BC261W	B/B	CVD	5 3/4
BC263W	B/B	CVD	7



**DUROTIP TC METZENBAUM
Wavecut Delicate Scissors**


Catalog Number	Tip	Blade	Length (in.)
BC271W	B/B	CVD	7



PLASTIC SURGERY SUPERCUT SCISSORS




BC903R




BC903R

SuperCut Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC903R	S/S	CVD	5



BC909R



BC909R

GOLDMANN-FOX SuperCut Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC909R	S/S	CVD	5 1/8






PLASTIC SURGERY SUPERCUT SCISSORS


GORNEY SuperCut Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
BC920R	B/B	STR	7 3/4
BC921R	B/B	CVD	7 3/4




KAYE SuperCut Scissors, Delicate

Catalog Number	Blade	Length (in.)
BC933R	CVD	7 3/4



FOMON SuperCut Scissors


Catalog Number	Tip	Blade	Length (in.)
BC939R	B/B	CVD	4 1/2




PLASTIC SURGERY

MILLIMETERS

PLASTIC SURGERY NOIR™ SUPREME CUT SCISSORS



Microtome Edge



Wavecut

Noir Supreme Cut Scissors

■ **Noir**

The combination of the renowned Aesculap Tungsten Carbide technology with the Aesculap Noir coating results in an extremely high surface hardness and therefore high wear resistance in daily hospital use.

■ **Microtome edge**

The extremely sharp cutting edge of the scissors blade facilitates gentle cutting of tissue with minimal use of force.

■ **Wavecut helps prevent tissue from slipping out during the cutting process.**

■ **The WB after the item number indicates the Noir scissors are the wavecut style, e.g. BC211WB.**

For more information see our brochure number **DOC637**.



BC210B-BC211WB

1/2



BC210B

1/1



BC211B, BC211WB

1/1

Noir Iris Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC210B	S/S	STR	4 1/4
BC211B	S/S	CVD	4 1/4
BC211WB	S/S	CVD	4 1/4



BC257B-BC257WB

1/2



BC257B-BC257WB

1/1

Noir Strabismus Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC257B	B/B	CVD	4 1/2
BC257WB	B/B	CVD	4 1/2



BC259B-BC259WB

1/2



BC259B-BC259WB

1/1


Noir Baby METZENBAUM Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC259B	B/B	CVD	5 3/4
BC259WB	B/B	CVD	5 3/4





PLASTIC SURGERY NOIR™ SUPREME CUT SCISSORS



Noir METZENBAUM Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC261B	B/B	CVD	5 3/4
BC261WB	B/B	CVD	5 3/4



Noir METZENBAUM Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC263B	B/B	CVD	7
BC263WB	B/B	CVD	7
BC271B	B/B	CVD, DEL	7
BC271WB	B/B	CVD, DEL	7



Noir TOENNIS-ADSON Dissecting Scissors

Catalog Number	Tip	Blade	Length (in.)
BC273B	B/B	CVD, DEL	6 7/8
BC273WB	B/B	CVD, DEL	6 7/8



PLASTIC SURGERY MICRO SCISSORS



FM010R-FM040R


1/2




1/1



1/1



FM010R



FM040R



FM011R-FM041R

1/2



1/1



1/1



1/1




1/1



1/1



FM011R




FM012R



FM013R



FM014R



FM041R




FD251R-FD252R

1/2




1/1



FD252R



1/1



FD251R

MUELLER Micro Scissors

Catalog Number	Tip	Handle Color	Blade	Length (in.)
FM010R	S/S	BLUE	STR	4 3/4
FM040R	S/S	BLUE	STR	6 3/4
FM011R	S/S	BLUE	CVD	4 3/4
FM012R	S/B	BLUE	CVD	4 3/4
FM013R	B/B	BLUE	CVD	4 3/4
FM014R	S/BALL-PT	BLUE	CVD	4 3/4
FM041R	S/S	BLUE	CVD	6 3/4

Micro Scissors, Golf Ball Handle




Catalog Number	Blade	Length (in.)
FD252R	SLIGHT CVD, 1/SERR	7
FD251R	CVD	7



PLASTIC SURGERY



PLASTIC SURGERY MICRO SCISSORS AND NEEDLE HOLDERS



Micro Scissors, Golf Ball Handle

Catalog Number	Blade	Length (in.)
FM042R	STR	7
FM043R	CVD	7

YASARGIL Micro Vessel Scissors

Catalog Number	Blade	Handle	Length (in.)
FD022R	STR	FLAT, SERR	7 1/4
FD023R	CVD	FLAT, SERR	7 1/4


HEGAR-OLSEN Needle Holder

Catalog Number	Tip	Length (in.)
MG521R	SERR	4 3/4



PLASTIC SURGERY DUROGRIP™ TC NEEDLE HOLDERS






Aesculap features Needle Holders with Tungsten Carbide (TC) inserts.
A smooth jaw and three different serration patterns are available, depending on the size of the suture needle.



We recommend carbide inserts as follows:

Jaw Surface Smooth 0.2 mm
pitch of serration 0.4 mm pitch
of serration 0.5 mm pitch of
serration

Suture Material
9/0 - 11/0
6/0 - 10/0
4/0 - 6/0
up to 3/0



**DUROGRIP TC CASTROVIEJO
Micro Needle Holder, Flat Handle**

Catalog Number	Jaw	Length (in.)
BM001R	SMOOTH	5 3/4
BM002R	SERR	5 3/4
BM004R	SERR	7
BM005R	SERR	8 1/2





PLASTIC SURGERY DUROGRIP™ TC NEEDLE HOLDERS


**DUROGRIP TC LEWIN Needle Holder/
Scissors**

Catalog Number	Jaw	Length (in.)
BM127R	SMOOTH	5



BM127R


1/2



BM127R

**DUROGRIP TC HEGAR-OLSON
Needle Holder/Scissors**

Catalog Number	Jaw	Length (in.)
MG005R	SMOOTH	4 3/4
MG006R	SERR	4 3/4
BM128R	SERR	5 3/4
BM124R	SERR	6 1/2



MG005R-MG006R, BM124R, BM128R




MG005R

1/1




MG006R

1/1



BM124R

0.4



BM128R

1/1


0.4



PLASTIC SURGERY DUROGRIP™ TC NEEDLE HOLDERS



BM129R



**DUROGRIP TC GILLIES Needle Holder/
Scissors**


Catalog Number	Jaw	Handle	Length (in.)
BM129R	SERR	OFF-SET	6 1/4




BM016R-BM021R

**DUROGRIP TC CRILE-WOOD
Needle Holder**


Catalog Number	Jaw	Length (in.)
BM016R	SERR	5 3/4
BM017R	SERR	7 1/4
BM019R	SERR	7 7/8
BM030R	SERR	10 5/8
BM021R	SERR	12




BM016R




BM019R




BM017R



BM021R



BM030R



BM030R



PLASTIC SURGERY



AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C




Certified Register

PLASTIC SURGERY DUROGRIP™ TC NEEDLE HOLDERS


**DUROGRIP TC CRILE-WOOD
Needle Holder**

Catalog Number	Jaw	Length (in.)
BM029R	SMOOTH	6



BM029R

1/2



PLASTIC SURGERY



PLASTIC SURGERY DRESSING AND TISSUE FORCEPS



MB295R



MB295R



MB294R



MB294R



BD505R



BD505R



BD329R



BD329R



BD330R



BD330R

CASTROVIEJO Suture Forceps

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
MB295R	0.12	1 x 2	4 1/8

CASTROVIEJO COLIBRI Type Corneal Forceps

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
MB294R	0.12	1 x 2	4 1/8

BONN Eye Forceps, W/ Platform

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
BD505R	0.2	1 x 2	3 3/4

Jewelers Forceps

Catalog Number	Tip Size (mm)	Tip	Length (in.)
BD329R	0.3	ANG	4 1/4

Micro Forceps #4 Jewelers

Catalog Number	Tip Size (mm)	Length (in.)
BD330R	0.2	4 1/4



PLASTIC SURGERY



PLASTIC SURGERY DRESSING AND TISSUE FORCEPS



BD345R



BD345R

Micro Jewelers Forceps

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
BD345R	0.4	1 x 2	5 1/4



BD336R



BD336R

VISE Micro Suture Forceps

Catalog Number	Tip Size (mm)	Tip	Length (in.)
BD336R	0.3	STR	5 1/4
BD337R	0.3	CVD	5 1/4



BD337R



BD337R

Micro Suture Forceps, W/ Platform

Catalog Number	Tip Size (mm)	Length (in.)
BD338R	0.3	4 3/4



BD338R



BD338R



BD340R



BD340R

AUSTIN Micro Suture Forceps, W/ Platform

Catalog Number	Tip Size (mm)	Length (in.)
BD340R	0.8	5 1/2

PLASTIC SURGERY DRESSING AND TISSUE FORCEPS



Micro Suture Forceps, Round Handle

Catalog Number	Tip Size (mm)	Tip	Length (in.)
FD284R	0.3	STR	6
FD286R	0.5	STR	7
FD285R	0.3	CVD	6
FD287R	0.5	CVD	7



FD284R




FD285R



FD286R




FD287R



MICRO-ADSON-BIEMER Forceps

Catalog Number	Teeth	Length (in.)
BD509R	1 x 2	4 3/4



BD509R






PLASTIC SURGERY DRESSING AND TISSUE FORCEPS



FM001R-FM006R

MUELLER Micro Forceps


Catalog Number	Tip Size (mm)	Handle Color	Teeth	Length (in.)
FM001R	0.15	BROWN	-	4 1/4
FM002R	0.20	BROWN	-	4 1/4
FM003R	0.30	BROWN	-	4 1/4
FM004R	0.50	BLACK	1 x 2	4 1/4
FM005R	1.20	GREY	RING	4 1/4
FM006R	1.20	GREY	TRI	4 1/4



BD260R

LAZAR Micro Suture Forceps


Catalog Number	Tip Size (mm)	Length (in.)
BD260R	0.5	6 1/8
BD261R	0.8	6 1/8



BD260R



BD261R



BD261R





PLASTIC SURGERY BONE AND TENDON FORCEPS


HERZ Tendon Seizing Forceps

Catalog Number	Jaw	Length (in.)
BF303R	CVD	6 1/8



JACKSON Tendon Seizing Forceps

Catalog Number	Teeth	Length (in.)
BF306R	4 x 5	6 1/8



F0940R, MB363C

BUNNELL Tendon Passing Probe

Catalog Number	Length (in.)
MB363C	7 3/8
F0940R	9



PLASTIC SURGERY




PLASTIC SURGERY BONE AND TENDON FORCEPS

Tendon Seizing Forceps, W/ 3 Pins


Catalog Number	Length (in.)
BF307R	4 3/4

BRAND Tendon Holding Forceps

Catalog Number	Jaw	Teeth	Length (in.)
BF324R	STR	1 x 2	6
BF326R	STR	1 x 2	7 1/4
BF315R	CVD	1 x 2	6
BF316R	CVD	1 x 2	7 1/4



BF307R



BF324R



BF326R



BF315R



BF316R





PLASTIC SURGERY

PLASTIC SURGERY BONE AND TENDON FORCEPS

KLEINERT-KUTZ Tendon Retriever

Catalog Number	Size	Shaft	Length (in.)
MG091R	-	FLEXIBLE	8
MB985R	SMALL	RIGID	8 1/4



CARROLL Tendon Pulling Forceps

Catalog Number	Length (in.)
MB474R	6






PLASTIC SURGERY BONE AND TENDON FORCEPS

Tendon Interlacing Forceps

Catalog Number	Jaw	Length (in.)
F0950R	STR	6
F0951R	CVD	6



F0950R-F0951R



F0950R



F0951R

Tendon Stripper

Catalog Number	Length (in.)
F0935R	11 3/4
F0939R	19 3/4



F0935R-F0939R



F0935R



1/1



F0939R



1/1



PLASTIC SURGERY SPLINTER FORCEPS



BD952R



BD952R



BD963R-BD965R



BD963R



BD964R



BD965R

Splinter Forceps

Catalog Number	Length (in.)
BD952R	3 1/8

FEILCHENFELD Splinter Forceps

Catalog Number	Length (in.)
BD963R	3
BD964R	3 3/4
BD965R	4 1/2

STIEGLITZ Splinter Forceps

Catalog Number	Jaw	Length (in.)
BD992R	STR	5 1/2
BD993R	CVD	5 1/2

PEET Splinter Forceps

Catalog Number	Length (in.)
BD995R	4 3/4



BD992R-BD993R



BD992R



BD993R



BD995R



BD995R




1/1





PLASTIC SURGERY HEMOSTATIC FORCEPS



BH100R




BH100R

CARREL Artery Forceps

Catalog Number	Length (in.)
BH100R	3 1/4



BH139R-BH143R



BH139R



BH143R

DANDY Forceps

Catalog Number	Jaw	Jaw Shape	Length (in.)
BH139R	SERR	CVD	5 1/2
BH143R	x 2 Teeth	LAT-CVD	5 1/2



PLASTIC SURGERY HOOKS AND RETRACTORS



OA332R-OA333R



OA332R



OA333R

GRAEFTE Hook

Catalog Number	Tip	Length (in.)
OA332R	SHARP	5
OA333R	BLUNT	5



OA334R-OA335R



OA334R




OA335R

GUTHRIE Double Hook


Catalog Number	Tip	Length (in.)
OA334R	SHARP	5
OA335R	SHARP	5 1/8




OA336R-OA338R



OA336R



OA337R



OA338R

ROLLET Hook Retractor

Catalog Number	Tip	Blade End	Length (in.)
OA336R	SHARP	2 PRONG	5 1/8
OA337R	SHARP	4 PRONG	5 1/4
OA338R	BLUNT	4 PRONG	5 1/4



OA340R



OA340R

PAUL Tendon Hook

Catalog Number	Tip	Blade End	Tip Size (mm)	Length (in.)
OA340R	BLUNT	3 PRONGS	11 x 8	4 3/4



OA351R



OA351R

JAEGER Hook Retractor, Fine

Catalog Number	Tip	Length (in.)
OA351R	SHARP	5 1/4





PLASTIC SURGERY HOOKS AND RETRACTORS



KLEINERT-KUTZ Skin Hook

Catalog Number	Size	Tip	Tip Size (mm)	Length (in.)
OL601R	1	SHARP	2.5	6 1/4
OL602R	2	SHARP	4.2	6 1/4
OL603R	3	SHARP	4.8	6 1/4

KILNER Skin Hook

Catalog Number	Tip	Length (in.)
OL610R	SHARP	6 1/4

GILLIES Skin Hook

Catalog Number	Hook End	Tip	Length (in.)
OL611R	SMALL	SHARP	7
OL612R	LARGE	SHARP	7

COTTLE Hook

Catalog Number	Tip	Length (in.)
OL619R	CVD, SHARP	6 1/4
OL617R	ANG, SHARP	6 1/4



PLASTIC SURGERY HOOKS AND RETRACTORS



OL620R, MA185R, OL622R-OL624R, MD581



JOSEPH Skin Hook

Catalog Number	Tip	Prongs	Tip Size (mm)	Length (in.)
OL620R	SHARP	1	-	6 1/2
MA185R	SHARP	2	2	6 1/4
OL622R	SHARP	2	3	6 1/2
OL623R	SHARP	2	6	6 1/2
OL624R	SHARP	2	8	6 1/2
MD581	SHARP	2	10	6 1/2



MG370R




MG370R

COTTLE-JOSEPH Double Hook

Catalog Number	Tip Size (mm)	Length (in.)
MG370R	12	5 3/4



OL621R



OL621R

FREER Double Hook


Catalog Number	Length (in.)
OL621R	6




OA357R-OA359R



OA357R



OA358R



OA359R

GRAEFE Strabismus Nerve/Vessel Hook

Catalog Number	Tip	Length (in.)
OA357R	SMALL	5 1/2
OA358R	MEDIUM	5 1/2
OA359R	LARGE	5 1/2











PLASTIC SURGERY HOOKS AND RETRACTORS

MANNERFELT Retractor

Catalog Number	Prongs	Tip	Blade Shape	Blade Size (mm)	Length (in.)
BT901R	1	SHARP	HOOK CVD	-	5 1/4
BT902R	1	SHARP	HOOK CVD	-	6 1/8
BT911R	1	SHARP	SMALL CVD	SMALL PRONG	6 1/8
BT921R	1	SHARP	SMALL CVD	LARGE PRONG	6 1/8
BT916R	1	BLUNT	SMALL CVD	SMALL PRONG	6 1/8
BT926R	1	BLUNT	SMALL CVD	LAGRGE PRONG	6 1/8
BT904R	2	SHARP	SMALL CVD	3.5 x 12	6 1/8
BT912R	2	SHARP	SMALL CVD	5 x 4	6 1/8
BT922R	2	SHARP	SMALL CVD	8 x 8	6 1/8
BT917R	2	BLUNT	SMALL CVD	5.5 x 4	6 1/8
BT927R	2	BLUNT	SMALL CVD	7 x 8	6 1/8
BT906R	2	SHORT ANG		8 x 12.5	6 1/8
BT913R	3	SHARP	SMALL CVD	5.5 x 7.2	6 1/8
BT923R	3	SHARP	SMALL CVD	7 x 11	6 1/8
BT918R	3	BLUNT	SMALL CVD	5.5 x 7.2	6 1/8
BT928R	3	BLUNT	SMALL CVD	7 x 11	6 1/8



PLASTIC SURGERY






PLASTIC SURGERY HOOKS AND RETRACTORS



LANGENBECK-MANNERFELT Retractor

Catalog Number	Blade Size (mm)	Length (in.)
BT943R	15 x 6	6 1/8
BT944R	20 x 6	6 1/8
BT945R	25 x 6	6 1/8
BT946R	30 x 8	6 1/8



BT943R



BT944R



BT945R




BT946R



MANNERFELT Saddle Hook

Catalog Number	Tip Size (mm)	Length (in.)
BT940R	16 x 16	6 1/8
BT948R	12 x 12	6 1/8
BT949R	12 x 14	6 1/8
BT951R	14 x 16	6 1/8
BT958R	12 x 12 FEN	6 1/8




BT940R



BT948R



BT949R



BT951R



BT958R



MANNERFELT Tendon Double Hook

Catalog Number	Length (in.)
BT964R	6 1/4



BT964R





PLASTIC SURGERY HOOKS AND RETRACTORS



BT966R



1/1



1/1


BT966R

MANNERFELT Double Spatula

Catalog Number	Length (in.)
BT966R	6 1/4



MB727R-MB733R



1/1

1/1

MB727R-MB733R

CRILE Retractor

Catalog Number	Size	Blade Size (mm)	Length (in.)
MB727R	1	3 x 11 / 6 x 21	4 1/4
MB728R	2	6 x 13 / 10 x 21	6 1/8
MB729R	3	16 x 25 / 17 x 35	6 1/4
MB730R	4	22 x 29 / 29 x 44	7
MB731R	5	29 x 41 / 35 x 51	7
MB732R	6	32 x 64 / 35 x 70	7 7/8
MB733R	7	41 x 76 / 51 x 79	8 7/8

TEBBETTS Breast Spatula, Semi-Malleable

Catalog Number	Length (in.)
MB384R	12 1/2




MB384R

1/2

PLASTIC SURGERY HOOKS AND RETRACTORS

TEBBETTS Breast Retractor

Catalog Number	Blade	Length (in.)
MB647R	DBL-ENDED	12



MB647R

MB647R

FREEMAN Facelift Retractor

Catalog Number	Length (in.)
MB934R	7



MB934R





PLASTIC SURGERY HOOKS AND RETRACTORS

COTTLE Alar (MASING) Retractor

Catalog Number	Blade Size (mm)	Length (in.)
OL211R	6 x 12	5 1/2
OL212R	10 x 12	5 1/2

DESMARRES Lid Retractor

Catalog Number	Blade Size (mm)	Length (in.)
OA310R	9 x 9	5 1/2
OA311R	9 x 12	5 1/2
OA312R	11 x 14	5 1/2
OA313R	12 x 16	5 1/2
OA314R	13 x 18	5 1/2



OL211R



OL211R



OL212R



OL212R



OA310R-OA314R



OA310R



OA311R




OA312R



OA313R



OA314R





PLASTIC SURGERY HOOKS AND RETRACTORS



BV001R

1/2




BV001R

1/1



BV002R

1/2




BV002R

1/1



BV003R

1/2



BV003R

1/1



BV005R

1/2



BV005R

1/1



BV008R, BV015R

1/2



BV015R

1/1



BV008R

1/1

FINSEN Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV001R	SHARP	2 x 3 PRONGS	2

FINSEN Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV002R	SHARP	3 x 4 PRONGS	2 3/4

FINSEN Modified Retractor


Catalog Number	Tip	Blade End	Length (in.)
BV003R	BLUNT	2 x 3 PRONGS	3

FINSEN Retractor

Catalog Number	Tip	Blade End	Length (in.)
BV005R	SHARP	3 x 4 PRONGS	2 3/4

HEISS Self-Retaining Retractor

Catalog Number	Tip	Prongs	Length (in.)
BV015R	SHARP	4 x 4	4
BV008R	BLUNT	4 x 4	4



PLASTIC SURGERY



PLASTIC SURGERY HOOKS AND RETRACTORS

Spring Wire Retractor

Catalog Number	Tip	Prongs	Length (in.)
BV006R	BLUNT	3 x 3	1 3/4
BV007R	BLUNT	3 x 3	2 1/4
MD471	SHARP	1 x 1	3

Spring Wire Retractor

Catalog Number	Length (in.)
MB194R	3 1/2
MB195R	4 1/2

Spring Wire Retractor

Catalog Number	Size	Length (in.)
MB777R	SMALL	5
MB778R	MEDIUM	5 1/8
MB196R	LARGE	9




BV006R-BV007R




BV006R




BV007R




MD471



MB194R-MB195R



MB777R-MB196R





PLASTIC SURGERY HOOKS AND RETRACTORS



BV010R-BV011R

ALM Retractor


Catalog Number	Tip	Prongs	Length (in.)
BV010R	SHARP	4 x 4	2 3/4
BV011R	SHARP	4 x 4	4



BV010R




BV011R




BV028R

ALLPORT Retractor

Catalog Number	Tip	Prongs	Length (in.)
BV028R	SHARP	4 x 4	3 3/4



BV028R








PLASTIC SURGERY ALUMINUM HAND

Aluminum Hand, Blue Silicone Coated

Catalog Number	Size (in.)
F0917A	10 3/4 x 7 1/2
F0918A	16 1/8 x 11



F0917A




F0918A



PLASTIC SURGERY KNIVES

Areola Cookie Cutter

Catalog Number	Diameter (mm)	Blade
AD123R	33.0	TOOTHED
AD125R	38.0	TOOTHED
AD127R	42.0	TOOTHED
AD129R	48.0	TOOTHED



SCHINK Knife, Adjustable From 0.1 To 2.0MM

Catalog Number	Length (in.)
BA706R	11
BA708R	3/4


BA706R Includes:
1 Spare Blade BA708R and
1 Skin Straightening Plate BA709R



PLASTIC SURGERY



PLASTIC SURGERY KNIVES



BA709R

Skin Straightening Plate

Catalog Number	
	BA709R



MD025

GOULIAN Skin Graft Knife Set

Catalog Number	Length (in.)
MD025	7 1/4
BA711	BOX/50

MD025 includes 3 grafting templates:
.20, .25, and .30 mm thick and blades



BA711R



MG575R

TEBBETTS Precision Knife Handle

Catalog Number	Blade Size (mm)	Length (in.)
MG575R	10 - 15	5 3/4



MG575R



PLASTIC SURGERY

MILLIMETERS



PLASTIC SURGERY KNIVES

Skin Graft Knife

Catalog Number	Item	Length (in.)
BA718R	SPARE BLADE, STERILE, 10/PKG	-
BA719R	KNIFE, W/ O BLADE	11 3/4



BA718R-BA719R

KEYES Dermal Punch Set

Catalog Number	OL670R
----------------	--------

OL670R includes:
5 interchangeable tips, handle and rack



OL670R

FOX Dermal Curettes

Catalog Number	Diameter (mm)	Length (in.)
MG150R	2	5 5/8
MG151R	3	5 5/8
MG152R	4	5 5/8
MG153R	5	5 5/8
MG154R	6	5 5/8



MG150R-MG154R



MG150R



MG151R



MG152R



MG153R



MG154R



PLASTIC SURGERY



PLASTIC SURGERY SUCTION INSTRUMENTS

FRAZIER Suction Tube


Catalog Number	Tip Size (FR)	Description	Working Length (in.)
MD124R	7		4 1/4
MD124300	7	SUCTION STYLET W/ OLIVE TIP FOR MD124R	
MD125R	7		4 1/4
MD125300	8	SUCTION STYLET W/ OLIVE TIP FOR MD125R	
MA501R	9		4 1/4
MA501300		SUCTION STYLET W/ OLIVE TIP FOR MA501R	
MD127R	10		4 1/4
MD127300	10	SUCTION STYLET W/ OLIVE TIP FOR MD127R	
MA499R	11		4 1/4
MA499300	11	SUCTION STYLET W/ OLIVE TIP FOR MA499R	
MD129R	12		4 1/4
MD129300	12	SUCTION STYLET W/ OLIVE TIP FOR MD129R	

ANDREWS Pynchon Suction Tube

Catalog Number	Description	Length (in.)
MD613	-	9
MD613001	TIP ONLY F/ MD613	-

GORNEY Suction Dissector, Elevator

Catalog Number	Tip Size (mm)	Length (in.)
OL168R	6-10	7
MA142R	3	8 1/2





PLASTIC SURGERY BONE INSTRUMENTS



D0711R

LUCAS Chisel, Angled

Catalog Number	Tip Size (mm)	Length (in.)
D0711R	4	6 1/4



D0711R



FL120R-FL128R

Chisel, Angled

Catalog Number	Tip Size (mm)	Length (in.)
FL120R	4	5 1/2
FL122R	6	5 1/2
FL124R	8	5 1/2
FL126R	10	5 1/2
FL128R	12	5 1/2



FL120R-FL128R



FL173R-FL176R

PARTSCH Chisel

Catalog Number	Tip Size (mm)	Length (in.)
FL173R	3	5 1/2
FL174R	4	5 1/2
FL175R	5	5 1/2
FL176R	6	5 1/2
FL183R	3	6 3/4
FL184R	4	6 3/4
FL185R	5	6 3/4
FL186R	6	6 3/4
FL187R	7	6 3/4



FL173R-FL176R



FL183R-FL187R



FL183R-FL187R





PLASTIC SURGERY BONE INSTRUMENTS



LUCAS Gouge

Catalog Number	Tip Size (mm)	Length (in.)
D0715R	3	6 1/4

Gouge

Catalog Number	Tip Size (mm)	Length (in.)
FL130R	4	5 1/2
FL132R	6	5 1/2
FL134R	8	5 1/2
FL136R	10	5 1/2

MANNERFELT Gouge

Catalog Number	Tip Size (mm)	Length (in.)
FL153R	3	6 1/8

PARTSCH Gouge

Catalog Number	Tip Size (mm)	Length (in.)
FL213R	3	5 1/2
FL214R	4	5 1/2
FL215R	5	5 1/2
FL216R	6	5 1/2
FL223R	3	6 3/4
FL224R	4	6 3/4
FL225R	5	6 3/4
FL226R	6	6 3/4
FL227R	7	6 3/4

ALEXANDER Gouge

Catalog Number	Tip Size (mm)	Length (in.)
FL343R	7	6 3/4
FL345R	9	6 3/4
FL347R	11	6 3/4



PLASTIC SURGERY BONE INSTRUMENTS



OL307R



OL307R

COTTLE Fishtail Chisel, W/ Protected Ends

Catalog Number	Tip Size (mm)	Length (in.)
OL307R	8	7



FL110R-FL118R



FL110R-FL118R

STRAND Osteotome

Catalog Number	Tip Size (mm)	Blade	Length (in.)
FL110R	4	STR	5 1/2
FL112R	6	STR	5 1/2
FL114R	8	STR	5 1/2
FL116R	10	STR	5 1/2
FL118R	12	STR	5 1/2



FL144R-FL148R



FL144R-FL148R

MANNERFELT Osteotome

Catalog Number	Tip Size (mm)	Length (in.)
FL144R	4	6 1/8
FL146R	6	6 1/8
FL148R	8	6 1/8

Ø 22,5 mm



DF112R

Mallet, Lead Filled, 80gr

Catalog Number	Diameter (mm)	Length (in.)
DF112R	22	6 1/2





PLASTIC SURGERY BONE INSTRUMENTS



F0900R-F0901R

1/2



F0901R

Finger Ring Saw

Catalog Number	Item	Length (in.)
F0900R	COMPLETE	6
F0901R	REPLACEMENT BLADE	-



F0903R

1/2



F0903R

Nail Extracting Forceps

Catalog Number	Length (in.)
F0903R	5 1/4



PLASTIC SURGERY

MILLIMETERS



PLASTIC SURGERY BONE INSTRUMENTS



D0810R

OBWEGESER Awl F/ Lower Jaw

Catalog Number	Length (in.)
D0810R	5 1/2



D0810R



D0812R

OBWEGESER Suture Wire Awl

Catalog Number	Length (in.)
D0812R	7 3/4
D0811R	9



D0812R



D0811R



D0811R



PLASTIC SURGERY

INCHES



PLASTIC SURGERY BONE INSTRUMENTS

Raspatory W/ PEEK Handle








Catalog Number	Tip	Tip Size (mm)	Blade	Length (in.)
FK318R	FLAT	6	STR	7 1/2
FK317R	FLAT	3	CVD	7 1/4
FK319R	FLAT	6	CVD	7 1/2
FK321R	RND	6	CVD	7 1/2

MOLT Raspatory and Elevator, Double Ended

Catalog Number	Length (in.)
DO543R	7

SELDIN Raspatory, Double Ended

Catalog Number	Size	Length (in.)
DO552R	SMALL	7 3/4
DO553R	LARGE	7 3/4




PLASTIC SURGERY BONE INSTRUMENTS






MILLER Bone Rasp

Catalog Number	Tip Size (mm)	Length (in.)
D0773R	6.0 / 5.0	7



MANNERFELT Modified Elevator

Catalog Number	Tip	Tip Size (mm)	Length (in.)
FK530R	SPOON	4	6 1/8
FK532R	SPOON	7	6 1/8
FK534R	SPOON	9	6 1/8
FK541R	CVD	3	6 1/8
FK543R	CVD	6	6 1/8
FK552R	SLIGHT-CVD, ROUND	8	6 1/8



MARBURG Needle Nose Pliers

Catalog Number	Tip	Length (in.)
DP031R	SERR	5 1/2





PLASTIC SURGERY BONE INSTRUMENTS

Universal Wire Cutter

Catalog Number	For Wire Size	Length (in.)
DP501R	SOFT, UP TO 1.5MM AND HARD, UP TO 1.0MM	6
DP505R	SOFT, UP TO 1.5MM AND HARD, UP TO 1.2MM	6 1/4

Wire Cutter, W/ Specially Hardened Cutting Edges

Catalog Number	For Wire Size	Blade End	Length (in.)
DP520R	SOFT, UP TO 1.5MM AND	X-HARD EDGE	5 3/4
DP530R	SOFT, UP TO 1.5MM AND	X-HARD EDGE	5 1/2







PLASTIC SURGERY BONE INSTRUMENTS







DUROTIP™ TC Wire Cutting Scissors

Catalog Number	For Wire Size	Blade	Length (in.)
DP552R	SOFT, UP TO 1.2MM AND HARD, UP TO 0.8 MM	1-SERR	4 3/4



BEEBE Crown Scissors

Catalog Number	Tip	Blade	Length (in.)
DP560R	SHARP	STR	4 1/4
DP561R	SHARP	CVD	4 1/4
DP563R	BROAD	CVD	4 1/4
DP569R	ONE-SERR	CVD	4 1/4





PLASTIC SURGERY MICRO VESSEL CLIPS AND CLAMPS

- The arterial clips have a slightly inwards curved jaw tip to prevent slipping, see Fig.1
- The venous clips have a parallel closing jaw, see Fig. 2



Fig. 1








Fig. 2

PLASTIC SURGERY



ACLAND Micro Vessel Clips

Catalog Number	Jaw Length (mm)	Size	Description	Vessel Size (mm)
FD870R	8	B-1	Artery	0.3-1.0
FD871R	8	B-1	Vein	0.3-1.0
FD872R	11	B-2	Artery	0.6-1.4
FD873R	11	B-2	Vein	0.6-1.4
FD874R	16	B-3	Artery	1.0-2.2
FD875R	16	B-3	Vein	1.0-2.2
FD876R	24	HD-S	Artery	1.5-3.3
FD877R	24	HD-S	Vein	1.5-3.3
FD878R	35	RD-S	Artery	2.0-5.0
FD879R	35	RD-S	Vein	2.0-5.0
FD880R	8	B-1	Artery	0.3-1.0
FD881R	8	B-1	Vein	0.3-1.0
FD882R	11	B-2	Artery	0.6-1.4
FD883R	11	B-2	Vein	0.6-1.4
FD884R	16	B-3	Artery	1.0-2.2
FD885R	16	B-3	Vein	1.0-2.2
FD886R	24	HD-D	Artery	1.5-3.3
FD887R	24	HD-D	Vein	1.5-3.3
FD888R	35	RD-D	Artery	2.0-5.0
FD889R	35	RD-D	Vein	2.0-5.0
FD890R	8	B-1	Artery	0.3-1.0
FD891R	8	B-1	Vein	0.3-1.0
FD892R	11	B-2	Artery	0.6-1.4
FD893R	11	B-2	Vein	0.6-1.4
FD894R	16	B-3	Artery	1.0-2.2
FD895R	16	B-3	Vein	1.0-2.2
FD896R	5 1/2	B-1, B-2	Clip Applying Forceps	-
FD897R	5 1/2	B-3, HD, RD	Clip Applying Forceps	-

All items are sold individually



PLASTIC SURGERY MICRO VESSEL CLIPS AND CLAMPS




FD560R-FD562R

BIEMER Vessel Clip


Catalog Number	Jaw Length (mm)	Jaw Opening (mm)	Closing Force
FD560R	6	4	30 - 40
FD561R	9	5	30 - 40
FD562R	9	5	20 - 25



FD560R



FD561R



FD562R



FD568R

**BIEMER Applying Forceps,
F/ FD560R-FD562R, FD611R-FD615R**

Catalog Number	Description	Length (in.)
FD568R	W/ O LOCK	5 1/2
FD569R	W/ LOCK	5 1/2




FD569R



FD571R

BIEMER Micro Approximator Only

Catalog Number	Jaw Length (mm)	Jaw Opening (mm)	Closing Force (gr)
FD571R	6	2	30 - 40



FD571R



FD572R

BIEMER Applying Forceps, F/ FD571R

Catalog Number	Length (in.)
FD572R	5 3/4



FD572R





PLASTIC SURGERY MICRO VESSEL CLIPS AND CLAMPS








Micro Approximator

Catalog Number	Description
FD575R	COMPLETE WITH KEY
FD576R	KEY ALONE

BIEMER-MÜLLER Approximator

Catalog Number	Blade	Jaw Length (mm)	Jaw Opening (mm)	Closing Force (gr)
FD613R	STR, DEL	9	4	25 - 30
FD615R	STR	4,5	3	25 - 30
FD611R	STR	9	3	25 - 30
FD614RANG, DEL		9	3	25 - 30
FD612R	ANG	9	3	25 - 30



BIEMER Approximator, F/ Nerves

Catalog Number	Description
FD576R	KEY ALONE
FD635R	APPROXIMATOR COMPLETE WITH KEY AND NEEDLES
FD636R	SPARE NEEDLES, PR