



Wholesale Supplier
Hospital Furniture Deluxe Quality
A to Z Hospital Equipment and
Surgical Products B to C



8275140877

9822688071

info@amramedtech.com



**Sr No.29/1/2, Adarsh Nagar,
Uruli Devachi, Hadapsar, Pune**
Branches :- All Maharashtra




amramedtech.com



Delivering Quality Hospital Equipments

HOSPITAL FURNITURE





AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C




Certified Register

FIVE FUNCTIONS ELECTRIC ICU BED

Description


- Frame work made of strong Rectangular CRCA pipe.
- Four section top made of perforated ABS structure
- Back Rest, Knee Rest, TB, RTB & Hi-Low Position operated with hand remote & nurse control panel.
- Heavy duty ABS curved Head & Foot Panel with attached buffer.
- Heavy duty ABS full length Tuck away safety side rails with angle indicator on both side.
- Nurse Control panel which operate all functions with electronic CPR & provision for locking positions of bed.
- Battery backup option for uncertain power cut.
- Bed mounted on 125mm Dia double tier Castors (2 with breaks)
- S.S telescopic IV Rod with 4 foldable hooks.
- Provision for Mattress guard & urine bag holder on both side.
- 4" Thick Medical Grade Mattress of 40 Density.
- Fibre base cover to protect from dust .
- Overall Size : 2200L X 1000W X 600-800H mm.
- Finish :- Epoxy powder coated .



FIVE FUNCTIONS ELECTRIC ICU BED

Description

- Frame work made of strong Rectangular CRCA pipe.
- Four section top made of CRCA sheet.
- Back Rest , Knee Rest , TB , RTB & Hi-Low Position operated with hand remote only.
- ABS Head & Foot Panel.
- ABS Tuck away safety side rails.
- Bed mounted on 125mm Dia double tier Castors (2 with breaks)
- S.S telescopic IV Rod with 4 foldable hooks.
- Provision for Mattress guard & urine bag holder on both side.
- 4" Thick Medical Grade Mattress of 40 Density.
- Fibre base cover to protect from dust .
- Overall Size : 2200L X 1000W X 600-800H mm.
- Finish :- Epoxy powder coated .



FIVE FUNCTIONS ELECTRIC ICU BED

Description

- Frame work made of strong Rectangular CRCA pipe.
- Four section top made of CRCA sheet.
- Back Rest , Knee Rest , TB , RTB & Hi-Low Position operated with hand remote only.
- Heavy duty ABS Head & Foot Panel.
- Heavy duty ABS Tuck away safety side rails with angle indicator.
- Bed mounted on 125mm Dia double tier Castors (2 with breaks)
- S.S telescopic IV Rod with 4 foldable hooks.
- Provision for Mattress guard & urine bag holder on both side.
- 4" Thick Medical Grade Mattress of 40 Density.
- Overall Size:-2200L X 1000W X 600-800H mm.
- Finish : Epoxy powder coated.



Optional :



Nurse Control
System(ACP)



Central Locking
Castors.



ICU BED ELECTRIC

Description



- Advance linear actuator system [04Pc Control box, 01Pc Handset, 01 Pc. Nurse Control System (ACP)]
- All positions can be controlled by a wired remote hand set placed on the side rail or by ACP (attendant control panel) at the foot end
- Easily detachable polymer head and foot bows with buffer
- Extra long 4pcs. Under side Occidental style side rails with inclination degree indicators
- Four section ABS top
- Framework made up of M.S tubes
- Bed Mounted on 5" Central Locking ABS Castors, Durable and quiet
- S.S. Telescopic I.V. Rod with 4 Location
- Holder for Urine Bag
- Sitting position adjustable
- Finish : Pre-treated & Epoxy Powder Coated
- Overall Approx. Size: 2060L x 900W x 470-810H mm
- C4 Section Mattress (MP-2001C)
- LED Lamp as reading lamp or under Bed lamp
- Positions Obtained by Hand Remote:
 - Hi-Low · Backrest elevation
 - Knee rest elevation
 - Trendelenburg / Reverse Trendelenburg
- Positions obtained by membrane
 - Hi-Low · Backrest elevation · Knee rest elevation
 - Trendelenburg/Reverse Trendelenburg
 - Electric CPR for zero positioning of bed
 - Fowler position
- Positions obtained/Functions on Nurse Control System (ACP)
 - Hi-Low with Locking function for handset
 - Backrest elevation with Locking function for handset
 - Kneerest elevation with locking function for handset
 - Trendelenburg / Reverse Trendelenburg with locking function for handset
 - One button Electric CPR Position and manual CPR Function for any Emergency case



Nurse Control




Central Locking
Castors.



ICU BED ELECTRIC

Description



- Advance linear actuator system [04Pc Control box, 01Pc Handset, 01 Pc. Nurse Control System (ACP)]
- All positions can be controlled by a wired hand set
- Extra long 4pcs. Under side Occidental style side rails with inclination degree indicators
- Heavy ABS Head and Foot Panel with locking and xorner buffer ABS top.
- Framework made up of M.S tubes
- Bed Mounted on 5" Central Locking ABS Castors, Durable and quiet
- S.S. Telescopic I.V. Rod with 4 Location
- Holder for Urine Bag
- Sitting position adjustable
- Finish : Pre-treated & Epoxy Powder Coated
- Overall Approx. Size: 2060L x 900W x 470-810H mm
- C4 Section Mattress (MP-2001C)
- LED Lamp as reading lamp or under Bed lamp.
- Positions Obtained by Hand Remote:
 - Hi-Low · Backrest elevation
 - Knee rest elevation
 - Trendelenburg / Reverse Trendelenburg
- Positions obtained/Functions on Control System (ACP)
 - Hi-Low with Locking function for handset.
 - Backrest elevation with Locking function for handset.
 - Knee rest elevation with Locking function for handset.
 - Trendelenburg/Reverse Trendelenburg
 - Trendelenburg / Reverse Trendelenburg with locking function for handset
 - One button Electric CPR Position and manual CPR Function for any Emergency case



Nurse Control



Central Locking Castors.



CHAIR POSITION ICU BED

Description

- Advance linear actuator system [6 Pcs Motor, 01 Pc Control box, 01Pc Handset, 01 Pc. Nurse Control System (ACP)]
- Extra Long four pieces under Side occidental style side rails with inclination degree indicator
- Heavy ABS Head and Foot Panel with locking and membrane system
- Four Section top
- Back section X-Ray Translucent
- Framework made of M.S Tubes
- Bed Mounted on 6" Central Locking ABS Castors, Durable and quiet
- S.S. Telescopic I.V. Rod with 4 Location
- Holder for Urine Bag
- Sitting position adjustable
- Leg Section Extendable
- Finish: Pre-treated & Epoxy Powder Coated
- Overall Approx Size: 1970-2120L x 900W X 460-740H MM
- C4 Section Mattress
- LED Lamp as reading lamp or under Bed lamp Positions
- Obtained by Hand Remote:
 - Hi-Low
 - Backrest elevation
 - Knee rest elevation
 - Trendelenburg/ Reverse Trendelenburg
- Positions obtained by membrane
 - Hi-Low
 - Backrest elevation
 - Knee rest elevation
 - Trendelenburg/ Reverse Trendelenburg Electric CPR for zero positioning of bed Fowler position
- Positions obtained/Functions on Nurse Control System (ACP)
 - Hi-Low with Locking function for handset
 - Backrest elevation with locking function for handset
 - Kneerest elevation with locking function for handset
 - Trendelenburg/ Reverse Trendelenburg with locking function for handset
 - One button Electric CPR Position and manual CPR Function for any Emergency case
 - One Touch Chair Position



Nurse Control



Central Locking



ICU BED ELECTRIC



Description

- Frame work made of strong Rectangular CRCA pipe.
- Four section top made of CRCA sheet.
- Back Rest , Knee Rest , TB , RTB & Hi-Low Position operated with hand remote only.
- Heavy duty ABS Head & Foot Panel.
- Heavy duty ABS Tuck away safety side rails with angle indicator.
- Bed mounted on 125mm Dia double tier Castors (2 with breaks)
- S.S telescopic IV Rod with 4 foldable hooks.
- Provision for Mattress guard & urine bag holder on both side.
- 4" Thick Medical Grade Mattress of 40 Density.
- Overall Size:-2200L X 1000W X 600-800H mm.
- Finish : Epoxy powder coated.



Description


- Advance linear actuator system [04Pcs Motor, 01Pc Control box, 01Pc Handset]
- Collapsible side railing
- Heavy ABS Head and Foot Panel with buffer
- 4 section top made up of oval perforated CRCA M.S. sheet.
- Framework made up of M.S. tubes
- Bed mounted on 5" ABS castors durable and quiet
- S.S. Telescopic I.V. Rod with 4 location
- Holder for Urine Bag
- Sitting position adjustable
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2090L x 900W x 550-800H mm
- Positions obtained by Hand Remote:
 - • Hi-Low
 - Backrest elevation
 - Kneerest elevation
 - Trendelenburg / Reverse Trendelenburg
- • **Optional :**
 - • C4 Section Mattress
 - S.S. Telescopic I.V. Rod with 4 location



FIVE FUNCTIONS MANUAL ICU BED

Description

- Frame work made of strong Rectangular CRCA pipe.
- Four section top made of CRCA sheet.
- Back Rest , Knee Rest, TB, RTB & Hi-Low Position operated with S.S Crank Handles.
- ABS Moulded Head & Foot Panel.
- Collapsible safety side railing
- S.S telescopic IV Rod with 4 foldable hooks.
- Provision for Mattress guard & urine bag holder on both side.
- Bed mounted on 125mm Dia double tier Castors (2 with breaks)
- 4" Thick Medical Grade Mattress of 40 Density.
- Overall Size : 2200L X 1000W X 600-800H mm.
- Finish : Epoxy powder coated




ICU BED MECHANICAL

Description

- Four pieces individual side-rails, ABS moulded.
- Heavy ABS Head and Foot Panel with buffer
- 4 section top made up of oval perforated CRCA M.S. sheet.
- Framework made up of M.S. tubes
- Bed mounted on 5" ABS castors 2 with brakes.
- Holder for Urine Bag
- Sitting position adjustable
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2090L x 900W x 550-750H mm
- Positions obtained by Smooth Crank System with foldable handle:
 - Hi-Low
 - Backrest elevation
 - Kneerest elevation
 - Trendelenburg / Reverse Trendelenburg
- **Optional :**
 - C4 Section Mattress
 - S.S. Telescopic I.V. Rod with 4 location.
 - Stainless Steel / Mild Steel Railing
 - Stainless Steel Head & Foot Bows



S. S. Head &
Foot Panel




ICU BED MECHANICAL

Description

- Two pieces individual detachable side-rails.
- Collapsible side railing - ABS Head and Foot Panel.
- 4 section top made up of oval perforated CRCA M.S. sheet.
- Framework made up of M.S. tubes
- Bed mounted on 5" castors 2 with brakes.
- Holder for Urine Bag
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2000L x 900W x 600H mm
- Positions obtained by Smooth Crank System with foldable handle :
- Backrest elevation
- Kneerest elevation
- Trendelenburg / Reverse Trendelenburg
- **Optional :**
- C4 Section Mattress
- S.S. Telescopic I.V. Rod with 4 location.
- Stainless Steel Head & Foot Bows



S. S. Head & Foot Panel




FOWLER BED MECHANICAL

Description

- Two pieces individual detachable collapsible side-rails.
- ABS moulded Head and Foot Panel.
- 4 section top made up of oval perforated CRCA M.S. sheet.
- Framework made up of M.S. tubes
- Bed mounted on 4" castors 2 with brakes.
- Four location for I.V. Rod.
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2000L x 900W x 600H mm
- Positions obtained by Smooth Crank System with foldable handle :
- Backrest elevation
- Kneerest elevation
- **Optional :**
- 4 Section Mattress
- S.S. Telescopic I.V. Rod
- Four pieces individual side rails, ABS moulded.
- Stainless Steel Head & Foot Panels



S. S. Head & Foot Panel




FOWLER BED ELECTRIC

Description

- Advance linear actuator system [04Pcs Motor, 01Pc Control box, 01Pc Handset]
- Four pieces individual side-rails, ABS moulded.
- ABS moulded Head and Foot Panel with locking & corner buffer
- 4 section top made up of oval perforated CRCA M.S. sheet.
- Framework made up of M.S. tubes
- Bed mounted on 4" castors 2 with Coated.
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2000L x 900W x 600H mm
- Positions obtained by Hand Remote:
- Backrest elevation
- Kneerest elevation
- **Optional :**
- C4 Section Mattress
- S.S. Telescopic I.V. Rod with 4 location.
- 2 pieces individual detachable collapsible side rails.



S. S. Head & Foot Panel




HOSPITAL FOWLER BED

Description

- Mild steel tubular Head & Foot Bows of unequal height with vertical tube support.
- 4 section top made up of oval perforated CRCA M.S. sheet.
- Framework made up of M.S. tubes
- Bed mounted on 4" castors 2 with brakes.
- Four location for I.V. Rod.
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2000L x 910W x 600H mm
- Positions obtained by Smooth Crank System with foldable handle :
- Backrest elevation
- Kneerest elevation
- **Optional :**
- 4 Section Mattress
- S.S. Telescopic I.V. Rod
- S.S. Collapsible Railing




Also available in wooden color



HOSPITAL SEMI FOWLER BED

Description

- ABS moulded Head & Foot Panel with locking & corner buffer.
- 2 section top made up of oval perforated CRCA M.S. sheet.
- Framework made up of M.S. tubes
- Legs fitted with PVC stumps.
- Four location for I.V. Rod.
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2000L x 900W x 600H mm
- Positions obtained by Smooth Crank System with foldable handle :
- Backrest elevation
- **Optional :**
- 2 Section Mattress
- S.S. Telescopic I.V. Rod
- S.S. Collapsible Railing
- Advance linear actuator system (01 Pc Motor, 01 Pc control box, 01 Handset)
- M.S or S.S. Head & Foot Panels.



Description


- Mild steel tubular Head & Foot Bows of unequal height with vertical tube support.
- Top made up of oval perforated CRCA M.S. sheet.
- Framework made up of M.S. tubes · Legs fitted with PVC stumps.
- Four location for I.V. Rod.
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2000L x 910W x 600H mm
- **Optional :**
- Plain Mattress
- S.S. Telescopic I.V. Rod
- S.S. Collapsible Railing



ATTENDED BED

Description

- Overall approx. size: 1800L x 600W x 460H mm
- Framework made up of Rectangular / Square M.S. tubes
- Also available in 5cm from upholstery top.
- Legs fitted with PVC stumps.
- Finish: Pre-treated & Epoxy Powder Coated
- **Optional :**
- Plain Mattress
- 5 cm foam padded fixed mattress covered with rexine.



PATIENT'S OPD COUCH/CHAIR

Description

Blood Transfusion Chair

- Dimension : 890 x 600 x 300 mm
- Net Wt. 26 kg
- Made of MS Powder Coated with Back Adjustable



Description

Attendant Chair Cum Bed

- Overall approx size : 1880mmL x 500mmW x 400mH.
- M.S. tubular frame work mounted on PVC Stumps
- Epoxy Powder Coated.



PEDIATRIC BED

Description

- Tubular Head and Foot Bows of unequal height with vertical tube support
- Frame work made up of M.S. tubes
- Top made up of oval perforated CRCA M.S. sheet
- Legs fitted with PVC stumps - Four location for I.V. Rod
- Full length drop side rails on both sides
- Finish : Pre-treated & Epoxy Powder Coated
- Overall approx. size: 1460L x 830W x 600H mm
- Bed frame approx size: 1370L x 760W mm
- Also available in electric

Options

- Plain Mattress
- S.S. Telescopic I.V. Rod
- Full length S.S. Side Rails



MOTORIZED BACKREST/ RECLINER BED

Description

- Headrest adjustable upto 75°
- Can be used on the existing bed without any modification just plug & play
- Size with mattress 1820 x 900 x 200 mm
- Recliner Frame can be attach over normal home bed.

MATTRESS FOR HOSPITAL BEDS

Mattress for Plain Beds



Two Section Mattress for
Semi-Fowler Beds



Four Section Mattress
for Fowler / I.C.U. Beds



8275140877



9822688071



info@amramedtech.com



amramedtech.com



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C

Certified Register



WARD EQUIPMENTS

Patient's Medicine Cupboard (DLX)

Description

- Overall size: 400L x 400W x 810H mm
- CRCA sheet constructions having one box and one drawer
- Jointless S.S. top with three sides raise edges
- Finish : Pre-treated and Epoxy Powder Coated

Optional

- Also available in Complete S.S.



Patient's Medicine Cupboard (STD)

Description

- Overall size: 400L x 400W x 810H mm
- Jointless S.S. top with three sides raise edges
- One box with door made of M.S. sheet
- Four MS legs fitted with PVC stumps
- Finish : Pre-treated and Epoxy Powder Coated

Optional

- Also available in Complete S.S.

Patient's Medicine Cupboard (ABS)

Description

- Overall size: 450L x 450W x 760H mm
- One drawer
- One towel rail
- Four 50mm Castors



8275140877



9822688071



info@amramedtech.com



amramedtech.com

Pune Maharashtra India



WARD EQUIPMENTS

Over Bed Table (Adjustable by Pneumatic Gas Spring)

Description

- Rectangular frame work mounted on 5cms castors
- 30"x 18" laminated /Membrane Pressed top adjustable by Pneumatic gas spring.
- Finish : Pre-treated and Epoxy Powder Coated.



Over Bed Table Mayo's Type (Adjustable Height with Gear Handle)

Description

- M.S. Rectangular frame work mounted on 5 cms castors.
- 30" x 16" Laminated top adjustable by gear handle.
- Epoxy Powder Coated



Over Bed Table Mayo's Type S.S. (Adjustable by Pneumatic Gas Spring)

Description

- S.S. Rectangular frame work mounted on 5 cms castors.
- 30" x 16" Laminated top adjustable by gear handle.



Over Bed Table (Manual)

Description

- Rectangular frame work mounted on 5cms castors
 - 30" x 16" laminated top adjustable by side knob.
 - Epoxy Powder Coated
- Also available in All Stainless Steel**



O.T. Mayo's Instruments Trolley with S.S. Tray

Description

- Tray S.S. size : 56 L x 40 W cms
- Rectangular CRC tubes mounted on 5 cms castors
- Height adjustable by side knob
- S.S. frame with S.S. top



O.T. Mayo's Instruments Trolley Over the O.T. Table (Mechanically)

Description

- Approx tray S.S. Size: 91 L x 50W cms
- Frame work of S.S. tubes mounted on 10 cms castors
- Height adjustable by gear handle
- S.S. frame with S.S. top



WARD EQUIPMENTS

Wash Basin Stand (Double)

Description

- Five Legs base mounted on 5 cms dia castors
- With two S.S. basin approx. 35 cms
- Vertical S.S. tube made of 25mm x 18G
- Basin holder S.S. tube made of 16mm x 18G
- Finish : Stainless Steel frame with Stainless Steel Basin

Optional

- Also available in Epoxy Powder Coated frame




Wash Basin Stand (Single)

Description

- Five Legs base mounted on 5 cms dia castors
- With two S.S. basin approx. 35 cms
- Vertical S.S. tube made of 25mm x 18G
- Basin holder S.S. tube made of 16mm x 18G
- Finish : Stainless Steel frame with Stainless Steel

Basin Optional


- Also available in Epoxy Powder Coated frame



I.V. Stand

Description

- Five Pronged base fitted with 5 cms Castors
- S. S. adjustable
- I.V. Rod with double / four hook
 - (A) S.S. Pipe with plastic base
 - (B) M.S. Pipe with plastic base
 - (C) S.S. Pipe with S.S. base
 - (D) M.S. Pipe with M.S. base



WARD EQUIPMENTS

Bed Side Screen (4 Panels)

Description

- Overall size: 1680H x 2450L mm
- M.S. tubular frame work made up of 19mm x 18G in four section
- Legs made of 25mm x 18G mounted on 5 cms castors with hooks and spring along-with curtains
- Each span is 610 mm wide
- Finish : Pre-treated and Epoxy Powder Coated



Bed Side Screen (3 Panels)

Description

- Overall size : 1680H x 2450L mm
- M.S. tubular frame work made of 19mm x 18G in three section
- Legs made up of 25mm x 18G mounted on 5 cms castors
- With hooks and spring along-with curtains
- Each span is 610mm wide
- Finish : Pre-treated and Epoxy Powder Coated



Kick Bowl / Bucket

Description

- Approx. S.S. Bowl of 350 mm dia/SS Bucket
- Frame work made up of M.S. tubes
- Mounted on 50 mm castors
- Finish : Stainless Steel

Also available in M.S. frame work and S.S.
Basin/Bucket



PATIENT TRANSFER TROLLEYS

Emergency & Recovery Trolley (Hydraulic)

Description

- Two section top height adjusted by foot operated imported hydraulic pump from 660 mm to 910 mm
- Gas spring assisted Trendelenburg/Reverse Trendelenburg positions
- Complete with corner buffers, synthetic rubber covered handles, sliding X-ray cassette, accessories tray, oxygen cylinder holder
- S.S. Telescopic IV rod and swing away S.S. side rails
- Finish : Pre-treated and Epoxy Powder Coated
- Four swivel castors, 5" dia, two with total lock
- Xray translucent removable stretcher, backrest raised on ratchet
- Overall approx size: 210L x 65W x 65-95H cms

Also available in complete S.S.



Stretcher on Trolley

Description

- Frame work made of 1.25" x 16G vertical and 1" x 18G horizontal CRC tubes
- Trolley mounted on 15 cms. dia castors, 2 with brakes
- Removable stretchers top made of M.S. Sheet with S.S. handle
- Overall size : 213L x 56W x 81H cms
- Finish : Pre-treated and Epoxy Powder Coated

Optional

- Mattress
- Side railings M.S./S.S.

Also available in complete stainless steel



Stretcher Trolley with Mattress

Description

- Overall approx size : 180L x 56W x 81H cms
- Frame work made of tubular M.S. pipe with S.S. handle
- Swing away type safety side railing
- Adjustable I.V. Rod, Utility tray, Cylinder cage foam padded top
- Trolley mounted on 15 cms Dia castors- 2 with brakes
- Finish: Pre-treated & Epoxy Powder Coated

Also available in complete Stainless Steel



PATIENT TRANSFER TROLLEYS

Emergency & Recovery Trolley (Hydraulic)

Description

- Indian Recovery Trolley Highly customizable for best patient care
- Excellent maneuverability to reduce the risks related with patient handling
- Easy to operate and maintain, also prevent injury to patient
- Easy to access and allow unlimited control of all the function

FACILITIES

- X-ray cassette can be slide under the translucent X-ray s platform
- Recovery trolley frame allows the undercarriage and wide range of height adjustment make the handling with C arm easier The Recovery Trolley allows treatment of patient easily and quickly
- The diagnostics is facilitated by possibility of X-ray examination and C-arm Systems directly on the stretcher reducing time for diagnostics

EASY TO MOVE

- Recovery Trolley is easy Maneuver & 360° Rotation on 5th wheel (Optional). Adjustment of height is operated by foot controls through hydraulic pumps. Trendelenburg & reverse position by gas spring.



Standard Accessories



Oxygen Cylinder Holder



Side Rails



Monitor/Ventilator Tray



I.V. Rod

Emergency & Recovery Trolley (Electric)

Description

- Frame work made of Rectangular M.S. Tube.
- 4 Section Top made of Perforated M.S. Sheet.
- Back Rest , Knee Rest , TB , RTB & Hi-Low Position operated with cable remote control
- S.S 304 solid tubular handles on head-foot side for easy movement
- S.S Collapsible side rails
- Nursing Control panel which operate all functions with CPR & lock positions of bed .
- Mounted on 125mm Dia Castors (2 with breaks)
- S.S IV Rod of 4 foldable hooks with four locations
- Overall Size:-2000L X 700W X 600-800H mm.
- PU made base cover to protect from dust .
- Finish : Epoxy powder coated .



8275140877



9822688071



info@amramedtech.com



amramedtech.com



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C

Certified Register



PATIENT TRANSFER TROLLEYS

Folding Stretcher

Description

- Aluminium folding lengthwise as well as width wise Fold, provided with bag
- Unfolded size: 221L x 53W x 15H cms
- Folded size: 110L x 18W x 11H cms
- Load bearing capacity 159 Kgs
- Gross wt. 6 Kg, Net wt. 5 Kg



Head
Immobilizer



Stair
Stretcher



Air
Splint

Spine Board

Description



- Product Size : 184L x 45W x 5H cm
- Load Bearing : < 159 kg
- Net Weight : 9
- Gross weight : 12 kg
- Dimension of the packing box : 187L x 48W x 8H cm

Scoop Stretcher



Description

- Product size (LxWxH) : 164x 43x8cm(min.) 201x43x8cm (Max)
- Net weight : 9kg
- The bearing pressure : < 159kg
- Dimension of the packing box : 123x45x11cm MP - 545
- Gross weight : 11k

Ambulance Stretcher

Description

- Utilizes high precision mechanical transmission
- Can be used automatically folded through the left & right hand control handle
- Made of high quality aluminium alloy tubes
- Dimensions :
 - A - High Position: 1900 x 560 x 900mm
 - B - Low Position: 1900 x 560 x 280mm
- Load weight up to 120 Kg
- Product weight: 40 Kg
- Diameter of the castors: 125 mm
- Applicable for various kind of hospitals, emergency centers, ambulances & battle fields



EMERGENCY TROLLEYS

Laparoscopic / Monitor Trolley

Description

- Overall approx. size : 820 mm L x 620 mm D x 1680 H
- Four shelves and one drawers
- Trolley having two handle
- The trolley mounted on 125 mm dia castors, 2 with brakes
- Switch & Socket for medical equipment
- Supplied with S.S. IV Rod and M.S. Oxygen Gas Cylinder Cage
- Finish : Pre-treated and Epoxy Powder Coated
- **Also available in Complete S.S.**



Crash Cart

Description

- Overall approx. size : 960 L x 500 W x 1545H mm
- Frame work made up of S.S. tubes
- Two light weight polystyrene boxes each with three drawers
- Upper drawer with medicine container of different sizes
- Provision to hold oxygen cylinder made up of M.S.
- Massage board for cardiac patient are provided
- Six number colour hand out small bins at eye level
- Total three S.S. shelves with railings to carry medical equipments
- Provided with S.S. I.V. Rod, corner buffers
- The cart mounted on 125 mm dia castors 2 with brakes
- **Also available in Epoxy Powder Coated Framework & S.S. To**



Emergency Medicine Trolley

Description

- Overall approx. size : 750 x 500 W x 810H mm
 - Mild steel frame work with five draws for medicine
 - Upper drawer with medicine container of different sizes
 - Provision to hold oxygen cylinder made up of M.S.
 - S.S. shelves with railings to carry medical equipments
 - Provided with S.S. I.V. Rod, corner buffers
 - The cart mounted on 125 mm dia castors 2 with brakes
- Also available in Complete S.S.



WARD EQUIPMENTS

Medicine Trolley with 4 Drawers

Description

- Overall approx. size: 76L x 46W x 81H cms
- Frame work made up of M.S. tubes mounted on 10 cms castors
- Two S.S. shelves with three side railing on top shelf
- Two drawer under the each shelf
- Finish : Epoxy Powder Coated



Utility Trolley Two Shelves

Description

- M.S. Frame work
- S.S. Shelves
- One M.S. Drawer
- 10 cms Castors
- **Also available in three shelves**

Dressing Trolley

Description

- Approx. size : 30"x18"x32"
- Frame work made up of S.S. tubes
- Two S.S. shelves up with S.S. railings on all four sides
- The trolley fitted S.S. bucket and S.S. bowl
- The trolley mounted on 100 mm dia castors
- S.S. frame, S.S. shelves with S.S. basin & bucket



Also available in Epoxy Powder Coated Frame & S.S Top



Instrument Trolley Two Shelves

Description

- Frame work made up of 1"x 18G vertical M.S. tubes
- Two 22G S.S. shelves with three sides railing
- The trolley mounted on 100 mm dia castors
- Finish : S.S. frame with S.S. shelves
- 30"x18"x32"
- 36"x24"x32"
- 24"x18"x32"

WARD EQUIPMENTS

Solid Linen Trolley With Canvas Bag

Description

- Frame work made up of CRC tube mounted on 10 cms castors
- Superior quality canvas bag
- Finish: Pre-treated and Epoxy Powder Coated

Also available in Complete S.S



Solid Linen Trolley With Canvas Bag

Description

- Push handle at one end
- Frame work made up of CRC tube on 10 cms castors
- Superior quality plastic bucket
- Pre-treated & Epoxy Powder Coated
- Also available in S.S. frame
- Complete in S.S. frame
- **Also available in Complete S.S.**



Foot Step Single Step

Description

- Overall approx. size : 505L x 305W x 230H mm
- Step made up of CRCA Sheet fitted with Aluminium Tread flats by pop rivets
- Frame made up of 1" x 18G CRCA tube fitted with PVC stumps
- Finish: Pre treated and Epoxy Powder Coated

Also available in Complete S.S




Foot Step Double Step

Description

- Overall approx. step size : 505 mm L x 305 mm W
- First step height 230 mm and second step size 450 mm
- Step made up of CRCA sheet fitted with Aluminium Tread flats by pop rivets
- Frame made up of 1"x18G CRCA tubes fitted with PVC stumps
- Finish: Pre-treated and Epoxy Powder Coated

Also available in Complete S.S



BIO-MEDICAL WASTE DUSTBINS

Waste Dustbins

Description

- Size : 60/30/20/10 Liters
- Material : PP/ CP Injection Moulded
- Available in colour coding options of Red,
Blue, Black, Yellow & Green



Waste Bin (Push Can)

Description

- Size : 20x40cm
- Size : 26x47cm
- Size : 31x57cm



Dustbin Push
Type S.S.

Waste Bin (Paddle Type)

Description

- Size : 5 Liters
- Size : 7.5 Liters
- Size : 12 Liters
- Size : 20 Liters
- Size : 35 Liters



Dustbin Pedestal S.S. With Plastic
Bucket

BIO-MEDICAL WASTE SEGREGATION TROLLEY

Description

- Size/Capacity : 5ltr — 60ltr
- Waste segregation Trolley with one, Two,
Three & Four Compartment
- Trolley available in M.S. and S.S.
- Foot operated system for touch less
opening off lid
- Bins can be interchanged as and when required
- Color available : Red, Blue, Yellow, Black & Green



WARD EQUIPMENTS

Trolley For Oxygen Cylinder

Description

- Frame work made up of CRC tube mounted on 10 cms/15cms castors
- Finish: Epoxy Powder Coated
- (A) For D type cylinder
- (B) For B type cylinder



Patient Revolving Stool S.S. Top

Description

- Four leg base fitted with PVC stumps
- Height adjustment by machine screw from 18" x 27"
- Rest ring support made up of M.S. tube
- S.S. top supported with M.S. sheet top
- Finish : Pre-treated and Epoxy Powder Coated

Also available in S.S. Frame work



Patient's Stool

Description

- Approx. size: 12"x 12"x 20"H
- Top made up of CRCA sheet double bent
- Legs made up of 25 mm x 18 G and horizontal support of 19 mm x 18 G CRCA tube fitted with PVC stumps
- Finish : Pre-treated and Epoxy Powder Coated



Surgeon Chair

Description

- Cushioned seat and back
- Pneumatic height adjustable
- Height adjustable and swinging back rest
- Legs chrome plated base with foot ring
- Base fitted with 50 mm dia castors





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C

Certified Register



Baby Cradles

Description

- Overall approx. size: 36" x 15" x 39"
- Stand frame work made up of 25 mm x 1.22 mm M.S. tubular pipe
- Legs fitted with PVC stumps
- Wire mesh hanging cradle made up of 19 mm x 1.22 mm
- CRC tube and wire mesh on 5 mm thick
- Finish: Pre-treated and Epoxy Powder Coated



Medical Clip Board

Description

- Dimension 340 x 230 mm
- A BS Plastic Material
- Colour Available : Blue, Grey & Pink



Baby Bassinet

Description

- Dimension : 870x520x780-980mm
- Hi-Low 780-980mm
- Head upward angle : 0-12 degree
- Height Adjustable by Gas Spring
- Easy for move & baby to communicate
- Aluminium Column
- Transparent and Hi-Intensity baby basin
- PP cart Basis
- Four Silent Wheels with cross brakes, dia 100 mm
- Colour Available : Green, Blue & Pink



Case History Trolley (Single & Double)

Description

- ABS Plastic Material, Double rows (60 Layers) Single rows (50 Layers)
- Noiseless Castors
- Available with 2 drawers in double & single with single drawer
- Double dimension 750 x 465 x 810 mm, Weight 15 Kg
- Single dimension 460 x 430 x 810 mm, Weight 9 Kg



Instruments Cabinet

Description

- Made from CRC sheet three side closed
- Front glass door and glass shelves
- Pre-treated and Epoxy Powder Coated
- Approx size : 63" x 24" x 78"

Also available in five shelves Double Door



DELIVERY TABLES

Labour Delivery Room Bed Hydraulic

Description

- Overall Size : 1880L x 660W x 600 - 800H mm
- Tubular frame work mounted on 125 mm dia central locking castors
- Three section top made up of Bakelite Sheet with 3 section mattress
- Trendelenburg / Reverse Trendelenburg positions adjustment through gas spring
- Height adjustment through hydraulic pump
- Back section adjustment through crank mechanism
- Leg end section can slide under the main section
- Provided with a pair of Upholstery Lithotomy positions
- S.S. Telescopic IV Rod
- Detachable S.S. Tray



Labour Delivery Room Bed Electric

Description


- Overall size: 1880L x 660W x 600-800H m
- Tubular frame work mounted on 125 mm dia central locking castors
- Three section top made up of Bakelite Sheet covered with 3 section mattress
- Trendelenburg / Reverse Trendelenburg positions, height adjustment, Back section adjustment through electric actuator system
- Leg end section can slide under the main section
- Provided with a pair of Upholstery Lithotomy positions
- S.S. Telescopic IV Rod
- Detachable S.S. Tray



Patient Shifter

Description

- Maintenance free
- Easy to use
- Easy to carry
- Foldable DELIVERY TABLES
- Unbreakable
- Size : L 72" x w 20"





DELIVERY TABLES

Obstetric Labour Table Telescopic Adjustable Height

Description

- Approx. dimensions: 72"L x 27"W x 24"- 32"H
- Tubular frame work mounted on swiveling castors
- Bakelite sheet top in three sections with 3 section mattress
- Trendelenburg / Reverse Trendelenburg positions, adjustable by gas spring from back end
- Hydraulically height adjustable by foot paddle from side
- Backrest section adjustable crank mechanism from back end
- Leg end section can slide under the main section
- Height adjustable by a pair of knee crutches
- Pre-treated & Epoxy Powder Coated

Also available in S.S. Frame & S.S Top



Labour Table

Description

- Overall approx. size : 72" x 27" x 30"
- Frame work made up of S.S. tube
- Legs fitted with PVC stumps
- 3 section top made up of Bakelite Sheet
- Back rest section adjustable ratchet mechanism from back end
- Leg end section can slide under the main section
- The entire table top gives trendelenburg positions by gas spring from back end
- Drop down S.S. safety side railing are provided on 3 sides of the table
- Height adjustable a pair of foam padded knee crutches
- Patient hand grips made up of S.S. tubes
- Supplied with waste collecting S.S. Bowl on S.S. Bracket and one S.S. I.V. Rod

Also available in M.S. Frame



Gynecological Couch

Description

- Made in 3 sections with independent controls and trendelenburg position
- Adjustable back rest and leg section with mechanism to varying positions capable of tilting for obstetric and gynecological use
- Size of assembled table : approx. 1800L x 550W x 800H mm
- Framework made up of chrome plated rectangular steel tube with upholstered top
- Supplied completely with S.S. bowl and two standard padded leg holders Knockdown construction



BLOOD TRANSFUSION CHAIR

Description

Specification

- **Dimensions**
 - Length: 1800mm
 - Width: 600mm
 - Height: 540mm-830mm
 - Back cushion: 390*220*100 mm
 - Caster: 125mm, PU Caster
- **Material**
 - Chair: Mild Steel
 - Outside: PVC leather
- **Adjustment**
 - Back adjustment: -129C-809C
 - Knee Adjustment
 - Chair Position



Description

- Castor wheels with locks for effortless manual shifting and parking at convenient locations
- High density cushion seating for firm supporting of patient.
- Foldable Handrest for easy patient entry and exit.
- Optional provisions for I-V Rod, Foot rest cushion etc.
- 200 Kg Patient weight capacity.

Specification

- **Dimensions**
 - Length: 1800mm
 - Width: 600mm
 - Height: 540mm-830mm
 - Back cushion: 390*220*100 mm
 - Caster: 125mm, PU Caster
- **Adjustment**
 - Back adjustment: -12 C - 80 C 0
 - Knee Adjustment
 - Chair Position
 - Trendelenburg /Reverse Trendelenburg


DIALYSIS CHAIR



EXAMINATION COUCH

Description

- Back section adjustable to 75 degree
- Powder coated steel frame, firm and anti-rust
- With a shelf to assemble non woven, easy to clean and more anti-infection
- High quantity leather, it is waterproof
- Size: 190 X 80 X 45cm
- Other size and material offered



Description

- Overall approx. dimension : 1800 X 500 X 810 mm
- Construction made from mild steel CRCA tubes and sheets
- Two section mild steel CRCA sheet top
- Headrest adjustable on gas spring
- Three drawers
- Three cabinets with lock
- Complete with B.P. apparatus tray, foam mattress covered with rexine and built in step stool pre-treated and powder coated finish



Description

- Overall approx. dimension : 1800 X 500 X 810 mm
- Construction made from mild steel CRCA tubes and sheets
- Two section mild steel CRCA sheet top
- Headrest adjustable on gas spring
- Three drawers
- Three cabinets with lock
- Complete with B.P. apparatus tray, foam mattress covered with rexine and built in step stool pre-treated and powder coated finish




WARD EQUIPMENTS

Waiting Area Chair

Description

- Perforated sheet
- Galvanized arm rest
- Chrome plated finish

Also available in full Stainless Steel



Wheelchair Fix / Folding

Description

- Folding wheelchair: 600 X 790 X 870 mm
- Folding wheelchair with cushion: 600 X 790 X 870 mm
- Sturdy structures
- High loadbearing
- Durable wheels
- Brake handles provided
- Available in Chrome and powder coated



Patient Transfer Trolley Imported Type

Description

- Frame is made of CRCA Sheet & pipes
- Trolley is mounted over 125 mm diagonally locking castor wheels
- Provision for saline stand
- Provide with swing down polymer moulded safety side rails
- Backrest can be adjusted via gas spring up to 65"
- High-Low can be adjusted via crank mechanism
- Dimension (overall) : 2025 x 760 x 560 - 825 H mm

Also available in ABS Top

