

Linked Page



Certified Register



Wholesale Supplier
Hospital Furniture Economical Quality
A to Z Hospital Equipment and
Surgical Products B to C



 **8275140877**

 **9822688071**

 **info@amramedtech.com**



**Sr No.29/1/2, Adarsh Nagar,
Uruli Devachi, Hadapsar, Pune**
Branches :- All Maharashtra



amramedtech.com

Hospital Beds

AMRA - 001

Hospital Fowler Bed (SS Bows)

Frame work of Strong & precise Steel tubes having 'H' types legs 10cms castors-two with brakes. Having four section top of perforated CRCA sheet duly machine press bent with mattress guards. Backrest knee-rest section manoeuvred by two separate screws from foot end. Laminated panels fitted S.S. head & foot bows. Provision for, four saline locations

Size : 195L X 90W X 60H cms.

Finish: Epoxy Powder Coated.



AMRA - 002

Recovery Bed with Swing Type Side Railings & Castors

1. 60mm x 30mm CRCA tubular frame work
2. Three section perforated CRCA sheet top.
3. Backrest, Trendeleberg & Reverse Trendelenberg operated by screw and lever mechanism with attached SS folding handles.
4. A pair of Swing away side railings.
5. Two locations for |.V. Road
6. Four heavy duty swivel castors of 125 mm dia (two with brake).
7. Pretreated and Epoxy powder coated.



AMRA - 003

Hospital Semi-Fowler Bed (SS Bows)

Frame work of Strong & precise steel tubes having 'H' type legs on PVC stumps. Having two section top of perforated CRCA sheets duly machine press bent, with mattress guards. Laminated panels fitted with S.S. heads & foot bows. Backrest section manoeuvred by crank mechanism from foot end. Provision for locations for |.V.Rod.

Size: 180L X 75W X 60H cms.

Finish: Pre-treated & Epoxy Powder Coated.



AMRA - 004

Hospital Bed Plain

1. Top made up of three piece CRCA perforated sheet double press at four sides with full length stipner support at bottom.
2. S.S.bows at both end and offset legs.
3. Four |.V.locations pre-treated & epoxy powder coated.
(Also available in 50mm X 25mm rectangular CRCA tubes frame work)

Size: 206L X 90W X 60H cms.

Finish: Pre-treated & Epoxy Powder Coated.





AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

Gynac care

AMRA - 005

Gynac Examination Table (Plain)

1. Overall Size: 180L x GOW x 46H cms.
2. Frame work made of Rectangular / Square M.S.Tube.
3. Legs Fitted with PVC Stumps.
4. 5cm Foam cushioned top covered with rexine with perieal cut.
5. Adjustable lithotomy bars with straps.
6. Finish : Epoxy Powder Coated.



AMRA - 006

Obstetric Labour Table Telescopic (Fixed Height)

1. Approx Dimensions :72"L x 27"W x 30"H
2. Tubular Frame work mounted on PVC stumps.
3. Stainless steel top in three sections.
4. Trendlenburg positions adjustable by screw mechanism.
5. Backrest section adjustable on ratchet.
6. Leg end section can slide under the main section.
7. Height adjustable a pair of knee crutches.
8. Pre Treated & Epoxy Powder Coated.
9. Also available in complete stainless Steel.
10. Also available in Gas Spring System.



AMRA - 007

Labour Table (Foot End Detachable)

1. Two section stainless steel top with head low/up position by screw mechanism.
2. Detachable & folding foot end.
3. Available in full stainless steel/stainless top & epoxy powder coated base
4. optional Drainage trolley.



8275140877

9822688071

info@amramedtech.com

amramedtech.com

Pune Maharashtra India

ICCU Beds



AMRA - 008

Height Adjustable ICU Bed

With Collapsible Side Railing.



AMRA - 009

Fixed Height ICU Bed

With Collapsible Side railing



AMRA - 0010

Fixed Height ICU Bed

With Collapsible Side railing

Features

1. 60mm x 30mm CRCA Tubular frame work
2. Four section perforated CRCA sheet top.
3. Polymer Head and Foot Boards.
4. A pair of new design collapsible side railings (1320mm length).
5. S.S Telescopic Saline Rod with four locations.
6. Four heavy duty swivel castors of 125mm Dia. (two with brake)
7. Four side corner buffers.
8. Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg operated by screw & lever mechanism with three separate attached S.S. folding handles.
9. pretreated & epoxy powder coated.
10. Accessories : Mattress |V. Rod, Pillow Rexine Cover and bedsheet



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

Ward Equipments



AMRA - 0011

Baby Cradle

1. (36"x 15"x 39")
2. M.S. Tubular Pipe
3. Wire mesh hanging cradle.
4. Pre-treated and Epoxy powder coated.



AMRA - 0012

Monitor Stand Wall Mounting

AMRA - 0013

Height Adjustable Bed Side Table



AMRA - 0014

Soiled Linen Trolley With Canvas bas

1. Frame work of CRC tube mounted on 100mm castors
2. Superior quality canvas bag
3. pre-treated and Epoxy powder Finish.



AMRA - 0015

Attendant Bed

1. Overall Size : 72"L x 24"W x 18"H
2. Frame work made of Rectangular CRC Tube.
3. 50mm Foam padded top covered with Rexine
4. Legs fitted with PVC Stumps.
5. Finish : Epoxy Powder Coated.



AMRA - 0016

Bed side Screen (3 panels)

1. Size : 16:H x 96:L
2. Tubular frame work mounted on 50mm castors with curtains.
3. Pre-treated and Epoxy Powder coated Finish



8275140877



9822688071



info@amramedtech.com



amramedtech.com

Pune Maharashtra India

Ward / OT Trolley



AMRA - 0017

Z Type trolley (Multipurpose)

1. Three Powder coated shelves
2. Mounted on 2"x1" Tabular Pipe Frame Optional: 50mm Dia Castors



AMRA - 0018

Revolving Stool with S.S Top

1. 1" Tabular Pipe Frame with Rubber Shoe fitted on base.
2. Height Adjustment 18" to 24".



AMRA - 0019

Bed Side Trolley (ECG)

1. Three Powder Coated Shelves .
2. Mounted on Tabular frame Optional : 50mm Dia Castors



AMRA - 0020

Instrument Trolley (2 Shelves)

1. Size 27"L x 18"W x 30"H. Two S.S. Shelves
2. Mounted on 50mm Casters.
3. With Bowl and Bucket Full S.S. Fitted on side
4. Optional : Full S.S Mounted on 100mm Casters.



AMRA - 0021

Instrument Trolley (2 Shelves)

1. Size 27°L x 18°W x 30°H. Two S.S Shelves
2. Mounted on 50mm Casters
3. Optional : S.S. Mounted on 100mm Casters

AMRA - 0022

Medicine Trolley with 4 Drawers

1. Overall Approx Size 76L x 46W x 81H cms.
2. Frame of MS Tubes mounted on 100mm castors.
3. Two SS Shelves with three side railing on top shelf.
4. Two drawer under each shelf.
5. Epoxy Powder Coated Finish.




AMRA - 0023

Cylinder Trolley

1. Tubular Frame work.
2. Two 10 cms castors
3. M.S Sheet base plate
4. Powder Coated Finish
5. For Small Cylinder.





AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

OPD/ Ward Equipments



AMRA - 0024
Bed Side Locker

1. (16"L x 16"W x 32"H)
2. CRCA sheet Constructions.
3. Having one Box & One Drawer.
4. Laminated or S.S Top with 3 sides raise edges.
5. Epoxy Powder Coated.



AMRA - 0025
Bed Side Locker

1. (16"L x 16"W x 32"H)
2. Tubular Pipe Construction fitted PVC stumps assembled to M.S. cupboard
3. Stainless Steel top.
4. Epoxy Powder Coated



AMRA - 0026

Visitor Stool

1. (12"L x 12"/W x 20"H)
2. M.S. Tubular frame fitted with PVC stumps
3. CRC sheet duly double press bent
4. Pre-treated and Epoxy powder Coated



AMRA - 0027

Crash Cart

1. Overall approx size.
2. 940L x 490mmW x 1535mmH.
3. Six removable bins & two polystyrene strong unit with three drawer each.
4. 125mm dia castors, corner Buffer, Oxygen Cylinder Holder, Electric Lamp, J.V Rod Cardiac Massage board & 3 powder coated shelves.
5. Pre-Treated & epoxy Powder coated
6. Also available in SS Frame & J Sheeted



AMRA - 0028

Emergency Drugs Trolley

1. Trolley designed for quick mobility and utility purpose
2. storage for emergency medicines
3. cabinet for I.V storage
4. O2 attachment optional
5. Dynamic design for better and smart utility



AMRA - 0029
Fix Stool S.S Top

1. Multipurpose stool
2. Round S.S Top
3. Top fitted with self rebids for future maintenance.
4. PVC material used in shoes for long term service.



AMRA - 0030

Foot Step Double

1. (18"L x 7"W x 13"H)
2. M.S Tubular Frame fitted
3. CRC Sheet duly double press bent.
4. Pre-treated and epoxy.
5. Powder Coated.



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register


Ward & Neonatal Equipments



AMRA - 0031

Saline Stand

1. M.S Tubular Pipe.
2. Five leg base with grommet fitted
3. S.S adjustable Rod from 137 to 240 cms with double hook
4. Pre Treated & Epoxy powder coated



AMRA - 0032

Phototherapy Unit



AMRA - 0033

Phototherapy Unit with baby bassinet (Double side)



AMRA - 0034

Waiting Bench (Powder Coated)

1. Option 6 feet, 4 feet & 3 feet.
2. With floor adjustment knob.
3. Perforated sheet fitted with screws.
4. Base structure of Tabular steel of 50mm x 25mm
5. Top of Round 25mm Tabular Steel.



AMRA - 0035

Over Bed Table Sunmica Top

1. Size: 42"L x 16"W x 39"H.
2. Tubular Frame Work mounted on 50mm castors.
3. Fixed Laminated top 30" x 16".
4. Pre-treated & Epoxy powder coated.



AMRA - 0036

Baby Resuscitation Trolley

1. Stainless Steel tray in head low position.
2. Three drawer modular system.
3. Adjustable .I.V. Rod.
4. Oxygen Cylinder Cage
5. Bottom tray for mounting suction appaatus
6. Trolley mounted on 100mm Dia
7. Pre-treated & Epoxy powder coated
8. Available in Stainless Steel / Powder coated.



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

Ward Equipments



AMRA - 0037

Kick Bucket Trolley

1. One Stainless Steel bowls Dia 360mm.
2. Available in Stainless Steel Powder Coated



AMRA - 0038

Kick Bucket Trolley

1. Two Stainless Steel bowls Dia 360mm each.
2. Available in Stainless Steel / Powder Coated



AMRA - 0039

Mayo's Trolley (Double Bar) with Height Adjustable

1. Trolley Mounted on 50mm dia. castors.
2. Available in Stainless steel / Powder Coated.



AMRA - 0040

Instrument Cabinet with Glass Shelves

1. Size Option : 5"H x 2"W x 16"D (4 Shelves) Size Option : 6"H x 2.5"W x 16"D (5 Shelves)
2. Rectangular tubular structure.
3. Sheet Panel on back side.
4. Shelves of 5mm thick Glass.
5. Door & two side glass.
6. Available in Stainless Steel / Powder Coated



AMRA - 0041

Mayo's Instrument trolley with S.S tray

1. Tray size S.S. : 56Lx 40 W cms.
2. Rectangular CRC tubes mounted on 5cms castors, height adjustable by side Knob.
3. Pre-treated and Epoxy Powder Coated finish.
4. Also available in S.S. frame & tray.



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C




Certified Register

Health Check

AMRA - 0042

Examination Couch

- Frame : CRCA Rectangular Tubes, CRCA sheet cabinet.
- Four cabinets - One drawer with lock.
- 50 mm thick reversible foam mattress.
- Pre-treated & Epoxy powder coated.
- OPTIONAL : Lock facility for Cabins & 100mm foam mattress



AMRA - 0043

Gynac Examination Table couch with cabinet

- Size : 72"L x 22"W x 33"H
- Optional Table side trolley



AMRA - 0044

Monitor Trolley

- Adjustable 5 shelves
- Mounted on 125 mm Dia. castors. Heavy duty castors 2 with brakes
- 4 point electrical panel board with socket & switch
- Optional Co2 Cylinder cage.
- Available in Stainless Steel / Powder Coated.





AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

Patient Transfer Trolleys



AMRA - 0045

Emergency & Recovery Trolley

- Two sections removable & X-Ray permeable stretcher top.
- Height adjustable & Trendelenburg / Reverse Trendelenburg position by screw & lever mechanism.
- Oxygen cylinder cage.
- LV. Rod.
- S.S. Swing away type side railings.
- Four heavy duty swivel castors of 125mm Dia. two with brakes.
- Buffers at four corners.
- Pre-treated & Epoxy powder coated.



AMRA - 0046

Stretcher Trolley

- Frame work of CRC tubes mounted on 100mm castors.
- Removable stretcher top with chrome plated handle.
- Pre-treated and Epoxy Powder coated finish without Mattress.
- Optional 150mm dia Heavy Duty Wheels.



AMRA - 0047

Stretcher Trolley

- Folding stainless steel top fixed on CRCA tube frame mounted on 100 mm dia heavy duty wheels.
- Pre-treated & epoxy powder coated.
- Also available in Stainless steel.
- Optional 150mm dia Heavy Duty Wheels



AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

Invalid Wheel Chairs & Attendant Equipments



AMRA - 0048

Invalid Wheel chair Folding
(Detachable Footrest)



AMRA - 0049

Invalid Wheel chair with
Camod (Detachable
Footrest)



AMRA - 0050

Recliner Attended Bed

- Multi Functioning Attendant Chair
- The Chair Can Be Used Multi Functioned As Chair or Bed
- Comfortable For Sitting And for Sleep
- Mattress is Mounted With Foam on Both The Sides (Dimensions:- 24"W/ 72"L/ 20H)



AMRA - 0051

Attendant Bed (Sliding)

- with Cushion Top



OT Equipments



AMRA - 0052

Wall Mounted Scrub

- Hand Operated with Option of Soap & Hand Dryer
- Available sizes: L22"/W48"/H48" & L22"/W24"/H48"



AMRA - 0053

Wall Mounted Scrub

- Knee Operated with Option of Soap & Hand Dryer
- Available Sizes: L22"/W48"/H48" & L22"/W24"/H48"



AMRA - 0054

SCRUB (Powder Coated)

- Foot Operated with Option of Soap & Hand Dryer
- Available Sizes : L22"/W48"/H48" & L22"/W24"/H48"



AMRA - 0055

SCRUB Stand Model (Powder Coated)

- Foot Operated with Option of Soap & Hand Dryer
- Available Sizes : L22"/W48"/H48" & L22"/W24"/H48"