

Linked Page



Certified Register



Wholesale Supplier
Operation LED Light Economical Quality
A to Z Hospital Equipment and
Surgical Products B to C



 **8275140877**

 **9822688071**

 **info@amramedtech.com**



**Sr No.29/1/2, Adarsh Nagar,
Uruli Devachi, Hadapsar, Pune**
Branches :- All Maharashtra



amramedtech.com

CONTENTS

UNIQUE FEATURES LED SURGICAL LIGHTS

- Crystal Series
- Star Series
- Moon Series
- Pixel Series

LED EXAMINATION LIGHTS

- Portable Series

1

2

6

9

11

12

UNIQUE FEATURES

Medical Grade LED

We use German medical grade LED in our each & every products. It provides the highest CRI rating, which means in this type of LED you will be able to see true natural colour.

Touch Controller

Touch control panel provides a user friendly interaction and allows to control features of dome such as Intensity, Endo mode & ON /OFF.

Single LED Replacement

In this type of system, we can replace or service each & every individual LED, which is very cost-effective.

Branded Medical Grade LENS

We use Branded European Lenses which are highest in efficiency and Its efficiency remains constant in future.

Aero Dynamic Design

Light Head design enables the ideal laminar airflow on the surgical site which minimize contamination of bacteria.

Seamless Dome


Dome is smooth & continuous, with no apparent Gaps or Spaces between one part & the next.

CRYSTAL

AMRA Crystal 6 Duo

Specification

Dome Diameter	650 mm
Lux Intensity (At 1 Meter)	200000 x 2 LUX
No. of LED	126 x 2 LED
Focus Diameter	150-300 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Focusing	Yes
Power Consumption	80 x 2 Watt
Power Supply	220V / 50HZAC



AMRA Crystal 6

Specification

Dome Diameter	650 mm
Lux Intensity (At 1 Meter)	200000 LUX
No. of LED	126 LED
Focus Diameter	150-300 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Focusing	Yes
Power Consumption	80 Watt
Power Supply	220V / 50HZAC



CRYSTAL

AMRA Crystal 5 DUO

Specification

Dome Diameter	630 mm
Lux Intensity (At 1 Meter)	170000 x 2 LUX
No. of LED	105 x 2 LED
Focus Diameter	150-300 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Focusing	Yes
Power Consumption	60 x 2 Watt
Power Supply	220V / 50HZAC



AMRA Crystal 5

Specification

Dome Diameter	630 mm
Lux Intensity (At 1 Meter)	170000 LUX
No. of LED	105 LED
Focus Diameter	150-300 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Focusing	Yes
Power Consumption	60 Watt
Power Supply	220V / 50HZAC





CRYSTAL

AMRA Crystal DUO 4

Specification

Dome Diameter	610 mm
Lux Intensity (At 1 Meter)	140000 x 2 LUX
No. of LED	84 x 2 LED
Focus Diameter	150-300 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Focusing	Yes
Power Consumption	50 x 2 Watt
Power Supply	220V / 50HZAC



AMRA Crystal 4

Specificatio

Dome Diameter	610 rnm
Lux Intensity (At 1 Meter)	140000 2 LUX
No. of LED	84 LED
Focus Diameter	150-300 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Focusing	Yes
Power Consumption	50 Watt
Power Supply	220V / 50HZAC






CRYSTAL

AMRA Crystal 3 DUO

Specification

Dome Diameter	580 mm
Lux Intensity (At 1 Meter)	120000 x 2 LUX
No. of LED	63 x 2 LED
Focus Diameter	150-300 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Focusing	Yes
Power Consumption	40 x 2 Watt
Power Supply	220V / 50HZAC



AMRA Crystal 3

Specification

Dome Diameter	580 mm
Lux Intensity (At 1 Meter)	120000 LUX
No. of LED	63 LED
Focus Diameter	150-300 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Focusing	Yes
Power Consumption	40 Watt
Power Supply	220V / 50HZAC



STAR

AMRA Star 5 DUO

Specification

Dome Diameter	650 mm
Lux Intensity (At 1 Meter)	200000 x 2 LUX
No. of LED	45 x 2 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Focusing	Yes
Power Consumption	100 x 2 Watt
Power Supply	220V / 50HZAC



AMRA Star 5

Specification

Dome Diameter	650 mm
Lux Intensity (At 1 Meter)	200000 x 2 LUX
No. of LED	45 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Lens Size	36mm (HD Lens)
Focusing	Yes
Power Consumption	100 Watt
Power Supply	220V / 50HZAC





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C




Certified Register

STAR

AMRA Star 4 DUO

Specification

Dome Diameter	610 mm
Lux Intensity (At 1 Meter)	170000 x 2 LUX
No. of LED	36 x 2 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Focusing	Yes
Power Consumption	80 x 2 Watt
Power Supply	220V / 50HZAC



AMRA Star 4

Specification

Dome Diameter	610 mm
Lux Intensity (At 1 Meter)	170000 x 2 LUX
No. of LED	36 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Lens Size	36mm (HD Lens)
Focusing	Yes
Power Consumption	80 Watt
Power Supply	220V / 50HZAC




STAR

AMRA Star 3 DUO

Specification

Dome Diameter	580 mm
Lux Intensity (At 1 Meter)	140000 x 2 LUX
No. of LED	27 x 2 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Focusing	Yes
Power Consumption	60 x 2 Watt
Power Supply	220V / 50HZAC



AMRA Star 3

Specificatio

Dome Diameter	580 mm
Lux Intensity (At 1 Meter)	140000 LUX
No. of LED	27 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Focusing	Yes
Power Consumption	60 Watt
Power Supply	220V / 50HZAC



MOON

AMRA Moon 28 DUO

Specification


Dome Diameter	500 mm
Lux Intensity (At 1 Meter)	140000 x 2 LUX
No. of LED	28 x 2 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Power Consumption	60 x 2 Watt
Power Supply	220V / 50HZAC



AMRA Moon 28

Specification

Dome Diameter	500 mm
Lux Intensity (At 1 Meter)	140000 LUX
No. of LED	28 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Power Consumption	60 Watt
Power Supply	220V / 50HZAC



MOON

AMRA Moon 16 DUO

Specification


Dome Diameter	410 mm
Lux Intensity (At 1 Meter)	100000 x 2 LUX
No. of LED	28 x 2 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Power Consumption	40 x 2 Watt
Power Supply	220V / 50HZAC



AMRA Moon 16

Specification

Dome Diameter	410 mm
Lux Intensity (At 1 Meter)	100000 LUX
No. of LED	16 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Power Consumption	40 Watt
Power Supply	220V / 50HZAC





AMRAMED
TECH Pvt. Ltd.

Wholesale Supplier

A to Z Hospital Equipment and
Surgical Products B to C



Certified Register

PIXEL

AMRA Pixel 300 DUO

Specification


Dome Diameter	575 mm
Lux Intensity (At 1 Meter)	130000 x 2 LUX
No. of LED	27 x 2 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Power Consumption	60 x 2 Watt
Power Supply	220V / 50HZAC



AMRA Pixek 300

Specification

Dome Diameter	575 mm
Lux Intensity (At 1 Meter)	130000 LUX
No. of LED	27 LED
Focus Diameter	100-250 mm
Focus Depth	10-15 cm
Colour Rendering Index	96 RA
Colour Temperature	3200 - 5500 K
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Power Consumption	60 Watt
Power Supply	220V / 50HZAC



CRYSTAL & STAR

AMRA Crystal 21

Specification

Dome Diameter	175 mm
Lux Intensity (At 1 Meter)	50000 LUX
No. of LED	21 LED
Focus Diameter	100-250 mm
Colour Rendering Index	94 RA
Colour Temperature	4500 Kelvin
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	3 Watt Per LED
Power Consumption	12 Watt
Power Supply	220V / 50HZAC

AMRA Star 9

Specification

Dome Diameter	200 mm
Lux Intensity (At 1 Meter)	60000 LUX
No. of LED	9 LED
Focus Diameter	100-250 mm
Colour Rendering Index	94 RA
Colour Temperature	4300 Kelvin
Dimming Range	0-100% (Digital)
LED Life	60000 Working Hours
LED Wattage	5 Watt Per LED
Lens Size	36mm (HD Lens)
Power Consumption	15 Watt
Power Supply	220V / 50HZAC