



TM

AMRAMED[®]
TECH Pvt. Ltd.

Wholesale Supplier

STERILIZERS/DRESSING DRUM/TRAY /FORMALIN CHAMBER/ VIEW BOX
A to Z Hospital Equipment and
Surgical Products **B to C**



8275140877

9822688071

info@amramedtech.com



Sr No.29/1/2, Adarsh Nagar,
Uruli Devachi, Hadapsar, Pune
Branches :- All Maharashtra



amramedtech.com

STERILIZERS

Oxygen Fine Adjustment Valve (With Rota Meter & Humidifier Bottle)

Description:

Electric Instrument Sterilizer is the perfect solution for eliminating bacteria and viruses from your surgical instruments. Constructed with the highest quality stainless steel, this sterilizer is incredibly durable and easy to use. Its advanced design ensures maximum efficiency, making it one of the most dependable tools in your clinic, hospital, nursing home. With its ability to quickly reduce germs and microorganisms, you can rest assured that your medical supplies are germ-free when using this Instrument Sterilizer.



Features:

Its body is made of Medical grade SS.
Electrically operated with immersion heater,
Power source 220v/50hz
With Handles for Lid and Base
Can be used for Medical and Dental both

Size (Inch)	Size (Inch)
8" x 5" x 5"	18" x 8" x 8"
10" x 5" x 5"	20" x 8" x 8"
12" x 6" x 6"	24" x 8" x 8"
14" x 6" x 6"	24" x 12" x 8"
16" x 6" x 6"	24" x 12" x 12"
16" x 8" x 8"	

TRAY

Size (Inch)	Size (Inch)
8' x 6"	12" x 10"
9" x 6"	14" x 10"
10" x 8"	12" x 15"
11" x 7"	18" x 12"
12" x 8"	18" x 12" x 3 Baby Tray

STAINLESS STEEL KIDNEY TRAY

Size (Inch)	Size (Inch)
6"	10"
8"	12"

STAINLESS STEEL CATHETER TRAY WITH COVER

Size (Inch)
8"x3"
16"x4"



DRESSING DRUM

Dressing Drums Seamless (Joint Less) & (Jointed)

Description:

Dressing Drums are used for sterilizing and storing various types of dressings, instruments and surgical requisites that are used by the patients in a disinfected way. These drums are made from high quality hard stainless steel with different layers such that more items can be kept inside it. It also has holes that allow exhausting the steam from outside. Dressing drum is a vital part in a hospital as we believe that only the germ-free environment can bring the patients into good conditions.

Features:

Dressing drums at Maxlife are available in different sizes.

Made of stainless steel 202 & 304 grade.

Products that have been mentioned as Super Deluxe are made up of heavy quality.

Available in both construction jointed and seamless.

Highly reliable and durable.

Can be used to store sterile bathrobes, dressings, and operational linen etc.

Dressing Drums Seamless (Joint Less)



Dressing Drums Seamless (Jointed)



Size D x H (Inch)	Size D x H (Inch)
6" x 4"	6" x 6"
9" x 5"	9" x 9"
11" x 5"	11" x 9"
14" x 5"	14" x 9"
15" x 12"	12" x 10"
9" x 18"	9" x 18"

Size D x H (Inch)	Size D x H (Inch)
6" x 4"	6" x 6"
11" x 5"	8" x 8"
12" x 10"	9" x 9"
14" x 10"	11" x 9"
9" x 6"	15" x 12"
9" x 18"	9" x 18"

FORMALIN CHAMBER

Stainless Steel

- | | |
|--------------------------|-------------|
| 1] Size: 8" X 4" X 4" | Single Tray |
| 2] Size:- 10" X 4" X 4" | Double Tray |
| 3] Size :- 12" X 4" X 4" | Double Tray |
| 4] Size :- 14" X 8" X 6" | Three Tray |
| 5] Size:- 18" X 8" X 6" | Three Tray |
| 6] Size :- 21" X 8" X 6" | Three Tray |
| 7] Size :- 24" X 8" X 8" | Three Tray |
| 7] Size :- 28" X 8" X 8" | Three Tray |



Acrylic

- | | |
|--------------------------|-------------|
| 1] Size: 8" X 4" X 4" | Single Tray |
| 2] Size:- 10" X 4" X 4" | Double Tray |
| 3] Size :- 12" X 4" X 4" | Double Tray |
| 4] Size :- 14" X 8" X 6" | Three Tray |
| 5] Size:- 18" X 8" X 6" | Three Tray |
| 6] Size :- 21" X 8" X 6" | Three Tray |
| 7] Size :- 24" X 8" X 8" | Three Tray |
| 7] Size :- 28" X 8" X 8" | Three Tray |



Foot suction SS

Foot suction PC

Slow suction ELC



X-RAY ILLUMINATOR

**X-Ray Viewer Single Film Compact, Solid
with Backlight of LED Based Panel**

Description

X- Ray Viewer Single Film

- Long life approximate LED's life 50,000 hrs
- Uniform light at the total surface area
- No heat emission
- Wall mounted
- Can be used for tracing on X-Ray
- Auto-sensor
- Overall size: 410 X 490 X 60 mm
- Screen size: 430 X 355 mm



X- Ray Viewer Double Film

- Long life approximate LED's life 50,000 hrs
- Uniform light at the total surface area
- No heat emission
- Wall mounted
- Can be used for tracing on X-Ray
- Auto-sensor
- Overall size: 490 X 765 X 60 mm
- Screen size: 430 X 710 mm

X- Ray Viewer Three Film

- Long life approximate LED's life 50,000 hrs
- Uniform light at the total surface area
- No heat emission
- Wall mounted
- Can be used for tracing on X-Ray
- Separate On/Off switch for each film
- Overall size: 490 X 1120 X 60 mm
- Screen size : 430 X 1060 mm



X- Ray Viewer Four Film

- Long life approximate LED's life 50,000 hrs
- Uniform light at the total surface area
- No heat emission
- Wall mounted
- Can be used for tracing on X-Ray
- Separate On/Off switch for each film
- Overall size: 490 X 1500 X 60 mm
- Screen size : 430 X 1425 mm